

## Aluminum Shower Chair and Commode with Casters

HCPCS: E0163



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- Easy-to-assemble frame
- Padded, open front, vinyl toilet seat (Figure A)
- Aluminum construction is lightweight and durable
- 3" non-skid, rust resistant, swivel casters (2 rear locking) (Figure B)
- Brass push button allows for easy assembly (Figure C)
- Removable, tool-free back
- Plastic arms for added comfort
- Comes with removable 12 qt commode bucket and cover (Figure D)

### PRODUCTS

### Description

11114KD-1	Shower Chair,Alum.w/Cstrs,KD,1/cs
-----------	-----------------------------------

### SPECIFICATIONS

### Description

Carton Shipping Weight	18 lbs.
Height	21"
Seat Dimensions	16" (W) x 16" (D)
Warranty	Limited Lifetime
Weight	14 lbs.
Weight Capacity	300 lbs
Width	18" between arms   21.5" Outside