



PPU For Home & Care Centers

Consumer Model

The "People Picker Upper" (PPU) safely lifts a fallen person from the floor, eliminating the need for a 911 call while transforming a stressful fall into a non-event.

It can raise individuals up to 300 lbs. (136 kg) to a seated position allowing the user to stand up or transfer to and/or from a wheelchair, bed, commode, exercise equipment or even the floor.

Can be self-operated or with the assistance of a caregiver.

Restores dignity and independence.



Self Operation



In the home, individuals can be lifted from the floor after a fall or when getting back up from exercise, therapy, gardening or playing with a family pet.

In Care Centers (assisted living, senior housing, senior centers, hospice and mobile caregivers), allows caregivers to easily & safely recover fallen residents.

Assisted Operation




Seated Transfers



Summary of Features

- Beveled seat-front for easy slide on transfer surface
- Wired remote for convenient operation
- Swinging/Removable Rise Assist Handles for easy, secure transfers
- Securement Belt for added safety
- Powder coated for easy cleaning
- Small footprint for tight spaces
- Lithium-Ion battery powered with automatic charger
- 9' (2.74m) power cord



Physical Specifications

Safe Working Load.....	300 lbs. (136 kg)
Maximum Height.....	50" (1270 mm)
External Width.....	15" (381 mm)
Maximum Seat Height.....	21" (533 mm)
Wheels.....	5" (127 mm)
Total Weight.....	65 lbs. (29.5 kg)

Electrical Specifications

Standard 2-prong 110 VAC Power Cord (USA)
Standard 2-prong 220 VAC Power Cord (International)
Operating Environment: >+5 °C <+40 °C (41° to 104°F)
Lithium-Ion Battery: 24 VDC



Contact@IndeeLift.com

www.IndeeLift.com 1-844-700-5438 International@IndeeLift.com

