

IndeeLift®

The People Picker Upper!

PPU-S For Home and Care Centers

Consumer Model

The "People Picker Upper" standing model (PPU-S) safely lifts a fallen person from the floor, eliminating the need for a call to emergency services while transforming a stressful fall into a non-event.

It can raise individuals up to 300 lbs. (136 kg) to a seated or standing position, and allows for easy transfers to and from a wheelchair, bed, commode, exercise equipment or even the floor.

The standing lift function allows users to be safely raised to access a high bed or an overhead cabinet - a step stool without the step!

Can be self-operated or with the assistance of a caregiver.

Restores dignity and independence.



Assisted Operation



Self Operation




**Standing
Function**



Summary of Features

- Beveled seat-front for easy slide on transfer surface
- Wired remote for convenient operation
- Swinging/Removable Rise-Assist Handles for easy secure transfer
- Dual safety belts for added safety
- Dual safety belts for added safety
- Powder coated for easy cleaning
- Small footprint for tight spaces
- Lithium-Ion battery powered with automatic charger
- 9' (2.74m) power cord



Physical Specifications

Safe Working Load	300 lbs. (136 kg)
Overall Depth	26.5" (673mm)
Maximum Overall Height (raised)	70.75" (1746mm)
Minimum Overall Height (lowered) ...	41.75" (1009mm)
External Width incl Rise Handles.....	23" (584mm)
Maximum Seat Height	30" (762mm)
Seat Width	14" (356mm)
Wheel Diameter.....	5" (127mm)
Unit Weight	78 lbs. (35kg)

Electrical Specifications

Standard 2-prong 110 VAC Power Cord (USA)
 Standard 2-prong 220 VAC Power Cord (International)
 Operating Environment: >+5° C <+40° C (41° to 104°F)
 Lithium-Ion Battery: 24 VDC

Maintains Privacy, Dignity and Independence

