

**UOM** 

1/cs







HCPCS K0824



## **Cirrus Plus HD**

Heavy-Duty Folding Power Wheelchair Rear-Wheel Drive

- Foldable, lightweight frame for easy portability. (Figure A)
- Steel double cross-brace frame.
- Adjustable-tension back upholstery.
- Armrests are padded, adjustable and removable.
- Adjustable-length controller mount.
- Adjustabl-tension back upholstery.
- Swing-away footrests with heel loops.
- Optional elevating leg rests (Item #LELR-TF) sold separately.

- Optional 2" seat depth extension (Model #STDS1818N) sold separately.
- Positioning belt included.
- Padded back upholstery with storage pouch and cushioned seat.
- Programmable controller.
- Flat-free, non-marking tires.
- Easy freewheel operation.
- Manual wheel locks.
- General foam seat cushion included.
- · Anti-tip wheels.

specifications	
HCPCS	K0824
Weight Capacity	400 lbs.
Top Speed	5 mph
Maximum Range	15 miles
Turning Radius	35"
Climbing Angle	6°
Ground Clearance	4"
Dimensions	42.5" (L) × 29" (W)
Floor-to-Seat Height	20"
Seat Dimensions	22" (W) x 16" (D)
Controller	PG VSi 50A
Motors (2)	24V x 320W x 4600 rpm

And tip wheels.	
Batteries (pair)	12V x 36AH (U1)
Charger	4A Offboard
Brakes	Electromagnetic
Caster Wheels	8" x 2" Solid
Drive Wheels	12.5" x 2.25" Flat-free
Base Weight	90 lbs.
Battery Weight (pair)	56 lbs.
Seat Weight	N/A
Total Weight	146 lbs.
Warranty on Frame	Lifetime
Warranty on Electronic Controller/Drive Train Components 14 Months	
Warranty on Batteries	6 Months