



* HCPCS CODING:
E2601 / E2602

General Use



CURVE


*see chart for dimensions

-covers-

STRETCH-AIR FABRIC — designed for air transmission and heat reduction. This two-ply, breathable laminate combines a high-stretch polyester top layer with a waterproof, breathable vapor barrier.

COMFORT-TEK — fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multidimensional stretch and helps alleviate pressure on bony prominences.

MILD ANTI-THRUST SHELF
PREVENTS FORWARD MIGRATION.



MEDIAL AND LATERAL
THIGH SUPPORTS PROMOTE
PROPER LEG ALIGNMENT.

HIGH DENSITY, MOLDED FOAM
INCREASES SURFACE AREA REDUCING
POTENTIAL HIGH-PRESSURE AREAS.

CUSHION STYLE: Non-Skid Bottom

FABRIC TYPE: Stretch-Air

MODEL NUMBER	W X L	FRONT	REAR	ABDUCTOR	LEG TROUGH	WT. LIMIT
CU-SN-1616	16 x 16	2.75"	2.75"	3.5"	2"	250 lbs.
CU-SN-1618	16 x 18	2.75"	2.75"	3.5"	2"	250 lbs.
CU-SN-1620	16 x 20	2.75"	2.75"	3.5"	2"	250 lbs.
CU-SN-1816	18 x 16	2.75"	2.75"	3.5"	2"	300 lbs.
CU-SN-1818	18 x 18	2.75"	2.75"	3.5"	2"	300 lbs.
CU-SN-1820	18 x 20	2.75"	2.75"	3.5"	2"	300 lbs.
CU-SN-2016	20 x 16	2.75"	2.75"	3.5"	2"	350 lbs.
CU-SN-2018	20 x 18	2.75"	2.75"	3.5"	2"	350 lbs.
CU-SN-2020	20 x 20	2.75"	2.75"	3.5"	2"	350 lbs.
CU-SN-2216	22 x 16	2.75"	2.75"	3.5"	2"	400 lbs.
CU-SN-2218	22 x 18	2.75"	2.75"	3.5"	2"	400 lbs.
CU-SN-2220	22 x 20	2.75"	2.75"	3.5"	2"	400 lbs.
CU-SN-2416	24 x 16	2.75"	2.75"	3.5"	2"	400 lbs.
CU-SN-2418	24 x 18	2.75"	2.75"	3.5"	2"	400 lbs.
CU-SN-2420	24 x 20	2.75"	2.75"	3.5"	2"	400 lbs.

CUSHION STYLE: Hook & Loop

FABRIC TYPE: Stretch-Air

MODEL NUMBER	W X L	FRONT	REAR	ABDUCTOR	LEG TROUGH	WT. LIMIT
CU-SV-1616	16 x 16	2.75"	2.75"	3.5"	2"	250 lbs.
CU-SV-1618	16 x 18	2.75"	2.75"	3.5"	2"	250 lbs.
CU-SV-1620	16 x 20	2.75"	2.75"	3.5"	2"	250 lbs.
CU-SV-1816	18 x 16	2.75"	2.75"	3.5"	2"	300 lbs.
CU-SV-1818	18 x 18	2.75"	2.75"	3.5"	2"	300 lbs.
CU-SV-1820	18 x 20	2.75"	2.75"	3.5"	2"	300 lbs.
CU-SV-2016	20 x 16	2.75"	2.75"	3.5"	2"	350 lbs.
CU-SV-2018	20 x 18	2.75"	2.75"	3.5"	2"	350 lbs.
CU-SV-2020	20 x 20	2.75"	2.75"	3.5"	2"	350 lbs.
CU-SV-2216	22 x 16	2.75"	2.75"	3.5"	2"	400 lbs.
CU-SV-2218	22 x 18	2.75"	2.75"	3.5"	2"	400 lbs.
CU-SV-2220	22 x 20	2.75"	2.75"	3.5"	2"	400 lbs.
CU-SV-2416	24 x 16	2.75"	2.75"	3.5"	2"	400 lbs.
CU-SV-2418	24 x 18	2.75"	2.75"	3.5"	2"	400 lbs.
CU-SV-2420	24 x 20	2.75"	2.75"	3.5"	2"	400 lbs.

CUSHION STYLE: Non-Skid Bottom

FABRIC TYPE: Comfort-Tek

MODEL NUMBER	W X L	FRONT	REAR	ABDUCTOR	LEG TROUGH	WT. LIMIT
CU-FN-1616	16 x 16	2.75"	2.75"	3.5"	2"	250 lbs.
CU-FN-1618	16 x 18	2.75"	2.75"	3.5"	2"	250 lbs.
CU-FN-1620	16 x 20	2.75"	2.75"	3.5"	2"	250 lbs.
CU-FN-1816	18 x 16	2.75"	2.75"	3.5"	2"	300 lbs.
CU-FN-1818	18 x 18	2.75"	2.75"	3.5"	2"	300 lbs.
CU-FN-1820	18 x 20	2.75"	2.75"	3.5"	2"	300 lbs.
CU-FN-2016	20 x 16	2.75"	2.75"	3.5"	2"	350 lbs.
CU-FN-2018	20 x 18	2.75"	2.75"	3.5"	2"	350 lbs.
CU-FN-2020	20 x 20	2.75"	2.75"	3.5"	2"	350 lbs.
CU-FN-2216	22 x 16	2.75"	2.75"	3.5"	2"	400 lbs.
CU-FN-2218	22 x 18	2.75"	2.75"	3.5"	2"	400 lbs.
CU-FN-2220	22 x 20	2.75"	2.75"	3.5"	2"	400 lbs.
CU-FN-2416	24 x 16	2.75"	2.75"	3.5"	2"	400 lbs.
CU-FN-2418	24 x 18	2.75"	2.75"	3.5"	2"	400 lbs.
CU-FN-2420	24 x 20	2.75"	2.75"	3.5"	2"	400 lbs.

CUSHION STYLE: Hook & Loop

FABRIC TYPE: Comfort-Tek

MODEL NUMBER	W X L	FRONT	REAR	ABDUCTOR	LEG TROUGH	WT. LIMIT
CU-FV-1616	16 x 16	2.75"	2.75"	3.5"	2"	250 lbs.
CU-FV-1618	16 x 18	2.75"	2.75"	3.5"	2"	250 lbs.
CU-FV-1620	16 x 20	2.75"	2.75"	3.5"	2"	250 lbs.
CU-FV-1816	18 x 16	2.75"	2.75"	3.5"	2"	300 lbs.
CU-FV-1818	18 x 18	2.75"	2.75"	3.5"	2"	300 lbs.
CU-FV-1820	18 x 20	2.75"	2.75"	3.5"	2"	300 lbs.
CU-FV-2016	20 x 16	2.75"	2.75"	3.5"	2"	350 lbs.
CU-FV-2018	20 x 18	2.75"	2.75"	3.5"	2"	350 lbs.
CU-FV-2020	20 x 20	2.75"	2.75"	3.5"	2"	350 lbs.
CU-FV-2216	22 x 16	2.75"	2.75"	3.5"	2"	400 lbs.
CU-FV-2218	22 x 18	2.75"	2.75"	3.5"	2"	400 lbs.
CU-FV-2220	22 x 20	2.75"	2.75"	3.5"	2"	400 lbs.
CU-FV-2416	24 x 16	2.75"	2.75"	3.5"	2"	400 lbs.
CU-FV-2418	24 x 18	2.75"	2.75"	3.5"	2"	400 lbs.
CU-FV-2420	24 x 20	2.75"	2.75"	3.5"	2"	400 lbs.

ALL CUSHIONS ARE AVAILABLE WITH RIGID INSERT AND EXTRA INCONTINENT PROOF LINER