

Jazzy[®]

600 ES

The Jazzy[®] 600 ES provides excellent performance and outstanding stability. Its compact design and standard Active-Trac[®] ATX Suspension allow it to excel indoors and out, making it a superb, all-around power chair choice.

Jazzy[®] 600 ES
shown with deluxe comfort
high-back seat



Pride[®]
Mobility Products Corp.

Specifications

Weight capacity	300 lbs.
Maximum speed	Up to 4 mph
Ground clearance	3"
Turning radius	22"
Overall length	41" without front riggings
Overall width	30.5"
Drive wheels	14" (run flat)
Front wheels	6" solid OMNI-casters
Rear wheels	6" solid OMNI-casters
Suspension	Active-Trac® ATX Suspension
Drivetrain	Two-motor, Mid-Wheel Drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	50 amp, Dynamic Shark
Specialty controls	N/A
Battery charger	Off-board, 5 amp
Per-charge range	Up to 17.1 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	NF-22
Battery weight	38 lbs. each
Base weight	138 lbs.
Standard seat weight	51.5 lbs. (high-back)
Max. seat dimensions	20" x 20"
Warranty	Lifetime limited on frame; 13-month limited on electronics; 13-month limited on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

Models

Jazzy 600 ES 2S-SS | (Code: K0822)
Jazzy 600 ES 2S-C | (Code: K0823)

Features

- OMNI-Casters (nylon, spherical-shaped casters) on front and rear to prevent wheel hang-ups
- Side-mounted, easily accessible freewheel levers
- Active-Trac® ATX Suspension for enhanced performance over more varied terrain
- Easy front access to batteries

Options

- Solid seat pan
- Cup holder
- Swing-away joystick
- Swing-away footrests
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Synergy® Seat
- TRU-Balance® 3 Static Seating

Recommended Lifts

- Backpacker® AVP
- Backpacker® MV
- Backpacker® Plus
- Outlander
- PrideBoom 250



Color Choices



Red