



## Bariatric Aluminum Rehab Shower Commode Chair

with Four Locking Casters

Item #	Color	UOM
NRS185008	White and Blue	1/ea

- Can be positioned over a standard toilet or used as a portable self-contained commode.
- Floor to seat height can be adjusted 4" without tools.
- White aluminum frame is attractive and rust resistant.
- Padded seat back and arms.
- Height adjustable.
- Comes with swing-away footrests.

specifications	
Overall Length (with Riggings)	35.04"
Overall Width	26"
Overall Height	40" - 42"
Seat to Floor Height	20.5" - 22.4"
Seat Dimensions	22" x 17.1"
Backrest Height	20.1"
Net Weight	36.8 lbs.
Gross Weight	43.4 lbs.
Weight Capacity	500 lbs.

accessories		
RTL12085	Biodegradable Sanitary Commode Liner	12 liners/box, 6bx/cs