



Quickie GP/GPV January 2014

EIGP70

Mark For:

Submitting for:

 Quote

 Order

Date: _____

PO#: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: Fax Email

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.

The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME STYLE

<input type="checkbox"/>	024WM23	K0005	Quickie GP-70°	\$1,695.00	<input type="checkbox"/>	024WM23	K0005	Quickie GPV-70° (NA 12"-13")	\$1,695.00
--------------------------	---------	-------	----------------	------------	--------------------------	---------	-------	------------------------------	------------

2. FRAME STYLE & SET UP

Frame Angle

<input type="checkbox"/>	024FA6	60°	NC
<input type="checkbox"/>	024FA7	70°	NC
<input type="checkbox"/>	024FA1	75° (NA with rigid footrest)	NC
<input type="checkbox"/>	024FA8	80° (NA with rigid footrest)	NC

Standard Frame Colors

<input type="checkbox"/>	024C140	Aztec Gold	NEW!	NC
<input type="checkbox"/>	024C2	Black (Not avail as base w/Zebra)		NC
<input type="checkbox"/>	024C141	Black Cherry	NEW!	NC
<input type="checkbox"/>	024C10	Black Opal		NC
<input type="checkbox"/>	024C38	Blue Opal		NC
<input type="checkbox"/>	024C36	Candy Blue		NC
<input type="checkbox"/>	024C14	Candy Purple		NC
<input type="checkbox"/>	024C35	Candy Red		NC
<input type="checkbox"/>	024C142	Electric Blue	NEW!	NC
<input type="checkbox"/>	024C143	Evergreen	NEW!	NC
<input type="checkbox"/>	024C126	Glow		NC
<input type="checkbox"/>	024C144	Green Applie	NEW!	NC
<input type="checkbox"/>	024C147	Hot Sparkle Pink	NEW!	NC
<input type="checkbox"/>	024C148	Mauve Pink	NEW!	NC
<input type="checkbox"/>	024C103	Rootbeer		NC
<input type="checkbox"/>	024C114	Sparkle Silver		NC
<input type="checkbox"/>	024C145	Sunrise Orange	NEW!	NC
<input type="checkbox"/>	024C113	Titanium Color Paint		NC
<input type="checkbox"/>	024C8	Yellow		NC

Matte Frame Colors

<input type="checkbox"/>	024C116	Matte Black		NC
<input type="checkbox"/>	024C151	Matte Black Cherry	NEW!	NC
<input type="checkbox"/>	024C150	Matte Electric Blue	NEW!	NC
<input type="checkbox"/>	024C149	Matte Evergreen	NEW!	NC
<input type="checkbox"/>	024C152	Matte Purple	NEW!	NC

Custom Frame Color

<input type="checkbox"/>	024C112	Desert Camouflage		\$295.00
<input type="checkbox"/>	024C43	Mossy Oak Camouflage		\$295.00
<input type="checkbox"/>	024C42	Stars & Stripes		\$295.00
<input type="checkbox"/>	024C108	Zebra		\$295.00

Xtender Power Assist

Rear seat height of 16.5" w/ Std and Small frames, 16.5"-17.5" with Tall

Higher seat heights available through Built for Me

<input type="checkbox"/>	024XT2	E0986	Xtender Power Assist 4X Multiplier NICKEL	\$6,595.00
--------------------------	--------	-------	---	------------

<input type="checkbox"/>	024XT3	E0986	Xtender Power Assist 4X Multiplier LITHIUM	\$7,595.00
--------------------------	--------	-------	--	------------

Please select one of the following handrim options

<input type="checkbox"/>	024XH3		Plastic Coated	NC
<input type="checkbox"/>	024XH4		Stainless Steel	NC

Xtender - Tire Option

<input type="checkbox"/>	024XW1	E2211 E2212	Pneumatic	NC
--------------------------	--------	----------------	-----------	----

Xtender - Other Options

<input type="checkbox"/>	024XA1		Spare NiMH Power Pack - NICKEL	\$1,395.00
<input type="checkbox"/>	024XA7		Spare NiMH Power Pack - LITHIUM	\$1,795.00

Frame Width

Seat Width 1/2" Narrower

<input type="checkbox"/>	024W3		12"	NC
<input type="checkbox"/>	024W4		13"	NC
<input type="checkbox"/>	024W5		14"	NC
<input type="checkbox"/>	024W6		15"	NC
<input type="checkbox"/>	024W7		16"	NC
<input type="checkbox"/>	024W8		17"	NC
<input type="checkbox"/>	024W9		18"	NC
<input type="checkbox"/>	024W10		19"	\$125.00
<input type="checkbox"/>	024W11	E2201	20"	\$450.00
<input type="checkbox"/>	024W12	E2201	21"	\$450.00
<input type="checkbox"/>	024W13	E2201	22"	\$450.00

Customer Service: 800-333-4000 Fax: 800-300-7502 Email: orders@sunmed.com www.sunrisemedical.com

Specifications are subject to change without notice. Copyright © 2012 Sunrise Medical (US) LLC. (9.09) Quickie GP/GPV 930910 REV. AA

2. FRAME STYLE & SET UP (CONT.)

Sling Depth			
<input type="checkbox"/>	024SD5	12"	NC
<input type="checkbox"/>	024SD6	13"	NC
<input type="checkbox"/>	024SD7	14"	NC
<input type="checkbox"/>	024SD8	15"	NC
<input type="checkbox"/>	024SD9	16"	NC
<input type="checkbox"/>	024SD10	17"	NC
<input type="checkbox"/>	024SD11	18"	NC
<input type="checkbox"/>	024SD72	19"	\$75.00
<input type="checkbox"/>	024SD13	E2203 20"	\$475.00
<input type="checkbox"/>	024SD14	E2203 21"	\$475.00
<input type="checkbox"/>	024SD33	E2204 22"	\$795.00

Frame Length			
<input type="checkbox"/>	024FL30	18-1/4" (Avail. on 12" Sling Depth)	NC
<input type="checkbox"/>	024FL31	19-1/4" (Avail. on 13" Sling Depth)	NC
<input type="checkbox"/>	024FL32	20-1/4" (Avail. on 14", 15" Sling Depth)	NC
<input type="checkbox"/>	024FL33	21-1/4" (Avail. on 14", 15", 16" Sling Depth)	NC
<input type="checkbox"/>	024FL34	22-1/4" (Avail. on 16", 17" Sling Depth)	NC
<input type="checkbox"/>	024FL35	23-1/4" (Avail. on 17", 18" Sling Depth)	NC
<input type="checkbox"/>	024FL36	24-1/4" (Avail. on 18" Sling Depth)	NC
<input type="checkbox"/>	024FL37	25-1/4" (Avail. on 19" Sling Depth)	NC
<input type="checkbox"/>	024FL38	26-1/4" (Avail. on 20" Sling Depth)	NC
<input type="checkbox"/>	024FL39	27-1/4" (Avail. on 21" Sling Depth)	NC
<input type="checkbox"/>	024FL40	28-1/4" (Avail. on 22" Sling Depth)	NC

Frame Height (See Chart)			
<input type="checkbox"/>	024FH2	Short (front 9.25"-rear 8.25") NA with Xtender Power Assist	\$110.00
<input type="checkbox"/>	024FH5	Standard (front 10" - rear 9")	NC
<input type="checkbox"/>	024FH7	Tall (front 11" - rear 10")	\$110.00

3. SEATING

Back Type			
<input type="checkbox"/>	024BT3	Non-Folding Backrest (Fixed) <i>NA with push handles or adj. upholstery</i>	
<input type="checkbox"/>	024BB1	8" Back Height	NC
<input type="checkbox"/>	024BB2	9" Back Height	NC
<input type="checkbox"/>	024BB3	10" Back Height	NC
<input type="checkbox"/>	024BB4	11" Back Height	NC
<input type="checkbox"/>	024BB5	12" Back Height	NC
<input type="checkbox"/>	024BB6	13" Back Height	NC
<input type="checkbox"/>	024BB7	14" Back Height	NC
<input type="checkbox"/>	024BB8	15" Back Height	NC
<input type="checkbox"/>	024BB9	16" Back Height	NC
<input type="checkbox"/>	024BT1	Folding Backrest (Adjustable)	
<input type="checkbox"/>	024BA1	11" - 14" Back Height	\$65.00
<input type="checkbox"/>	024BA2	12" - 16" Back Height	\$65.00
<input type="checkbox"/>	024BA3	14" - 18" Back Height	\$65.00
<input type="checkbox"/>	024BA4	16" - 20" Back Height	\$65.00

Back Upholstery			
<input type="checkbox"/>	024BU1	Standard	NC
<input type="checkbox"/>	024BU8	Adjustable	\$175.00
<input type="checkbox"/>	024BU16	Reinforced (NA w/ adj. upholstery)	\$30.00
<input type="checkbox"/>	024BU14	Omit Back Upholstery	-\$15.00

Backrest Options			
<input type="checkbox"/>	024BT33	Folding w/ Integral Push Handles	NC
<input type="checkbox"/>	024A85	Bolt-On Push Handle (Folding only) NA with Xtender Power Assist	\$115.00
<input type="checkbox"/>	024SR2	Regular Stroller Handle NA with amputee axle plate or Xtender Power Assist	\$125.00
<input type="checkbox"/>	024SR3	Tall Stroller Handle NA with amputee axle plate or Xtender Power Assist	\$125.00

JAY Back (Optional)			
<input type="checkbox"/>	024JB01	J3 Back (See J3 Order Form)	Price List

Cushion			
<input type="checkbox"/>	024CU9	Reinforced Seat Sling	\$30.00
<input type="checkbox"/>	024CU1	K0669 2" Cushion	\$95.00
<input type="checkbox"/>	024CU2	Omit Cushion	NC
<input type="checkbox"/>	024CU4	Omit Cushion and Seat Sling	-\$15.00

JAY Bundled Cushions - Optional			
<i>To order any JAY cushions on-chair, please see JAY Price List. Cushion name or part number must be filled out with desired width & depth.</i>			
<input type="checkbox"/>	024JC20	JAY Bundled Cushion Cushion Name or Part #: _____ Width: _____ Depth: _____	Price List

4. ARMREST

Single Post				
<input type="checkbox"/>	024AA9	E0973	Single Post Ht. Adjust. Low <i>Must be removed to fold back down</i>	\$265.00
<input type="checkbox"/>	024AM1		Desk Length - Classic (10")	NC
<input type="checkbox"/>	024AM2		Full Length - Contour (14")	NC
<input type="checkbox"/>	024AA8	E0973	Single Post Ht. Adjust. Std <i>Must be removed to fold back down</i>	\$265.00
<input type="checkbox"/>	024AM1		Desk Length - Classic (10")	NC
<input type="checkbox"/>	024AM2		Full Length - Contour (14")	NC

Padded Swing Away Armrest				
<i>All Padded Swing-Away NA with amputee axle plate or Xtender Power Assist</i>				
<input type="checkbox"/>	024AA19	E0973	Padded S-A Medium	\$195.00
<input type="checkbox"/>	024AA16	E0973	Padded S-A Adult	\$195.00

5. FOOTREST

Footplates

Heel loop or leg strap not included, See Footplate Options Section

1. GP FOOTREST OPTIONS

<input type="checkbox"/>	024F56	K0040	Rigid NA w/ Ext. Tubes; NA with 75° or 80° frame angle	NC
<input type="checkbox"/>	024F50	K0040	Platform NA w/ Ext. Tubes	\$95.00
<input type="checkbox"/>	024F1	K0040	Composite Footplate	NC
<input type="checkbox"/>	024F3	K0040	Adult Angle Adj.	\$185.00
<input type="checkbox"/>	024F48	K0040	Platform Flip Up	\$190.00

2. GPV FOOTREST OPTIONS

<input type="checkbox"/>	024F56	K0040	Rigid NA w/ Ext. Tubes; NA with 75° or 80° frame angle	NC
<input type="checkbox"/>	024F50	K0040	Platform NA w/ Ext. Tubes	\$95.00
<input type="checkbox"/>	024F48	K0040	Platform Flip Up	\$190.00

Footrest Extension Tubes

<input type="checkbox"/>	024E1		Short (14" to 16-1/2")	NC
<input type="checkbox"/>	024E2		Medium (16-1/2" to 19")	NC
<input type="checkbox"/>	024E5		Long (19" to 21-1/2")	NC
<input type="checkbox"/>	024E16	K0037	High Mount Only avail. with adult angle adjustable footplates	\$55.00

Footplate Options (Ext Tube & Direct Mount)

<input type="checkbox"/>	024HL1	E0951	Heel Loop	\$30.00
<input type="checkbox"/>	024HL5	K0038	Leg Strap	\$35.00
<input type="checkbox"/>	024HL4		Omit Heel Loop	NC
<input type="checkbox"/>	024HL6		Omit Leg Strap	NC
<input type="checkbox"/>	024F7		Omit Footplate & Ext Tube	-\$15.00
<input type="checkbox"/>	024IG2		Impact Guard Neoprene	\$85.00

6. CASTERS

Caster Wheels

<input type="checkbox"/>	024CW1		3" Micro	NC
<input type="checkbox"/>	024CW3		4" Low Profile Polyurethane	NC
<input type="checkbox"/>	024CW31		4" Micro Lighted	\$70.00
<input type="checkbox"/>	024CW6		5" Soft Roll Semi-Pneumatic	\$85.00
<input type="checkbox"/>	024CW5		5" Low Profile Polyurethane	NC
<input type="checkbox"/>	024CW7		5" x 1.5" Semi-Pneumatic only avail. For use in bottom fork holes	\$85.00
<input type="checkbox"/>	024CW33	E2214 E2215	6" Pneumatic	\$75.00
<input type="checkbox"/>	024CW34		6" Polyurethane	NC
<input type="checkbox"/>	024CW35		6" x 1.5" Semi-Pneumatic only avail. For use in bottom fork holes	\$85.00
<input type="checkbox"/>	024CW21	E2214 E2215	8" Pneumatic Only available on 60° frame option	\$75.00

Caster Wheels - Continued

<input type="checkbox"/>	024CW32		8" Polyurethane Only available on 60° frame option	NC
<input type="checkbox"/>	024CW36	E2214 E2215	8" x 2" Pneumatic Only available on 60° frame option	\$145.00
<input type="checkbox"/>	024CW37	E2214 K0108	8" x 2" Pneumatic w/Airless Insert Only available on 60° frame option	\$225.00

Caster Fork Options

<input type="checkbox"/>	024CS5		3/4" Longer Fork Stem Bolt	NC
<input type="checkbox"/>	024CS9		1-1/2" Longer Fork Stem Bolt	NC
<input type="checkbox"/>	024CS3	K0073	Caster Pin Locks	\$195.00
<input type="checkbox"/>	024CW12		Omit Casters	NC

7. AXLE & WHEEL OPTIONS

Axle Plate

N/A w/ Power Assist

<input type="checkbox"/>	024AL2		Standard Adjustable	NC
<input type="checkbox"/>	024AL13		Double Slotted Adjustable	\$60.00
<input type="checkbox"/>	024AL6	E0959	Amputee Adjustable NA with padded swing-away armrests or stroller handles	\$175.00

Axle

<input type="checkbox"/>	024AX1		Quick Release Steel	NC
<input type="checkbox"/>	024AX11		Quad Quick Release	\$90.00

Rear Anti-Tip

<input type="checkbox"/>	024AT1	E0971	Rear Anti-Tip	\$135.00
--------------------------	--------	-------	---------------	----------

Rear Wheel Type

Spinergy NA w/ Xtender, spoke guards. Contact Built 4 Me for Spinergy Color Options

<input type="checkbox"/>	024RW8		22" Mag (NA with Xtender Power Assist)	NC
<input type="checkbox"/>	024RW9		22" Spoke (NA with Xtender Power Assist)	\$200.00
<input type="checkbox"/>	024RW1		24" Mag	NC
<input type="checkbox"/>	024RW2		24" Spoke	NC
<input type="checkbox"/>	024RW11		26" Spoke (NA with Xtender Power Assist)	\$65.00
<input type="checkbox"/>	024RW42		24" Spinergy (540) SPOX Black	\$695.00
<input type="checkbox"/>	024RR25		Omit Rear Wheel and Axle Intended Size for Set Up :	-\$55.00

Rear Wheel Tire

<input type="checkbox"/>	024RT1	E2211 E2212	Pneumatic	NC
<input type="checkbox"/>	024RT2	E2211 E2213	Pneumatic Airless Insert NA with 26" wheels	\$80.00

Rear Wheel Tire (Continued)

<input type="checkbox"/>	024RT3	E2218	Full Poly (24" Only) NA with 26" wheels	NC
<input type="checkbox"/>	024RT6	E2211 E2212	High-Pressure Clincher NA w/ 26", N/C w/ Performance	\$95.00

Hand Rim

<input type="checkbox"/>	024HR1		Aluminum Anodized	NC
<input type="checkbox"/>	024HR4		Plastic Coated -NA on all rims	\$95.00
<input type="checkbox"/>	024HR7		Omit Hand Rim	NC

Projection Hand Rim

Projections NA w/ plastic coated handrims

<input type="checkbox"/>	024HR8		8 Oblique (24"/26")	\$240.00
--------------------------	--------	--	---------------------	----------

Wheel Lock

<input type="checkbox"/>	024WL1		High Mount Push	NC
<input type="checkbox"/>	024WL2		High Mount Pull	NC
<input type="checkbox"/>	024WL5		High Mount Scissor NA w/ extension handles, grade aides	NC
<input type="checkbox"/>	024WL12		Do Not Mount Wheel Locks	NC
<input type="checkbox"/>	024WL10		Omit	NC

Extension Handle

<input type="checkbox"/>	024EH1	E0961	6" Pair	\$75.00
--------------------------	--------	-------	---------	---------

Grade Aid

<input type="checkbox"/>	024GA1	E0974	Grade Aid NA w/ poly, high-pressure clincher tires, or Xtender	\$150.00
--------------------------	--------	-------	--	----------

8. KEY OPTIONS & ACCESSORIES

Backpack			
<input type="checkbox"/>	024A1		Quickie Backpack Black \$55.00
<input type="checkbox"/>	024A2		Quickie Kids Backpack Black \$50.00

Seat Pouch			
<input type="checkbox"/>	024A8		Quickie Seat Pouch Black \$45.00

Crutch Holder			
<input type="checkbox"/>	024A53	E2207	Left Side \$195.00
<input type="checkbox"/>	024A52	E2207	Right Side \$195.00

Extra Leg Strap			
<input type="checkbox"/>	024A55	K0038	Extra Leg Strap Black \$30.00

Miscellaneous			
<input type="checkbox"/>	024A83		Tool Kit \$45.00

Positioning Belts			
<input type="checkbox"/>	024A142	E0978	Belt 1.5" Auto Buckle w/Clamps \$70.00
<input type="checkbox"/>	024A12	E0978	Belt 2" Aircraft Buckle w/Clamps \$110.00
<input type="checkbox"/>	024A13	E0978	Belt 2" Aircraft Buckle Pad w/Clamp \$120.00

Spoke Guards			
<i>NA w/ Xtender Power Assist; Only available w/ Mag and Lite Spoke wheels</i>			
<input type="checkbox"/>	024A23	K0065	Black Spoke Guards \$115.00
<input type="checkbox"/>	024A24	K0065	Clear Spoke Guards \$115.00

Tray Table			
<input type="checkbox"/>	024A31	E0950	Tray Table \$250.00

Side Guard			
<input type="checkbox"/>	024SG3	K0108	Plastic Regular \$125.00
<input type="checkbox"/>	024SG4	K0108	Plastic Kids \$125.00

9. MISCELLANEOUS

GP Series Rear STF Chart

**All measurements given are an estimate and have been rounded to the nearest .5"*

SHORT Frame (9 1/4" Front / 8 1/4" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	14.5	15	15.5	16	16.5	17
22"	Mag / Spoke	Pneum / Insert	14.5	15	16	16	17	17.5
24"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17	17.5
24"	Spoke / Spinerger	Clincher	15	16	16.5	16.5	17	18
24"	Mag / Spoke / Spinerger	Pneum / Insert	15.5	16	17	17	17.5	18
26"	Spoke / Spinerger	Clincher	16	17	17.5	17.5	18	19
26"	Spoke / Spinerger	Pneumatic	16.5	17	17.5	18	18.5	19

TALL Frame (11" Front / 10" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	16	16.5	17.5	17.5	18	19
22"	Mag / Spoke	Pneum / Insert	16.5	17	17.5	18	18.5	19
24"	Mag / Spoke	Full Poly	17	17.5	18	18.5	18.5	19.5
24"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
24"	Mag / Spoke / Spinerger	Pneum / Insert	17	17.5	18.5	19	19	20
26"	Spoke / Spinerger	Clincher	18	18.5	19	19.5	20	20.5
26"	Spoke / Spinerger	Pneumatic	18	19	19.5	19.5	20	21

STANDARD Frame (10" Front / 9" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
22"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17.5	18
22"	Mag / Spoke	Pneum / Insert	15.5	16	17	17	17.5	18
24"	Mag / Spoke	Full Poly	16	16.5	17	17	18	18.5
24"	Spoke / Spinerger	Clincher	16	16.5	17	17.5	18	18.5
24"	Mag / Spoke / Spinerger	Pneum / Insert	16	16.5	17.5	17.5	18	18.5
26"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
26"	Spoke / Spinerger	Pneumatic	17	18	18.5	18.5	19.5	20

Maximum Dump 2 3/4"
 *** Custom Frame Height Through Configurator - REAR 6 1/2 - 13" (Max frame squeeze = 3 1/2") Calculate Accordingly

GP Series Front STF Chart

**All measurements given are an estimate and have been rounded to the nearest .5"*

SHORT Frame (9 1/4" Front / 8 1/4" Rear)										
Fork	Spacer	Fork Hole	3"	4" Poly	5" Poly	6" Poly	6" Pne	8" Poly	8" Pne	8x2
			4" Lighted	5" Soft Roll	6" x 1.5					
4 1/4	STD	Bottom	15.50	16.00	16.50	17.00	17.00	x	x	x
4 1/4	3/4	Bottom	16.50	17.00	17.00	17.50	18.00	x	x	x
4 1/4	1 1/2	Bottom	17.00	17.50	18.00	18.50	18.50	x	x	x
5 1/4	STD	Top	x	x	16.50	17.00	17.00	x	x	x
5 1/4	STD	Bottom	16.50	17.00	17.5	18.00	18.00	19.00	20.5	x
5 1/4	3/4	Top	x	x	17.50	18.00	18.00	x	x	x
5 1/4	3/4	Bottom	17.50	18.00	18.50	18.50	19.00	19.50	20.00	x
5 1/4	1 1/2	Top	x	x	18.00	18.50	18.50	x	x	x
5 1/4	1 1/2	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50	x
6	STD	Top	x	x	17.50	18.00	18.00	18.50	19.00	x
6	STD	Bottom	17.50	18.00	18.50	19.00	19.00	19.50	20.00	x
6	3/4	Top	x	x	18.00	18.50	18.50	19.50	19.50	x
6	3/4	Bottom	18.00	18.50	19.00	19.50	19.50	20.50	20.50	x
6	1 1/2	Top	x	x	19.00	19.50	19.50	20.00	20.50	x
6	1 1/2	Bottom	19.00	19.50	20.00	20.50	20.50	21.00	21.50	x
7	STD	Top	x	x	18.50	19.00	19.00	19.50	20.00	x
7	STD	Bottom	18.50	19.00	19.50	20.00	20.00	20.50	21.00	x
7	3/4	Top	x	x	19.00	19.50	19.50	20.50	20.50	x
7	3/4	Bottom	19.00	19.50	20.00	20.50	20.50	21.50	21.50	x
7	1 1/2	Top	x	x	20.00	20.50	20.50	21.00	21.50	x
7	1 1/2	Bottom	20.00	20.50	21.00	21.00	21.50	22.00	22.50	x
8x2	STD	Top	x	x	x	x	x	x	x	19.00
8x2	STD	Bottom	x	x	x	x	x	x	x	20.00
8x2	3/4	Top	x	x	x	x	x	x	x	19.50
8x2	3/4	Bottom	x	x	x	x	x	x	x	20.50
8x2	1 1/2	Top	x	x	x	x	x	x	x	20.50
8x2	1 1/2	Bottom	x	x	x	x	x	x	x	21.00

STANDARD Frame (10" Front / 9" Rear) ADD 3/4" to each measurement above

TALL Frame (11" Front / 10" Rear) ADD 1 3/4" to each measurement above

*** Custom Frame Height through Configurator - Front 9 - 13" - Calculate Accordingly (Minimum Front STF = 15 1/2" / Max = 26")

Frame Dimensions

