






Pride[®]'s Classic Collection delivers superb value and excellent comfort with features like an innovative arm design and overstuffed fiber-filled backs. In addition, Classic Collection Lift Chairs feature patented, quiet and smooth operation and stylish designs to add to any home décor.





1-800-800-8586 (U.S.) 1-888-570-1113 (Canada)

COLLECTION

www.pridemobility.com





CLASSIC COLLECTION MODELS



2-Position Medium to Large Button-back Design 5'2" - 5'10"

CL-15 3-Position Medium to Large Button-back Design

5'2" - 5'10"

2-Position Medium to Large Sewn Pillow-back Design 5'4" - 6'0" **CL-30**

3-Position Medium to Large Chaise Lounger Split-back Design 5'4" - 6'0"

SPECIFICATIONS

	CL-10	CL-15	CL-20	CL-30
Seat To Floor:	20"	20"	20"	20"
Seat Depth:	20"	20"	20"	20"
Seat Width:	20"	20"	20"	22"
Top of Back to Seat:	25"	27"	24.5"	26"
Distance From Wall:	12"	18"	11"	12"
Overall Chair Width:	33.5"	33.5"	32.5"	34"
Chair Weight:	110 lbs.	119 lbs.	114 lbs.	118 lbs.

FABRIC SELECTION STANDARD FABRICS



BENEFITS

- · Engineered furniture grade laminate/ hardwood frame
- Quiet & smooth lift system
- Heavy-duty lift actuator & scissor mechanism
- · Standard head and arm covers
- Fire retardant (California specifications; bulletin 117)
- · Standard, large dual pockets
- Removable back no tools required
- 375 lbs. weight capacity on CL-20 & CL-30; 325 lbs. weight capacity on CL-10 & CL-15
- Integrated battery backup in case of power failure¹
- Transformer features self-diagnostic electronics with 12' cord
- Lifetime warranty on steel lift frame
- 2-year limited warranty

OPTIONS*

- Heat & Massage (CL-20 & CL-30 only)
- Head Pillow
- Lumbar pillow

* Additional charges apply for all options

1. Batteries not included.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice



Additional lead-time may apply for some options.

www.pridemobility.com