



REBEL

NEW!



CARTON SIZE MAKES THE REBEL  
PERFECT FOR RETAIL!



PACKAGED IN SMALL CARTON THAT  
CAN BE SHIPPED UPS!

QUICK  
RELEASE  
REAR AXLE



FOLD  
DOWN  
BACK



TOOL FREE  
REMOVABLE  
FRONT  
CASTERS



## Rebel Wheelchair

Single Axle

### FEATURES:

- Quick release rear wheels allow for easy storage or transport.
- Fold-down back and folding frame allows transport in almost any trunk.
- Easily removable swing away footrests are also convenient for storage.
- Flip up desk style armrests make transfers easy and accessibility better.
- Lightweight and extremely maneuverable, easy to navigate tight areas.
- Large accommodating seating surface.
- Limited Lifetime Warranty.

1-Pack Available for All Wheelchair Models

Due to product improvements and changes, all images in this catalog are subject to change.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

**Rebel – SINGLE AXLE**

|                |  |      |
|----------------|--|------|
| RTLREB18DDA-SF | 18", Flip Back Desk Arm, Swing-away Footrest | 1/cs |
|----------------|--|------|

**REPLACEMENT FRONT RIGGINGS**

|           |                     |        |
|-----------|---------------------|--------|
| RTLREB-12 | Swing-Away Footrest | 1pr/cs |
|-----------|---------------------|--------|

**ARMS**

|           |                           |      |
|-----------|---------------------------|------|
| RTLREB-20 | Right, Flip Back Desk Arm | 1/cs |
|-----------|---------------------------|------|

|           |                          |      |
|-----------|--------------------------|------|
| RTLREB-19 | Left, Flip Back Desk Arm | 1/cs |
|-----------|--------------------------|------|

**WHEEL LOCKS (Push to Lock)**

|            |  |      |
|------------|--|------|
| STDS8M241R | Right Wheel Lock, Push to Lock, Detachable Arm | 1/ea |
|------------|--|------|

|            |   |      |
|------------|---|------|
| STDS8M241L | Left Wheel Lock, Push to Lock, Detachable Arm | 1/ea |
|------------|---|------|

**SEAT BELTS**

|                           |      |
|---------------------------|------|
| STDS850 (Auto Clasp Type) | 1/bx |
|---------------------------|------|

|                                |      |
|--------------------------------|------|
| STDS851 (Velcro® Type Closure) | 1/bx |
|--------------------------------|------|

**BARIATRIC SEAT BELTS**

|                           |      |
|---------------------------|------|
| STDS855 (Auto Clasp Type) | 1/bx |
|---------------------------|------|

|                                |      |
|--------------------------------|------|
| STDS856 (Velcro® Type Closure) | 1/bx |
|--------------------------------|------|

**UNIVERSAL OXYGEN "E" CYLINDER HOLDER**

|         |  |
|---------|--|
| STDS803 |  |
|---------|--|

|                                   |  |
|-----------------------------------|--|
| STDS804 with I.V. Pole Attachment |  |
|-----------------------------------|--|

**WHEEL LOCK EXTENSIONS**

|         |      |
|---------|------|
| STDS801 | 1/ea |
|---------|------|

**ANTI FOLD BAR**

|         |      |
|---------|------|
| STDS806 | 1/cs |
|---------|------|

**WHEELCHAIR CANE/CRUTCH HOLDER**

|          |      |
|----------|------|
| STDS1034 | 1/bx |
|----------|------|

**HEEL LOOPS**

|         |        |
|---------|--------|
| STDS831 | 1pr/bx |
|---------|--------|

**ACCESSORIES ON PAGES 57-78**

1-Pack Available for All Wheelchair Models

**SPECIFICATIONS**

| Seat Width | Overall Width (open) | Overall Width (closed) | Overall Height | Seat To-Floor Height | Back Height | Width Seat | Depth Seat | Overall Length (with riggings) | Product Weight | Carton Shipping Weight | Weight Limit |
|------------|----------------------|------------------------|----------------|----------------------|-------------|------------|------------|--------------------------------|----------------|------------------------|--------------|
| 18"        | 25.5"                | 10.5"                  | 35"            | 19.5"                | 17"         | 18"        | 16"        | 43"                            | 37 lbs.        | 42 lbs.                | 250 lbs.     |

Due to product improvements and changes, all images in this catalog are subject to change.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.