Profile Mid-Wheel Drive Lift AL580 / AL580XL





| Specifications | AL580 | AL580XL Extra Large |
|----------------------|--|--|
| Accommodates | Mid-wheel power chairs with 10" drive wheels | Mid-wheel power chairs with 14" drive wheels |
| Platform | 16.5" x 27" | 17.5" x 27" |
| Lifting Capacity | 350 lbs. (160 kg) | 350 lbs. (160 kg) |
| Installed Weight | 63 lbs. (31 kg) | 63 lbs. (31 kg) |
| Securement | Padded hold-down arm & wheel cradle | Padded hold-down arm & wheel cradle |
| Hitch | Class II or greater | Class II or greater |
| Included Accessories | AL123 Class III hitch adapter / AL107 hitch post | AL123 Class III hitch adapter / AL107 hitch post |
| Options | Battery pack, Swing-Away, Power chair covers | Battery pack, Swing-Away, Power chair covers |
| Warranty | 3-year transferable | 3-year transferable |







