

## Viper Wheelchair

Deluxe High-Strength, Lightweight, Dual-Axle



HCPCS  
K0004

- Viper weighs 32 lbs. (Excluding front rigging)
- Easier to propel and transport.
- Frame style eliminates seat guides and allows for custom backs and accessories.
- Back height adjustable from 18"– 20".
- Carbon steel black powder-coated frame double coated for an attractive, chip-proof, maintainable finish.
- Durable, flame retardant nylon upholstery resists mildew and bacteria.
- Composite Mag-style wheels with aluminum hand rims are lightweight and maintenance free.
- Urethane rear tires and casters offer superior performance, smoother ride and are lightweight.
- 8" Front casters have three height adjustments.
- Nylon, padded armrests provide added patient comfort.
- Flip-back arms allow for easy transfer.
- Composite steel lined front receiver is angled so the arm will lock into place every time.
- The front mounted release lever can be flipped up or down to unlock flip-up arm, depending on the patient's needs. (Figure A - B)
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability.
- Dual axle provides easy transition of seat height to hemi level.
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings.
- Composite foot plates with heel loops are durable and lightweight.
- Comes standard with flip-up anti-tippers.
- Comes with push-to-lock wheel locks.
- Available in 12" and 14" seat widths.



Item #	Seat Width	Armrest	Product Front Riggings	UOM
L412DDA-SF	12"	Flip-Back Desk Arm	Swing-away Footrests	1/cs
L412DDA-ELR	12"	Flip-Back Desk Arm	Swing-away, Elevating Leg rests	1/cs
L414DDA-SF	14"	Flip-Back Desk Arm	Swing-away Footrest	1/cs
L414DDA-ELR	14"	Flip-Back Desk Arm	Swing-away, Elevating Leg rests	1/cs
L416DDA-SF	16"	Flip-Back Desk Arm	Swing-away Footrests	1/cs
L416DDA-ELR	16"	Flip-Back Desk Arm	Swing-away, Elevating Leg rests	1/cs
L416DFA-SF	16"	Flip-Back Full Arm	Swing-away Footrests	1/cs
L416DFA-ELR	16"	Flip-Back Full Arm	Swing-away, Elevating Leg rests	1/cs
L416ADDA-SF	16"	Flip-Back/Adj. Height Desk Arm	Swing-away Footrests	1/cs
L416ADDA-ELR	16"	Flip-Back/Adj. Height Desk Arm	Swing-away, Elevating Leg rests	1/cs
L418DDA-SF	18"	Flip-Back Desk Arm	Swing-away Footrests	1/cs
L418DDA-ELR	18"	Flip-Back Desk Arm	Swing-away, Elevating Leg rests	1/cs
L418DFA-SF	18"	Flip-Back Full Arm	Swing-away Footrests	1/cs
L418DFA-ELR	18"	Flip-Back Full Arm	Swing-away, Elevating Leg rests	1/cs
L418ADDA-SF	18"	Flip-Back/Adj. Height Desk Arm	Swing-away Footrests	1/cs
L418ADDA-ELR	18"	Flip-Back/Adj. Height Desk Arm	Swing-away, Elevating Leg rests	1/cs
L420DDA-SF	20"	Flip-Back Desk Arm	Swing-away Footrests	1/cs
L420DDA-ELR	20"	Flip-Back Desk Arm	Swing-away, Elevating Legrests	1/cs
L420DFA-SF	20"	Flip-Back Full Arm	Swing-away Footrests	1/cs
L420DFA-ELR	20"	Flip-Back Full Arm	Swing-away, Elevating Leg rests	1/cs

**specifications**

Width Open (12" Seat   14" Seat)	20"   22.5"	Length with Riggings with 12" and 14" Seats	36"
Width Open (16" Seat   18" Seat)	22.5"   24.5"	Length with Riggings with 16", 18" and 20" Seat	42"
Width Open with 20" Seat	28"	Product Weight with 12" Seat	30 lbs.
Width Closed	12"	Product Weight with 14" Seat	31 lbs.
Height with 12" and 14" Seats	32"	Product Weight with 16" Seat	31.5 lbs.
Height with 16", 18" and 20" Seat	36"	Product Weight with 18" Seat	32 lbs.
Seat-To-Floor Height (12" - 14") (16" - 20")	13.5" - 15.5"   17.5" - 19.5"	Product Weight with 20" Seat	34 lbs.
Back Height	16" - 18"	Carton Shipping Weight	39 lbs. or less
Seat Depth with 12" Seat	12"	Weight Limit (with 12" and 14" seat)	250 lbs.
Seat Depth with 14" Seat	14"	Weight Limit (with 16", 18" and 20" seat)	300 lbs.
Seat Depth with 16", 18" and 20" Seat	16"	Warranty	Limited Lifetime

**accessories**

**Amputee Adapter**

STDS826	Amputee Adapter	1/pr
STDS826C	Amputee Adapter	1/pr

**Universal Forward Stabilizer**

STDS860	Universal Forward Stabilizer	1/pr
---------	------------------------------	------

**Anti-Tippers**

STDS807	Anti-Tippers With Wheels	1/pr
STDS814	Anti-Tippers With Wheels	1/pr

**6" Wheel Lock Extension**

STDS801	Wheel Lock Extension	1/ea
---------	----------------------	------

**Seat Belts**

STDS850	Seat Belts, Auto-Clasp Type Closure	1/ea
STDS851	Seat Belts, Velcro® Type Closure	1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea
STDS856	Seat Belts, Velcro® Type Closure, Bariatric	1/ea

**Oblique Projection Handrim**

STDS1022	Oblique Projection Handrim	1/ea
----------	----------------------------	------

**Wheelchair Cane/Crutch Holder**

STDS1034	Wheelchair Cane/Crutch Holder	1/ea
----------	-------------------------------	------

**Limb Support**

WASR	Limb Support, Right	1/ea
WASL	Limb Support, Left	1/ea

**Heel Loops**

STDS831	Heel Loops	1/pr
---------	------------	------

**One-Arm Drive**

STDS4055L	Left (For use with 18" and 20")	1/ea
STDS4055R	Right (For use with 18" and 20")	1/ea

**Telescoping IV Pole Attachment**

STDS820	Telescoping IV Pole Attachment	1/ea
---------	--------------------------------	------

**Universal Oxygen "E" Cylinder Holder With IV Pole**

STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment	1/ea
---------	--	------

**Seat Rail Extension**

STDS810	Seat Rail Extensions	1/pr
STDS1818N	Black Extension	1/ea
STDS1618N	Black Extension	1/ea
STDS2018N	Seat Rail Extension Kit	1/e

**Overhead Anti-Theft Device**

STDS823	Overhead Anti-Theft Device	1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks	1/ea
STDS834	Overhead Anti-Theft Device Single Pole	1/ea

**Chart Carry Pocket**

STDS835	16", 18" and 20" Wheelchairs w/ Nylon Upholstery	1/ea
---------	--	------