


Geo-Matt® Wedge Cushion

High resiliency, high density foam wedge has a 4" (front) to 2" (back) taper that helps keep users' hips from sliding forward. Facilitates restraint-free positioning without compromising skin integrity, thanks to its tissue-friendly Geo-Matt® surface geometry. Cover consists of anti-slip top surface and non-skid bottom. Built-in sleeve in bottom of cover allows use of molded Sling-Fill base (not included-see page 33); 18-month warranty. Weight limit 250 lbs.

Medicare Code:
General Use – E2601.



Item #	Description	Qty/Cs
GMW1616-96	Geo-Matt Wedge w/ anti-slip cover, 16"W x16"L	6
GMW1616-96-1	Geo-Matt Wedge w/ anti-slip cover, 16"W x16"L	1
GMW1816-96	Geo-Matt Wedge w/ anti-slip cover, 18"W x16"L	6
GMW1816-96-1	Geo-Matt Wedge w/ anti-slip cover, 18"W x16"L	1
GMW1818-96	Geo-Matt Wedge w/ anti-slip cover, 18"W x18"L	6
GMW1818-96-1	Geo-Matt Wedge w/ anti-slip cover, 18"W x18"L	1
GMW2018-96	Geo-Matt Wedge w/ anti-slip cover, 20"W x18"L	6
GMW2018-96-1	Geo-Matt Wedge w/ anti-slip cover, 20"W x 18"L	1

Fire barrier versions can be ordered to conform with 2000 NFPA 101 (Life Safety Code) and Cal TB# 133. To order, replace the "-96" suffix with "-29" in item number.

QUICK FACTS

Use for:

- Comfort
- Prevention of breakdown
- Assistance with positioning

