



Yes, you can.®

Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

*Saving and emailing the eForms is available with Acrobat Professional software

If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

Please note: The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit www.adobe.com.



Yes, you can.®

TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____

Company Name _____ Phone # _____

Address _____ City _____ State _____ Zip Code _____ Country _____

Contact Name _____ Comments _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

FRAME SELECTIONS

FRAME TYPE

- PROX4S X4 Folding Frame w/SWNG Front..... \$2,043.00 K0005
- PROX4F70 70 Degree Fixed Front Frame \$2,314.00
Not available with Heavy Duty, TRRO, or TRBKTS.
- PROX4F80 80 Degree Fixed Front Frame \$2,314.00
Not available with Heavy Duty, TRRO, or TRBKTS.

FRAME WIDTH

- W14 14" Wide No Charge A9900
Not available with Heavy Duty Option.
- W15 15" Wide No Charge
Not available with Heavy Duty Option.
- W16 16" Wide No Charge
Not available with Heavy Duty Option.
- W17 17" Wide No Charge
Not available with Heavy Duty Option.
- W18 18" Wide No Charge
- W19 19" Wide No Charge A9900
- W20 20" Wide \$300.00 E2201 A9900
- W22 22" Wide \$300.00 E2201 K0108
Must order Heavy-Duty Option.

SEAT DEPTH

- D14 14" Seat Depth No Charge A9900
Not available with Heavy-Duty Option.
- D15 15" Seat Depth No Charge
Not available with Heavy-Duty Option.
- D16 16" Seat Depth No Charge
Not available with Heavy-Duty Option.
- D17 17" Seat Depth No Charge
Not available with Heavy-Duty Option.
- D18 18" Seat Depth No Charge
- D19 19" Seat Depth No Charge A9900
- D20 20" Seat Depth \$300.00 E2203
- D21 21" Seat Depth \$300.00 E2203 K0108
Must order Heavy-Duty Option.
- D22 22" Seat Depth \$300.00 E2204 K0108
Must order Heavy-Duty Option.

HEAVY DUTY OPTION

Only available on model PROX4S and not available with TRRO or TRBKTS.

- HD Heavy Duty Option \$400.00 K0108

FRONT SEAT - TO - FLOOR

- FSTF155 Front STF Height 15 1/2" \$100.00 K0056 A9900
- FSTF16 Front STF Height 16" \$100.00 K0056 A9900
- FSTF165 Front STF Height 16 1/2" \$100.00 K0056 A9900
- FSTF17 Front STF Height 17" No Charge
- FSTF175 Front STF Height 17 1/2" No Charge
- FSTF18 Front STF Height 18" No Charge
- FSTF185 Front STF Height 18 1/2" No Charge
- FSTF19 Front STF Height 19" No Charge
- FSTF195 Front STF Height 19 1/2" No Charge
- FSTF20 Front STF Height 20" No Charge
- FSTF205 Front STF Height 20 1/2" No Charge
- FSTF21 Front STF Height 21" No Charge
- FSTF215 Front STF Height 21 1/2" \$100.00 K0056 A9900

| Width and Depth Matrix | | | | | | | | | |
|------------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| Depth Width | 14" | 15" | 16" | 17" | 18" | 19" | 20" | 21" | 22" |
| 14" | X | X | X | X | X | S | S | | |
| 15" | X | X | X | X | X | S | S | | |
| 16" | X | X | X | X | X | S | S | | |
| 17" | X | X | X | X | X | S | S | | |
| 18" | X | X | X | X | X* | X* | X* | S | S |
| 19" | S | S | X | X | X* | X* | X* | S | S |
| 20" | S | S | X | X | X* | X* | X* | * | * |
| 22" | | | S | S | S | S | S | S | * |

X - Sizes available
* - Sizes available with HD S - Sizes available through Specials

All prices subject to change without notice. All prices in U.S. dollars. All orders F.O.B. Pinellas Park, Florida USA

For the most current pricing information go to www.invacare.com.

Invacare® ProSpin™ X4™ Ultralightweight Active Performance Folder

| Front Seat - To - Floor Heights | | | | | |
|---------------------------------|--------|----------------|---------------------------|-----------------------|---------------------------|
| Fork | Spacer | Caster | | | |
| | | 3" AT900 | 4" 4SPLT 4ALUM2 | 5" 51DC 5ALUM2 | 6" 61DC 6ALUM2 |
| BRG1 (Short) | None | 15 1/2" 16" | 15 1/2" 16" 16 1/2" | 16 1/2" 17" | 17 1/2" |
| | 1" | 16 1/2" 17" | 16 1/2" 17" 17 1/2" | 17 1/2" 18" | 18 1/2" |
| | 2" | 17 1/2" 18" | 17 1/2" 18" 18 1/2" | 18 1/2" 19" | 19 1/2" |
| BRG2 (Long) | None | 17 1/2" 18" | 17 1/2" 18" 18 1/2" | 18" 18 1/2" 19" | 18 1/2" 19" 19 1/2" |
| | 1" | 18 1/2" 19" | 18 1/2" 19" 19 1/2" | 19" 19 1/2" 20" | 19 1/2" 20" 20 1/2" |
| | 2" | 19 1/2" 20" | 19 1/2" 20" 20 1/2" | 20" 20 1/2" 21" | 20 1/2" 21" 21 1/2" |
| SOFR (Suspension) | None | 16 1/2" 17" | 16 1/2" 17" 17 1/2" | 17 1/2" 18" | 18 1/2" |
| | 1" | 17 1/2" 18" | 17 1/2" 18" 18 1/2" | 18 1/2" 19" | 19 1/2" |
| | 2" | N/A | N/A | N/A | N/A |

Transport Tie Down*

TRRO includes four factory-installed transport brackets and wheelchair anchored pelvic belt. TRRO has been crash-tested in accordance with ANSI / RESNA WC Vol 1. Section 19 Frontal Impact Test requirements for wheelchairs with a 168 lbs. crash dummy, which corresponds to a person with a weight of 114 lbs. to 209 lbs. TRBKTS includes four factory - installed wheelchair transport brackets. TRBKTS has not been crash tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport.

As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

ARM SELECTION

ARM TYPE (PAIR)

- AT99 Swingaway Padded Arms No Charge
Not available with Heavy Duty Option or Angle Adjustable Backs (AAB1216 & AAB1720) or Transport Brackets (TRRO/TRBKTS)
- CAD Conventional Adj Hgt Desk Arms \$195.00 E0973
CAD and CAF are flip back and removable.
- CAF Conventional Adj Hgt Full Arms \$195.00 E0973
CAD and CAF are flip back and removable.
- TAD T - Arms, Adj. Height Desk Length \$225.00 E0973
Rigid sideguards included with the T-Arms. Bolt-on attachment.
- TAF T - Arms, Adj. Height Full Length \$225.00 E0973
Rigid sideguards included with the T-Arms. Bolt-on attachment.
- LESSARMS Less Arms No Charge

T - ARM HEIGHT

- L6-10 Low Range 6" - 10" (T - Arm Height) No Charge
- H8-14 High Range 8" - 14" (T - Arm Height) No Charge

| Arm Heights | | |
|---|------------------------|--|
| AT99 | 8" - 10" | Adjustable in 1" Increments |
| CAD CAF | 10" - 14" 10" - 14" | Adjustable in 1" Increments |
| TAD (Low Range) TAF (High Range) | 6" - 10" 8" - 14" | Adjustable in 1/2" Increments Adjustable in 1" Increments |
| All measurements from top of seat rail to top of armrest. On AT99 add 1" for upholstered arm pads. | | |

SUSPENSION SELECTIONS

TRANSPORT TIE DOWN

- Clear
- TRRO Transport Ready Option \$300.00
- TRBKTS Wheelchair Transport Brackets \$250.00

FRONT SUSPENSION

Not available with Heavy Duty Option.

- SOFR Frog Legs™ Front Suspension Forks \$315.00 E1015

FRONT SUSPENSION SETTING

- Clear
- SET160 Soft Setting (Max. User Weight 160 lbs.) No Charge
- SET300 Firm Setting (Max. User Weight 250 lbs.) Standard

All prices subject to change without notice. All prices in U.S. dollars. All orders F.O.B. Pinellas Park, Florida USA

For the most current pricing information go to www.invacare.com.

Invacare® ProSpin™ X4™ Ultralightweight Active Performance Folder

BACK SELECTION

BACK TYPE

- BACK1014 Adj. Back without Push Handles (10" - 14") No Charge
- BACK1216 Adj. Back without Push Handles (12" - 16") No Charge
- BACK1418 Adj. Back without Push Handles (14" - 18") No Charge
- BACK1620 Adj. Back without Push Handles (16" - 20") No Charge
- PH1216 Adj. Back with Push Handles (12" - 16") \$30.00
- PH1620 Adj. Back with Push Handles (16" - 20") \$30.00
- A1418 Adj. Back 10° Bend with Push Handles (14" - 18") \$30.00
- A1620 Adj. Back 10° Bend with Push Handles (16" - 20") \$30.00
- AAB1216 Adj. Angle Back with Push Handles (12" - 16") \$165.00
Not available with TRRO, TRBKTS or ADJTENS.
- AAB1720 Adj. Angle Back with Push Handles (17" - 20") \$165.00
Not available with TRRO, TRBKTS or ADJTENS.

BACK HEIGHT

- ABH10 10" High No Charge
- ABH11 11" High No Charge
- ABH12 12" High No Charge
- ABH13 13" High No Charge
- ABH14 14" High No Charge
- ABH15 15" High No Charge
- ABH16 16" High No Charge
- ABH17 17" High No Charge
- ABH18 18" High No Charge
- ABH19 19" High No Charge
- ABH20 20" High No Charge

BACK UPHOLSTERY TYPE

- SLIPBACK Slip On Back Upholstery No Charge
- ADJTENS Adjustable Tension Back Upholstery \$125.00 K0108
Not available with Heavy Duty Option.
- L8601 Omit Back Upholstery Less \$10.00

CASTER SELECTION

CASTER OPTIONS (PAIR)

- AT900 3" Roller Casters \$38.00 A9900
Not available with Heavy-Duty Option.
- 4ALUM2 4" Alum Soft Roll Casters \$92.00 A9900
- 4SPLT 4" Comp Urethane Light Up Casters No Charge
Not available with Heavy-Duty Option.
- 51DC 5" x 1" Urethane Casters No Charge
Not available with Heavy-Duty Option.
- 5ALUM2 5" Aluminum Soft Roll Casters \$92.00 A9900
- 61DC 6" x 1" Urethane Casters No Charge
Not available with Heavy-Duty Option.
- 6ALUM2 6" Aluminum Soft Roll Casters \$92.00 A9900

HARDWARE ANODIZING PACKAGE

Forks, spacers, caster housing and axle plate are included in the hardware anodizing package. **Not available with PROX4F70 or PROX4F80.**

- BLKANO Bold Black Anodized Accent Pkg. No Charge
- SVRANO Classic Silver Anodized Pkg. Standard

FORK TYPE (NON - SUSPENSION) (PAIR) Clear

- BRG1 Fork (3 - Hole Short Fork) No Charge
- BRG2 Fork (3 - Hole Long Fork) No Charge

FORK SPACER

- SPCR0 No Fork Spacer No Charge
- SPCR1 1" Fork Spacer No Charge
- SPCR2 2" Fork Spacer No Charge
Not available with SOFR.

REAR WHEELS

| Rear Seat-to-Floor Height Ranges | |
|---|----------------------|
| Rear Wheel | Seat-to-Floor Height |
| 20COM, 20SP | 13 1/2" - 18" |
| 22COM, 22SP | 14 1/2" - 19" |
| 24COM7, 24SP, 24HSP, 24SPIN | 15 1/2" - 20" |
| 25SPIN, 559TS, 26SP, 26HPS, 26SPIN | 16" - 20 1/2" |
| Rear Seat-to-Floor is factory set at 1" - 1 1/2" lower than Front Seat-to-Floor. | |

WHEEL OPTIONS (PAIR)

Spoke wheels come standard with dome, quick-release button axle.

- 20COM 20" Composite Wheels No Charge
Available with: PNEU, PNEUFF, or URETH tires.
- 22COM 22" Composite Wheels No Charge
Available with: PNEU, PNEUFF, or URETH tires.
- 24COM7 24" Composite-7 Spoke Wheels No Charge
Available with PNEU, PNEUFF, URETH, or KIKB tires.
- 20ITESPK20" Lightweight Spoke Wheel No Charge
**Not available with Heavy-Duty Option.
Available with PNEU or PNEUFF tires.**
- 22ITESPK22" Lightweight Spoke Wheel No Charge
**Not available with Heavy-Duty Option.
Available with PNEU or PNEUFF tires.**
- 24ITESPK24" Lightweight Spoke Wheel No Charge
**Not available with Heavy-Duty Option.
Available with PNEU, PNEUFF, PRIMO, or KIKB tires.**
- 26ITESPK26" Lightweight Spoke Wheel No Charge A9900
**Not available with Heavy-Duty Option.
Available with PNEU or PRIMO tires.**
- 559TS 559 Turbo Spoke \$117.00 A9900
**Not available with Heavy-Duty Option.
Available with PNEU or PRIMOTRD tires.**
- 24HPS 24" High Performance Spoke Wheels \$117.00 A9900
**Not available with Heavy-Duty Option.
Available with: PNEU or PNEUFF tires.**
- 26HPS 26" High Performance Spoke Wheels \$117.00 A9900
**Not available with Heavy-Duty Option.
Available with PNEU or PRIMO tires.**
- FN1624 24" Sun® Fusion 16 Spoke Wheel \$325.00 A9900
Not available with Heavy Duty Option.
- FN1625 25" Sun® Fusion 16 Spoke Wheel \$325.00 A9900
Not available with Heavy Duty Option.
- FN1626 26" Sun® Fusion 16 Spoke Wheel \$325.00 A9900
Not available with Heavy Duty Option.
- SPLX24 24" Spinergy® LX Spoke Wheels \$895.00 A9900
Choose Spoke color below.
- SPLX25 25" Spinergy® LX Spoke Wheels \$895.00 A9900
Choose Spoke color below.
- SPIN24 24" Spinergy® Spoke Wheels \$695.00 A9900
**Not available with Heavy Duty Option.
Available with PNEU tires. Choose Spoke color below.**
- SPIN25 25" Spinergy® Spoke Wheels \$695.00 A9900
**Not available with Heavy-Duty Option.
Available with PNEU or PRIMO tires. Choose Spoke color below.**
- SPIN26 26" Spinergy® Spoke Wheels \$695.00 A9900
**Not available with Heavy-Duty Option.
Available with PNEU or PRIMOTRD tires. Choose Spoke color below.**
- L8701 Omit Rear Wheels Less \$50.00
**Not available with Heavy-Duty Option.
When choosing omit rear wheels the frame will be constructed for 24" rear wheels.**

All prices subject to change without notice. All prices in U.S. dollars. All orders F.O.B. Pinellas Park, Florida USA

For the most current pricing information go to www.invacare.com.

Invacare® ProSpin™ X4™ Ultralightweight Active Performance Folder

REAR WHEELS continued

SPINERGY SPOKE COLOR

- Clear
- If color is not chosen, spokes will be black**
- SPCY Yellow No Charge
 - SPCBK Black Standard
 - SPCBW White No Charge
 - SPCR Red No Charge
 - SPCBL Blue No Charge

TIRE OPTIONS (PAIR)

- PNEU Pneumatic Tires \$60.00 ^{E2211} A9900
- PNEUFF Flat free Inserts \$65.00 ^{E2213}
- URETH Urethane (treaded) Tires No Charge
- PRIMO Primo Tires \$60.00 A9900
Not available with Heavy-Duty Option.
- PRIMOTRD Primo 559 Treaded Pneumatic \$60.00 A9900
Not available with Heavy-Duty Option.
- KIKB KIK® Black Tires \$60.00 A9900
- CRF123 Schwalbe® 1 in Blk Flat Resistant Non Marking Tires \$90.00
Only available with 24SP,26SP,24HPS,26HPS,SPIN24, SPIN25,SPIN26

HANDRIM OPTIONS (PAIR)

- ALUMPC Aluminum Anodized Handrim No Charge
- TIPC Titanium Handrims \$400.00 A9900
Not available with Heavy-Duty Option.
Not available with 20COM, 22COM, 24COM7, 20SP, 22SP, 559TS.
- A1512 Black Plastic Coated Handrims \$72.00 ^{K0059} A9900
- NFSDLT Natural-Fit™ LT Handrims Standard Grip \$325.00 A9900
Not available with 20COM, 22COM, 20SP, 22SP, 559TS.
- NFSPLT Natural-Fit LT Handrims Super Grip \$370.00 A9900
Not available with 20COM, 22COM, 20SP, 22SP, 559TS.
- NFITSD Natural-Fit Handrims Standard Grip \$295.00 A9900
Not available with 20COM, 22COM, 20SP, 22SP, 559TS.
- NFITSP Natural-Fit Handrims Super Grip \$340.00 A9900
Not available with 20COM, 22COM, 20SP, 22SP, 559TS.
- L8702 Omit Handrims Less 15.00
Not recommended with Heavy Duty Option.

MOUNTING HANDRIMS (SPOKE WHEEL ONLY) Clear

- Not available with Heavy-Duty Option.**
- IN Mounting Inside Handrim No Charge
 - OUT Mount Outside Handrim No Charge

AXLES (Pair)

- A1577 Quick Release Axles with Dome No Charge
Not available with HDPG.
- PAXLE Permanent Axles No Charge
- A1573 Quad Release Axles \$89.00
Not available with HDPG.

CAMBER ANGLE

- C00 0° Camber No Charge
- C03 3° Camber No Charge
- C06 6° Camber No Charge

FRAME FINISH

FRAME COLORS

- 24P Wet Black No Charge
- 29P White No Charge
- 60P Silver Metallic No Charge
- 115P Black Prism No Charge
- 105P Midnight Blue No Charge
- 30P Sunny Yellow No Charge
- 61P Electric Red No Charge
- 119P Bubblegum Pink No Charge
- 120P Grape Madness No Charge
- 116P Sunkissed Orange No Charge
- 117P Tangerine Red No Charge
- 62P Electric Blue No Charge
- 127P Grasshopper Green No Charge

Discontinued frame colors are available through our specials department.
C-070109030544

High Performance Wheel/Tire AVAILABLE Combinations

| Wheel Size (inches) | Spinergy LX & SPOX | SUN FUSION |
|---------------------|-------------------------|-------------------------------|
| 24" | PRIMO, CRF123 | PNEU, PRIMO, KIKB, CRF123 |
| 25" | PRIMO, PRIMOTRD, CRF123 | PRIMO, PRIMOTRD, CRF123, KIKB |
| 26" | PRIMO, CRF123 | PNEU, PRIMO, CRF123 |

POSITIONING STRAPS AND ACCESSORIES

SEAT POSITIONING STRAPS (EACH) Clear

- A1321BK Hook and Loop Seat Positioning Strap \$40.00 E0978
- A1515 Push Button Seat Positioning Strap \$50.00 E0978

BACK ACCESSORIES

- BP02 Back Pack \$50.00 A9270

SEAT ACCESSORIES Clear

- SP01 Seat Pouch \$35.00 A9270

ARM ACCESSORIES (PAIR)

- Only available with AT99.**
- A1571 Fabric Clothing Protectors \$46.00
 - 1671-SM Rigid Side Guard, Small \$122.00
 - 1671-LG Rigid Side Guard, Large \$122.00
 - 1671-SMH Quad Rigid Side Guard, Small \$122.00
 - 1671-LGH Quad Rigid Side Guard, Large \$122.00

PERFORMANCE DETAIL PACKAGES

Sport Package includes impact guards, rear frame protectors, and clear coat back cane brackets.

- SPORT ProSpin Sport Performance Detail Package \$65.00 A9270

All prices subject to change without notice. All prices in U.S. dollars. All orders F.O.B. Pinellas Park, Florida USA

For the most current pricing information go to www.invacare.com.

Invacare® ProSpin™ X4™ Ultralightweight Active Performance Folder

POSITIONING STRAPS AND ACCESSORIES continued

MISCELLANEOUS ACCESSORIES

- A1307 Spoke Protector (For 20", 22", 24", 26")(pr) \$90.00 K0065
 - 1209 Tire Pump \$35.00 A9270
 - A1580 Rear Anti - Tippers (pair) \$60.00 E0971
 - A1541 Caster Swivel Locks \$80.00
 - EXTPT Extended Amputee Axle Plate \$150.00 E0959 K0108
- Not available with Heavy-Duty Options.
Only available with 0° of camber. Not available with HUMV undermount wheel locks, cannot mount the wheellock far enough back to engage wheels. Not available with A1580.
- 6WLEXT 6" Wheel Lock Extensions \$39.00 E0961

WHEEL LOCK SELECTIONS

WHEEL LOCKS

- A1598 Hill Holder Clamp - On Wheel Locks \$170.00 E0974
Not available with URETH tires.
- WLPULL Pull To Lock Clamp - On Wheel Locks No Charge
- WLPUSH Push To Lock Wheel Locks No Charge
- HUMV Hiway Undermount Wheel Locks \$65.00 A9900
Not available with Heavy-Duty Option.
- L8704 Omit Wheel Locks Less \$10.00
Not available with 20COM, 22COM, 20SP, 22SP, 559TS

FRONT RIGGING OPTIONS

FRONT RIGGING (PAIR) Clear

- 60SSNTAPER 60° Swing In / Swing Out Footrests No Charge
- 70SSNTAPER 70° Non - Tapered Swing In / Out Ftrsts No Charge
- 70SSTAPER 70° Tapered Swing In / Swing Out Ftrsts No Charge
- ELRSS Elevating Legrests \$255.00 E0990

FOOTPLATE TYPE (PROX4S) (PAIR) Clear

- 1651 Flip - Up Composite Footplates No Charge
- AT5543 Adjustable Angle Flip - Up Footplates \$165.00 K0040
Doris/Plantar angle adjustable.
Not available with Elevating legrests.
- TUBULAR Individual Tubular Footplates \$59.00 A9900
Not available with Elevating legrests.

FOOTPLATE TYPE (PROX4F70 AND PROX4F80) Clear

- AT5545 One Pc. Adjustable Angle Flip Up Footplates \$165.00
Not available with Heavy-Duty Option.

FOOTREST / LEGREST ACCESSORIES Clear

- A1337 Calf Strap (each) \$24.00 K0038
- FRMPRBK ChairArmor Black Neoprene Frame Protectors \$65.00 A9270
- FRMPRR ChairArmor Red Neoprene Frame Protectors \$65.00 A9270
- FRMPRB ChairArmor Blue Neoprene Frame Protectors \$65.00 A9270
- FRMPRS ChairArmor Silver Neoprene Frame Protectors .. \$65.00 A9270
- ALPT Longer Pivot and Slide Tube (pair) \$89.00
- A9342 "H" Strap (each) \$51.00 K0039
- HLOOP Heel Loops without Ankle Straps (pair) \$51.00 E0951

| Knee - To - Heel Range (PROX4S) | | |
|--|---------|---------|
| Front Rigging Type | Min. | Max. |
| ELRSS | 15 1/4" | 21 1/4" |
| 60° | 14" | 18 3/4" |
| 70° | 13 3/4" | 18 3/4" |
| 70° Tapered | 14 1/2" | 19 1/2" |
| Range adjust in 1" increments Measurement taken from top of seat rail to back of footrests. | | |

| Knee - To - Heel Range (PROX4F70 and PROX4F80) | | |
|---|--------|--------|
| Front Rigging Type | Min. | Max. |
| Adjustable Angle Flip Up | 14.00" | 18.00" |
| 1 Pc. Tubular Flip Up | 14.00" | 18.00" |
| Range adjusts in 1/4" increments Measurement taken from top of seat rail to back of footrests. | | |

SEATING

SEATING CUSHION SELECTION Clear

- For additional information (size availability, additional seating options and pricing) See form 07-027 Rehab Seating and Positioning Products price list.
- EC Absolute Cushion \$104.00 E2601
Available in even sizes only. See form 07-027 for additional information.
 - SD Single Density Cushion \$170.00 E2602 E2601
 - ITSG InTouch™ Stabilite™ Gentle Contour Cushion \$320.00 E2607 E2608
 - ITSM InTouch Stabilite Max Contour Cushion \$320.00 E2605
 - ITFG InTouch Flovair™ Gentle Contour Cushion \$390.00 E2607 E2608
 - ITFM InTouch Flovair Max Contour Cushion \$390.00 E2607 E2608

SEATING BACK SELECTION

- For additional information (size availability, additional seating options and pricing) see form # 07-027 Rehab Seating and Positioning Products price list.
- PBR Personal Back 10™ Reg. Height \$500.00
 - PBT Personal Back 10 Tall Height \$500.00
 - PPR Personal Back 10 Plus Reg. Height \$600.00
 - PPT Personal Back 10 Plus Tall Height \$600.00
 - ITPR InTouch Propel™ Back \$490.00
 - ITPCSR InTouch PCS Regular \$645.00
 - ITPCST InTouch PCS Tall \$645.00

INTOUCH PCS™ ACCESSORIES

- PCSPS Pelvic Stabilizers, 3" (Pair).....\$230.00 E0956
- PCSTL Thoracic Lateral Supports, 4" (Pair).....\$230.00 E0956

For additional seating needs please refer to the following:
Invacare Matrix Seating and Positioning Products Price List Form
#11-036.
Please contact Invacare Customer Service to ship Invacare Matrix
Seating options with your chair order.

All prices subject to change without notice. All prices in U.S. dollars. All orders F.O.B. Pinellas Park, Florida USA

For the most current pricing information go to www.invacare.com.

Invacare® ProSpin™ X4™ Ultralightweight Active Performance Folder

PROSPIN X4 WEIGHT CHARTS

| FRAME WEIGHT SWING - AWAY | | | | | | | | |
|---------------------------|----------|--------|------|------|--------|--------|--------|-----------|
| W | 14" | 15" | 16" | 17" | 18" | 19" | 20" | HD |
| D | | | | | | | | |
| 14" | 7.25lbs. | 7.31 | 7.37 | 7.44 | 7.5-s | 7.56-s | 7.62-s | |
| 15" | 7.31 | 7.37 | 7.44 | 7.50 | 7.56-s | 7.62-s | 7.68-s | |
| 16" | 7.37 | 7.44 | 7.50 | 7.56 | 7.62-s | 7.68-s | 7.74-s | |
| 17" | 7.43 | 7.50 | 7.56 | 7.62 | 7.68-s | 7.74-s | 7.80-s | |
| 18" | 7.49 | 7.57 | 7.63 | 7.69 | 7.75 | 7.80 | 7.86 | |
| 19" | 7.55-s | 7.63-s | 7.69 | 7.75 | 7.81 | 7.87 | 7.93 | |
| 20" | 7.61-s | 7.69-s | 7.75 | 7.81 | 7.87 | 7.93 | 8.00 | |
| HD | | | | | | | | 12.42lbs. |

S = Sizes available through Specials Dept.

| FRAME WEIGHT FIXED FRONT | | | | | | | |
|--------------------------|----------|--------|------|------|--------|--------|--------|
| W | 14" | 15" | 16" | 17" | 18" | 19" | 20" |
| D | | | | | | | |
| 14" | 6.25lbs. | 6.31 | 6.37 | 6.44 | 6.5-s | 6.56-s | 6.62-s |
| 15" | 6.31 | 6.37 | 6.44 | 6.50 | 6.56-s | 6.62-s | 6.68-s |
| 16" | 6.37 | 6.44 | 6.50 | 6.56 | 6.62-s | 6.68-s | 6.74-s |
| 17" | 6.43 | 6.50 | 6.56 | 6.62 | 6.68-s | 6.74-s | 6.80-s |
| 18" | 6.49 | 6.57 | 6.63 | 6.69 | 6.75 | 6.80 | 6.86 |
| 19" | 6.55-s | 6.63-s | 6.69 | 6.75 | 6.81 | 6.87 | 6.93 |
| 20" | 6.61-s | 6.69-s | 6.75 | 6.81 | 6.87 | 6.93 | 7.00 |

S = Sizes available through Specials Dept.

| Component Weights | |
|---|-------|
| Arm Type (pair) | |
| AT99 Swing away padded tubular | +2.64 |
| Back Type (pair) | |
| Adjustable height back without push handles (10 - 14) | +1.60 |
| Adjustable height back without push handles (16 - 20) | +1.70 |
| Adjustable height back with push handles (12 - 14) | +1.67 |
| Adjustable height back with push handles (16 - 20) | +1.77 |
| Adjustable height 10° Bend back with push handles (14 - 18) | +1.00 |
| Adjustable height 10° Bend back with push handles (16 - 20) | +1.00 |
| Adjustable angle back with push handles (12 - 16) | +1.86 |
| Adjustable angle back with push handles (17 - 20) | +1.86 |
| Back Upholstery | |
| Slip Back upholstery | +.78 |
| Adjustable tension back upholstery | +97 |
| Front Casters (pair) | |
| 3" Roller Blade Casters | +486 |
| 4" Semi pneumatic | +686 |
| 4" Frog Legs® Aluminum Soft Roll Casters | +686 |
| 4" Composite Urethane Light Up Casters | +888 |
| 5" Semi Pneumatic Casters | +910 |
| 5" Frog Legs Aluminum Soft Roll Casters | +954 |
| 6" Composite Urethane 3 - Spoke Casters | +1.02 |
| 6" Frog Legs Aluminum Soft Roll Casters | +1.09 |
| Fork Type (pair) | |
| 3-Hole Short Fork | +1.63 |
| 3-Hole Short Fork with 1" Spacer | +1.79 |
| 3-Hole Short Fork with 2" Spacer | +1.94 |
| 3-Hole Long Fork Spacer | +1.72 |
| 3-Hole Long Fork with 1" Spacer | +1.87 |
| 3-Hole Long Fork with 2" Spacer | +2.07 |
| Frog Legs Front Suspension Forks | +2.03 |
| Frog Legs Front Suspension Forks with 1" Spacer | +2.18 |
| Rear Wheel Selections (pair) | |
| 20" Composite Wheels | +4.13 |
| 22" Composite Wheels | +4.70 |
| 24" Composite 7 Spoke Wheels | +5.28 |
| 24" High performance spoke | +4.01 |
| 25" High performance spoke | +4.17 |
| 26" High performance spoke | +4.34 |
| 24" Spinergy® | +3.18 |
| 25" Spinergy | +3.31 |
| 26" Spinergy | +3.49 |
| Tire Options (pair) | |
| Pneumatic Tires | +2.87 |
| Flat Free Inserts | +5.42 |
| Urethane (treaded) Tires | +3.00 |
| Primo Tires | +1.54 |
| Primo 559 Treaded Pneumatic Tires | +1.94 |
| KIK Black Tires | +4.28 |
| Handrims (pair) | |
| Aluminum | +1.80 |
| Titanium | +1.65 |
| Plastic Coated | +1.80 |
| Natural-Fit Standard Grip | +1.68 |
| Natural-Fit Supergrip | +1.70 |
| Wheel Locks (pair) | |
| Pull - to - Lock | +1.92 |
| Push - to - Lock | +1.92 |
| Sigeguards (pair) | |
| Rigid - small | +1.43 |
| Rigid - Large | +1.47 |
| Rigid Quad - small | +1.33 |
| Rigid Quad - large | +1.35 |
| Accessories | |
| Frame protectors | +1.10 |
| Calf strap | +1.10 |
| Seat pouch | +3.35 |
| Spoke guards | +1.60 |
| Anti - tip bars | +1.63 |
| Footrest (Swing - away with composite footplates-pair) | |
| 60° Swing - In Swing - Out Footrest (pair) | +2.34 |
| 70° Swing - In Swing - Out Footrests (pair) | +2.45 |
| Elevating Legrest (pair) | +7.45 |
| Footrest (Fixed Front version) | |
| Tubular | +2.34 |
| Adjustable Angle One Piece Footplate (width 14 - 15) | +2.41 |
| Adjustable Angle One Piece Footplate (width 16 - 17) | +2.52 |
| Adjustable Angle One Piece Footplate (width 18 - 20) | +2.65 |

Weights are approximate +/- .50 lbs.

All prices subject to change without notice. All prices in U.S. dollars. All orders F.O.B. Pinellas Park, Florida USA
 For the most current pricing information go to www.invacare.com.

Invacare® ProSpin™ X4™ Ultralightweight Active Performance Folder

Basic Chair Includes:

- Folding Curved Frame with Swingaway Front
- Seat Width (W14" - W19")
- Seat Depth (D14" - D19")
- Swingaway Padded Arms
- Front Seat To Floor Height Selection (17" - 21")
- Adjustable Height Back Selection w/o Push Handles
- Slip On Back Upholstery
- Powder Coated Frame Finish
- Caster Selection (SPLT4, 51DC, 61DC)
- Fork Type (BRG1 or BRG2)
- Fork Spacer Type (SPCR0, SPCR1, SPCR2)
- Wheel Options (20SP, 22SP, 24SP, 20COM, 22COM, 24COM7)
- Urethane (treaded) Tire (URETH)
- Aluminum Anodized Handrims (ALUMPC)
- Axle Type (PAXLE or A1577)
- Wheel Lock Selections (Push or Pull to Lock)
- Camber Angle Selection
- Handrim Options (AHDRM or CHDRM)
- Swingaway Front Rigging
- Flip Up Composite Footplates

Available Through Our Specials Department

| | |
|--|-----------------|
| Oxygen Holder | SM-010227014237 |
| 80° Swingaway Front Riggings | SM-200411121756 |
| Attendant Operated Wheel Locks | SM-071114092615 |
| One Arm Drive | SM-060314091819 |
| All specials are subject to final approval. Some special quotes are not compatible for all frame configurations. For your customization requests, please contact Customer Service for details and prices. | |

Specifications

| | |
|--|------------------------------------|
| Frame Style | Folding with Swingaway Front End |
| Seat Widths | 14" - 20" and 22" |
| Seat Depths | 14" - 22" |
| Back Heights | 10" up to 20" |
| Adjustable Angle Back | 0°, 5°, 10°, 15° |
| Camber Angles | 0°, 3°, 6° |
| Tubing O.D. | 1.0" |
| Weight of Chair | 15 lbs.* 30 lbs.** |
| Weight Limit | 250 lbs. |
| Heavy Duty | 350 lbs. |
| Warranty | Lifetime on frames and crossbraces |
| *16" x 16", BACK1216, SLIPBACK, 4ALUM2, BRG1, SPCR1, no rear wheels, no footrests, no anti - tippers, no wheel locks, no arms | |
| **16" x 18", AT99, BACK1216, SLIPBACK, 4ALUM2, BRG1, SPCR1, SPIN24, PRIMO, ALUMPC, A1577, 70SNTAPER, 1651, HLOOP, IMPACTG, 1671LG, WLPULL. | |

Weight Limitation 250 lbs.

Invacare recommends all Invacare wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

Comments

NOTE: All specifications and dimensions are approximate.
For the most current pricing information, go to www.invacare.com/homecare.

Invacare Corporation

www.invacare.com

USA

One Invacare Way
Elyria, OH USA
44035-4190
800-333-6900

©2010 Invacare Corporation

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. KIK is a registered trademark of KIK Technology, Inc. Spinergy is a registered trademark of Spinergy, Inc. Frog Legs is registered trademark of Frog Legs, Inc. Specifications and prices are subject to change without notification.

Page 8 of 8

