


<b>ORDER FORM</b>	
<b>Quickie® Match Point</b>	
	Failure to use an order form for a Quickie chair order may result in delay or wrong shipment



Please send all orders to:  
**Fax:** +44 (0)845 605 66 89  
**Email:** orders@sunmed.co.uk


Please ask for details on how to order online @  
[www.sunscript.co.uk](http://www.sunscript.co.uk)

**Sunrise Medical Ltd**  
 Thorns Road, Brierley Hill, West Midlands DY5 2LD  
[www.sunrisemedical.co.uk](http://www.sunrisemedical.co.uk)


**Standard features:**

- ◆ Base frame single integral anti-tip
- ◆ Choice of 16 frame colours
- ◆ Choice of 3 frame lengths
- ◆ Choice of seat widths ranging from 12" to 20"
- ◆ Fully adjustable frame
- ◆ Patented centre of mass adjustment
- ◆ Choice of seat depths ranging from 10" to 20"
- ◆ 24" spoked rear wheels
- ◆ Choice of folding and non-folding back
- ◆ 3" Micro castors
- ◆ Quick release axles
- ◆ Choice of 16° or 20° camber
- ◆ Maximum user weight limit 112kg (18st)


**How to measure:**




**Seat Depth**  
 From front of seat sling to back of seat sling



**Seat Width**  
 From outside seat rail tube to outside seat rail tube



**Frame Length**  
 From outside rear frame to outside front frame



**Back Height**  
 From seat sling to top of backrest sling

**Prices:**

All prices are in Pounds Sterling/Euros, exclude VAT. No responsibility held for technical changes or printing failures. Please consult Sunrise Medical with any queries you may have.

# ORDER FORM

## Quickie® Match Point



Failure to use an order form for a Quickie chair order may result in delay or wrong shipment

Please send all orders to:  
**Fax:** +44 (0)845 605 66 89  
**Email:** orders@sunmed.co.uk

Please ask for details on how to order online @  
[www.sunscript.co.uk](http://www.sunscript.co.uk)

Invoice to:  
 Name:  
 Street:  
 Town:  
 Postal Code:  
 Country:  
 Telephone: \_\_\_\_\_ Fax: \_\_\_\_\_

Ship to:  
 Name:  
 Street:  
 Town:  
 Postal Code:  
 Country:  
 Telephone: \_\_\_\_\_ Fax: \_\_\_\_\_

Order Date:  
 Purchase Order (Repeat on all pages):  
 Number Required:

Marked for:  
 Sunrise Order No:

Please fill boxes, do not tick       = Standard       = Option      N/C = No Cost Option      Std = Standard

		£			£	TOTAL
<b>MODEL (Standard)</b>						
156WM16	<input type="checkbox"/> Match Point (QPEISP3)	1,974				
<b>BASE FRAME (Standard)</b>						
156BF02	<input checked="" type="radio"/> Base frame single integral anti-tip	Std				
<b>FRAME COLOUR (Option)</b>						
156C140	<input type="radio"/> Aztec Gold	Std	156C114	<input type="radio"/> Sparkle Silver	Std	
156C2	<input type="radio"/> Black (not available as base with Zebra)	Std	156C145	<input type="radio"/> Sunrise Orange	Std	
156C141	<input type="radio"/> Black Cherry	Std	156C113	<input type="radio"/> Titanium	Std	
156C10	<input type="radio"/> Black Opal	Std	156C8	<input type="radio"/> Yellow	Std	
156C38	<input type="radio"/> Blue Opal	Std	156C116	<input type="radio"/> Matt Black	Std	
156C36	<input type="radio"/> Candy Blue	Std	156C151	<input type="radio"/> Matt Black Cherry	Std	
156C14	<input type="radio"/> Candy Purple	Std	156C150	<input type="radio"/> Matt Electric Blue	Std	
156C35	<input type="radio"/> Candy Red	Std	156C149	<input type="radio"/> Matt Evergreen	Std	
156C142	<input type="radio"/> Electric Blue	Std	156C152	<input type="radio"/> Matt Purple	Std	
156C143	<input type="radio"/> Evergreen	Std	<b>Custom Colours</b>			
156C126	<input type="radio"/> Glow	Std	156C112	<input type="checkbox"/> Desert Camouflage	209	
156C144	<input type="radio"/> Green Apple	Std	156C43	<input type="checkbox"/> Mossy Oak Camouflage	209	
156C147	<input type="radio"/> Hot Sparkle Pink	Std	156C42	<input type="checkbox"/> Stars & Stripes	209	
156C148	<input type="radio"/> Mauve Pink	Std	156C108	<input type="checkbox"/> Zebra (Choose base colour above)	209	
156C103	<input type="radio"/> Rootbeer	Std				
<b>FRAME LENGTH (Mark one box)</b>						
156FL50	<input type="radio"/> Frame length short 15" (38cm) (3)	Std	156FL52	<input type="radio"/> Frame length long 19" (43cm)	Std	
156FL51	<input type="radio"/> Frame length regular 17" (43cm) (2)	Std				
<b>SEAT WIDTH (See seat matrix) (Mark one box)</b>						
156S3	<input type="radio"/> Seat width 12" (30.5cm)	Std	156S9	<input type="radio"/> Seat width 18" (45.5cm)	Std	
156S4	<input type="radio"/> Seat width 13" (33cm)	Std	156S10	<input type="radio"/> Seat width 19" (48cm)	Std	
156S5	<input type="radio"/> Seat width 14" (35.5cm)	Std	156S47	<input type="radio"/> Seat width 20" (51cm)	Std	
156S6	<input type="radio"/> Seat width 15" (38cm)	Std	<b>Optional Seat Width:</b>			
156S7	<input type="radio"/> Seat width 16" (40.5cm)	Std	156OW1	<input type="radio"/> Additional 1/2" seat width (1)	N/C	
156S8	<input type="radio"/> Seat width 17" (43cm)	Std				
<b>SEAT DEPTH (See seat matrix) (Mark one box)</b>						
156SD3	<input type="radio"/> Seat depth 10" (25.5cm)	Std	156SD9	<input type="radio"/> Seat depth 16" (40.5cm)	Std	
156SD4	<input type="radio"/> Seat depth 11" (28cm)	Std	156SD10	<input type="radio"/> Seat depth 17" (43cm)	Std	
156SD5	<input type="radio"/> Seat depth 12" (30.5cm)	Std	156SD11	<input type="radio"/> Seat depth 18" (45.5cm)	Std	
156SD6	<input type="radio"/> Seat depth 13" (33cm)	Std	156SD12	<input type="radio"/> Seat depth 19" (48cm)	Std	
156SD7	<input type="radio"/> Seat depth 14" (35.5cm)	Std	156SD73	<input type="radio"/> Seat depth 20" (51cm)	Std	
156SD8	<input type="radio"/> Seat depth 15" (38cm)	Std				
<b>SEAT SLING</b>						
<b>Seat Sling:</b>			<b>Seat Pan:</b>			
156SL1	<input checked="" type="radio"/> Standard seat sling	Std	156PN5	<input type="checkbox"/> Seat pan half deep (2)	59	
<b>SEAT CUSHION (Option)</b>						
156CU1	<input type="checkbox"/> Cushion 2"	47	<input type="radio"/> Omit cushion		N/C	
<b>REAR SEAT HEIGHT (Mark one box)</b>						
156RH11	<input type="radio"/> Seat height 13" (33cm)	Std	156RH20	<input type="radio"/> Seat height 18" (45.5cm)	Std	
156RH13	<input type="radio"/> Seat height 14" (35.5cm)	Std	156RH21	<input type="radio"/> Seat height 19" (48cm)	Std	
156RH15	<input type="radio"/> Seat height 15" (38cm)	Std	156RH22	<input type="radio"/> Seat height 20" (51cm)	Std	
156RH17	<input type="radio"/> Seat height 16" (40.5cm)	Std	156RH23	<input type="radio"/> Seat height 21" (53.5cm)	Std	
156RH19	<input type="radio"/> Seat height 17" (43cm)	Std				
<b>FRONT SEAT HEIGHT (Mark one box)</b>						
156SH14	<input type="radio"/> Seat height 16" (40.5cm)	Std	156SH21	<input type="radio"/> Seat height 19" (48cm)	Std	
156SH16	<input type="radio"/> Seat height 17" (43cm)	Std	156SH24	<input type="radio"/> Seat height 20" (51cm)	Std	
156SH18	<input type="radio"/> Seat height 18" (45.5cm)	Std	156SH27	<input type="radio"/> Seat height 21" (53.5cm)	Std	

All prices are in Pounds Sterling/Euros. No responsibility held for technical changes or printing failures. Please consult Sunrise Medical with any queries you may have.

# Quickie® Match Point

Purchase Order (Repeat on all pages):

Please ask for details on how to order online @

[www.sunscript.co.uk](http://www.sunscript.co.uk)

		£			£	TOTAL
<b>BACKREST (Mark one type and choices underneath)</b>						
<b>Type:</b>				<b>Back Rest Upholstery:</b>		
156BT1	<input type="radio"/> Folding backrest	N/C	156BU1	<input type="radio"/> Standard upholstery	N/C	
156BT3	<input type="radio"/> Non-folding backrest	N/C	156BU13	<input type="checkbox"/> Adjustable sports upholstery	126	
<b>Back Height:</b>			156BU14	<input type="radio"/> Omit back upholstery	N/C	
156B128	<input type="radio"/> 9" (23cm)	N/C				
156B133	<input type="radio"/> 10" (25.5cm) - 13" (33cm)	N/C				
156B83	<input type="radio"/> 12" (30.5cm) - 16" (40.5cm)	N/C				
156B134	<input type="radio"/> 15" (38cm) - 18" (45.5cm)	N/C				
156B135	<input type="radio"/> 17" (43cm) - 20" (51cm)	N/C				
<b>FOOTRESTS (Mark one box)</b>						
156F234	<input type="radio"/> Platform	N/C	156F39	<input type="radio"/> High mount platform	N/C	
156IG20	<input type="radio"/> Impact guard neoprene					
<b>ADAPTIVE PLATFORM SPACER (Mark one box)</b>						
156AD1	<input type="checkbox"/> 2" Spacer (4)	79	156AD3	<input type="checkbox"/> 6" Spacer (4)	79	
156AD2	<input type="checkbox"/> 4" Spacer (4)	79				
<b>CASTORS (Mark one box)</b>						
156CW1	<input type="radio"/> 3" Micro	N/C				
<b>REAR WHEEL SIZE (Mark one box)</b>						
156RW53	<input type="radio"/> 24" X/S Spoke wheel (540)	N/C	156RW63	<input type="radio"/> Omit rear wheel and axle (540)	N/C	
156RW69	<input type="checkbox"/> 24" X/S QP Heat treated tennis wheel (559)	259	156RW64	<input type="radio"/> Omit rear wheel and axle (559)	N/C	
156RW57	<input type="checkbox"/> 24" Spinerig Sport Wheel (540) Black	492	156RW65	<input type="radio"/> Omit rear wheel and axle (590)	N/C	
156RW59	<input type="checkbox"/> 24" Spinerig Sport Wheel (559) Black	492				
156RW61	<input type="checkbox"/> 26" Spinerig Sport Wheel (590) Black	492				
<b>REAR WHEEL TYRES (Mark one box)</b>						
156RT31	<input type="radio"/> High pressure clincher	N/C				
156RT5	<input type="radio"/> Turbo tread	N/C				
<b>HANDRIMS (Mark one box)</b>						
156HR1	<input type="radio"/> Aluminium anodised	Std	156HR7	<input type="checkbox"/> Omit handrims	-20	
156HR4	<input type="checkbox"/> Plastic coated	67	<b>Handrim Tabs:</b>			
				156TB1	<input type="radio"/> Regular tab	N/C
<b>AXLES (Mark one box)</b>						
156AX4	<input type="radio"/> Quick-release axles	N/C				
156AX11	<input type="checkbox"/> Quad quick-release axles	77				
<b>CAMBER (Mark one box)</b>						
156CI10	<input type="radio"/> 16°	N/C	156CI11	<input type="radio"/> 20°	N/C	
<b>REAR ANTI-TIP (Mark one box)</b>						
156AT39	<input type="radio"/> Rear anti-tip integral	N/C				
156AT40	<input type="checkbox"/> Rear anti-tip quick release	57				
<b>ACCESSORIES (Optional)</b>						
<b>Back Pack:</b>			<b>Positioning Belt:</b>			
156A1	<input type="checkbox"/> Back pack adult	46	156A17	<input type="checkbox"/> Hook and loop closure	33	
156A2	<input type="checkbox"/> Back pack kids	46	156A18	<input type="checkbox"/> Nylon buckle closure	33	
<b>Seat Pouch:</b>			156A190	<input type="checkbox"/> Ratchet positioning belt (Large)	63	
156A8	<input type="checkbox"/> Seat pouch	46	156A191	<input type="checkbox"/> Ratchet positioning belt (Extra large)	63	
<b>Spoke Guard:</b>			<b>Miscellaneous:</b>			
156A24	<input type="checkbox"/> Spoke guard (Clear)	86	156A103	<input type="checkbox"/> Chest strap 2" (length 48")	17	
156A23	<input type="checkbox"/> Spoke guard (Black)	86	156A106	<input type="checkbox"/> Chest strap 4" (length 48")	18	
156A25	<input type="checkbox"/> Spoke guard sport (High profile)	101	156A90	<input type="checkbox"/> Pads for backrest crossbar	17	
156A26	<input type="checkbox"/> Spoke guard (Low profile)	101	156A98	<input type="checkbox"/> Tyre pump	67	
<b>Side Guard:</b>			156A180	<input type="checkbox"/> Toe clips	28	
156SG3	<input type="checkbox"/> Plastic (Regular)	127	156A181	<input type="checkbox"/> 5th Wheel support	63	
156SG4	<input type="checkbox"/> Plastic (Kids)	127				
156SG5	<input type="checkbox"/> Aluminium (Regular) (6)	175				
				<b>Base Price:</b>	£	
				<b>Options:</b>	£	
				<b>Total Price:</b>	£	

All prices are in Pounds Sterling/Euros, exclude VAT. No responsibility held for technical changes or printing failures. Please consult Sunrise Medical with any queries you may have.

# Quickie Match Point

Please ask for details on how to order online @

[www.sunscript.co.uk](http://www.sunscript.co.uk)

## FOOTNOTES:

- 1 1/2" seat width is built by selecting the Aluminium side guards
- 2 Not available with 10" to 11" seat depths
- 3 Not available with 17" to 20" seat depths
- 4 Not available with High mount platform footrest
- 5 Not available with Spinergy wheels
- 6 Aluminium side guards required for 1/2" seat width option

## Seat Matrix:

	DEPTH	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"
WIDTH												
12"		X	X	X	X	X	X	X				
13"			X	X	X	X	X	X	X			
14"				X	X	X	X	X	X			
15"					X	X	X	X	X	X		
16"					X	X	X	X	X	X		
17"						X	X	X	X	X	X	X
18"						X	X	X	X	X	X	X
19"						X	X	X	X	X	X	X
20"						X	X	X	X	X	X	X

All prices are in Pounds Sterling/Euros, exclude VAT. No responsibility held for technical changes or printing failures. Please consult Sunrise Medical with any queries you may have.