



YOUR
RESIDENTS
DESERVE A
❖ TRANQUIL ❖
NIGHT'S SLEEP



ALTEIRA 1385 HIGH-LOW BED

∞ QUIET ∞

STRENGTH

SAFETY

≡ ≡
< POSITIONING >
≡ ≡



QUIET

Top-of-the-line parallel electronics system provides exceptionally smooth and **QUIET** operation.

Perfectly synchronized high-low motors ensure the bed deck remains level even after frequent motion and long-term use.





STRENGTH

With a 500-lb. weight capacity, the Alterra 1385 has the **STRENGTH** to accommodate many of your residents.

Unique grid/rib deck prevents the structure from sagging.





SAFETY

Patent-pending folding footlock helps improve **SAFETY** by swinging away when not in use to reduce tripping.

The bed meets the FDA's guidelines to reduce the risk of entrapment.





< POSITIONING >

Multiple **POSITIONING** options are made possible by the bed's motors, which control the bed's height, knee hinge and reclining head deck.

The bed offers a full range of heights from a high height of 29.88" to a low height of 7.75".



ALTERRA 1385

HIGH-LOW BED



FOR A TRANQUIL NIGHT'S SLEEP

ORDERING INFORMATION





ITEM NO.	DESCRIPTION
FCE1385B	Alterra 1385 Bed *Shown with optional Soft-Touch Assist Rails and custom head and footboards.

[PRODUCT PAGE](#)

