



Yes, you can:

CS7™

Bed featured with ThinkSoft® Positioning Device, Traditional bed ends in Solar Oak and Solace® mattress



The innovative Carroll CS7 has already become the industry standard for long-term care beds. With an Auto Contour feature that simultaneously raises both the head and knees, residents stay comfortably in place while being better positioned for reading, watching TV or talking to family members. The simultaneous movement also helps prevent painful skin shearing that might otherwise occur during repositioning.

CS7™



Tool-less Assembly
Fast connect system allows for faster setup and easier alterations.



Hand Control Pendant
Carroll CS beds feature a rugged Hand Control, with bright, universal symbols, allowing residents and staff to easily operate the comfort controls of the bed.



Mobility-On-Demand
Floor lock system secures the bed to the floor at all height ranges but allows for mobility-on-demand.



Slat Deck
Exclusive arched slat deck design allows for improved infection control and reduces wear and tear on mattresses.

CS7 Specifications

- 7" – 30" Height Range
- 500 lbs. Weight Capacity
- Auto Contour
- One-Touch Floor Lock
- I.V. Rod Holders – standard
- Meets UL60601-1
- Meets CAN/CSA 22.2 No. 601.1-M-90
- Meets CAN/CSA 22.2 No. 601.1-2-38:03
- 8 Button Hand Pendant
- Ergonomically Designed, Rugged, Water Resistant Pendant
- Universal Deck Length (76" or 80")
- Mattress Keepers
- Optional Attendant Control Panel
- Optional Battery Backup
- Ratchet Heel Lift

CS7 Accessories

- Trapeze
- Hi-Low Lockout
- CPR Release
- 4" Bed Extender

Invacare Continuing Care

www.invacare-ccg.com

info_iccg@invacare.com

USA
One Invacare Way
Elyria, Ohio 44035
USA
800-668-2337

Canada
994 Hargrieve Road
London, Ontario
N6E 1P5 Canada
800-668-2337

Trademarks are identified by the symbol™ and®.
All trademarks are owned or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

© 2009 Invacare Corporation. Form 08.30.1/09



Yes, you can.®