


* HCPCS CODING:
E2601 / E2602

General Use



QUADRAGEL PACK HAS FOUR COMPARTMENTS DESIGNED TO PREVENT GEL MIGRATION. IT HELPS SKIN AND RETAINS BODY TEMPERATURE TO PROVIDE MAXIMUM COMFORT.



**Non-SKID
FLAT BOTTOM**

**SINGLE-LAYER FOAM
CONSTRUCTION** CONTOURS TO USERS BODY, REDUCING HIGH PRESSURE AREAS.

AVAILABLE IN BOTH 2" AND 3" THICKNESSES

-COVERS-

STRETCH-AIR FABRIC — designed for air transmission and heat reduction. This two-ply, breathable laminate combines a high-stretch polyester top layer with a waterproof, breathable vapor barrier.

COMFORT-TEK — fluid proof exterior has supple leathery properties that provide an easily cleaned surface and protection for the inner cushion. Comfort-Tek has multidimensional stretch and helps alleviate pressure on bony prominences.

CUSHION STYLE: 2" Thick w/ Gel	FABRIC TYPE: Stretch-Air	
MODEL NUMBER	W X L	WEIGHT LIMIT
46GS-1616-B	16 x 16	250 lbs.
46GS-1618-B	16 x 18	250 lbs.
46GS-1620-B	16 x 20	250 lbs.
46GS-1816-B	18 x 16	300 lbs.
46GS-1818-B	18 x 18	300 lbs.
46GS-1820-B	18 x 20	300 lbs.
46GS-2016-B	20 x 16	350 lbs.
46GS-2018-B	20 x 18	350 lbs.
46GS-2216-B	22 X 16	350 lbs.
46GS-2218-B	22 x 18	400 lbs.
46GS-2220-B	22 x 20	400 lbs.
46GS-2416-B	24 x 16	400 lbs.
46GS-2418-B	24 x 18	400 lbs.
46GS-2420-B	24 x 20	400 lbs.

ALL CUSHIONS AVAILABLE WITH RIGID INSERT AND EXTRA INCONTINENT PROOF LINER (PAGE 11)

CUSHION STYLE: 2" Thick w/ Gel	FABRIC TYPE: Comfort-Tek	
MODEL NUMBER	W X L	WEIGHT LIMIT
46G-1616-B	16 x 16	250 lbs.
46G-1618-B	16 x 18	250 lbs.
46G-1620-B	16 x 20	250 lbs.
46G-1816-B	18 x 16	300 lbs.
46G-1818-B	18 x 18	300 lbs.
46G-1820-B	18 x 20	300 lbs.
46G-2016-B	20 x 16	350 lbs.
46G-2018-B	20 x 18	350 lbs.
46G-2216-B	22 X 16	350 lbs.
46G-2218-B	22 x 18	400 lbs.
46G-2220-B	22 x 20	400 lbs.
46G-2416-B	24 x 16	400 lbs.
46G-2418-B	24 x 18	400 lbs.
46G-2420-B	24 x 20	400 lbs.

ALL CUSHIONS AVAILABLE WITH RIGID INSERT AND EXTRA INCONTINENT PROOF LINER (PAGE 11)

CUSHION STYLE: 3" Thick w/ Gel	FABRIC TYPE: Stretch-Air	
MODEL NUMBER	W X L	WEIGHT LIMIT
463GS-1616-B	16 x 16	250 lbs.
463GS-1618-B	16 x 18	250 lbs.
463GS-1620-B	16 x 20	250 lbs.
463GS-1816-B	18 x 16	300 lbs.
463GS-1818-B	18 x 18	300 lbs.
463GS-1820-B	18 x 20	300 lbs.
463GS-2016-B	20 x 16	350 lbs.
463GS-2018-B	20 x 18	350 lbs.
463GS-2216-B	22 X 16	350 lbs.
463GS-2218-B	22 x 18	400 lbs.
463GS-2220-B	22 x 20	400 lbs.
463GS-2416-B	24 x 16	400 lbs.
463GS-2418-B	24 x 18	400 lbs.
463GS-2420-B	24 x 20	400 lbs.

ALL CUSHIONS AVAILABLE WITH RIGID INSERT AND EXTRA INCONTINENT PROOF LINER (PAGE 11)

CUSHION STYLE: 3" Thick w/ Gel	FABRIC TYPE: Comfort-Tek	
MODEL NUMBER	W X L	WEIGHT LIMIT
463G-1616-B	16 x 16	250 lbs.
463G-1618-B	16 x 18	250 lbs.
463G-1620-B	16 x 20	250 lbs.
463G-1816-B	18 x 16	300 lbs.
463G-1818-B	18 x 18	300 lbs.
463G-1820-B	18 x 20	300 lbs.
463G-2016-B	20 x 16	350 lbs.
463G-2018-B	20 x 18	350 lbs.
463G-2216-B	22 X 16	350 lbs.
463G-2218-B	22 x 18	400 lbs.
463G-2220-B	22 x 20	400 lbs.
463G-2416-B	24 x 16	400 lbs.
463G-2418-B	24 x 18	400 lbs.
463G-2420-B	24 x 20	400 lbs.

ALL CUSHIONS AVAILABLE WITH RIGID INSERT AND EXTRA INCONTINENT PROOF LINER (PAGE 11)

