



ErgoFlight ***Ultra Lightweight Ergonomic Wheelchair***

- 14.5 lb Aluminum Frame. This Manual Wheelchair Competes Against Titanium
- Patented S-Style Ergonomic Seat Frame
- 6061 T-6 Aircraft-grade Aluminum Only 19.8 lbs!
- Companion Brakes Located on Handles for Added Safety and Comfort.
- Built in Silver Aegis Anti-bacterial Seat and Back Cushion
- Folding Backrest/Folding Seat for Transporting in Vehicle or Travel
- Fixed Armrests w/Ergonomic Arm Pads
- Removable Swing-in & Out Footrests
- 20” Rear Polyurethane, High Tread, Flat Free Wheels
- 6” x 1” Polyurethane Front Casters
- 18” x 16”, 16” x 16” Seat Width available
- Available Quick Release Rear Wheels
- Frame Color: Pearl Silver
- Weight Capacity: 220 lbs
- Limited Lifetime Warranty on Frame

Taking Aluminum to Titanium Heights

Titanium Weight at Aluminum price!

***U.K. ISO 7176/19 Dynamic
Crash Test Passed***

| Model No. | seat width (in.) | seat depth (in.) | armrest height (in.) | seat height (in.) | back height (in.) | overall height (in.) | overall open width (in.) | folded width (in.) | overall length (in.) | weight w/o riggings (lb.) | Wt. Cap. (lb.) |
|-------------|------------------|------------------|----------------------|-------------------|-------------------|----------------------|--------------------------|--------------------|----------------------|---------------------------|----------------|
| S-2512F18SS | 18 | 17 | 8 | 18.5 | 17.5 | 36 | 27 | 15 | 38 | 19 | 220 |
| S-2512F16SS | 16 | 17 | 8 | 18.5 | 17.5 | 36 | 25 | 15 | 38 | 19 | 220 |
| S-2512Q18SS | 18 | 17 | 8 | 18.5 | 17.5 | 36 | 27 | 15 | 38 | 19 | 220 |
| S-2512Q16SS | 16 | 17 | 8 | 18.5 | 17.5 | 36 | 25 | 15 | 38 | 19 | 220 |

****Available Options ~ Quick Release wheels, Anti-Tipper, Seat belt, Extra Upholstery, Memory Foam cushions, UV Resistant Canopy, Rigidifying Push Bar, Foldable/Adjustable Headrest***

karma KARMAN HEALTHCARE Inc.

19255 San Jose Avenue, City of Industry, CA 91748-1418

Order. 1-800-805-2762 / Tel. 626-581-2235 / Fax. 626-581-2335