

# GOCHAIR®

**Pride®**  
Mobility Products Corp.  
— LIVE YOUR BEST™ —

The first travel chair that combines super portability, maneuverability and style at a low, low cost!



## Specifications

<b>Drive Wheels</b>	8" (solid)
<b>Front Casters</b>	5" (solid)
<b>Rear Anti-tips</b>	3" (solid)
<b>Maximum Speed</b>	3.5 mph
<b>Ground Clearance</b>	1.75"
<b>Turning Radius</b>	26.75" with foot platform
<b>Length</b>	29" without foot platform
<b>Width</b>	19"
<b>Seat Type</b>	Foldable, compact, lightweight
<b>Batteries</b>	(2) 12V, 12 amp
<b>Range</b>	Up to 8 miles
<b>Charger</b>	Off-board, 2 amp
<b>Electronics</b>	40 amp, PG GC 2 (non-programmable)
<b>Weight Capacity</b>	250 lbs.
<b>Base Weight</b>	57 lbs. (front & rear)
<b>Standard Seat Weight</b>	25 lbs.
<b>Battery Pack Weight</b>	20 lbs.
<b>Suspension</b>	Limited
<b>Heaviest Piece (rear)</b>	34.5 lbs.

## Portability

- Easy-to-remove battery pack
- Conveniently disassembles into four manageable pieces
- Fits in most vehicle trunks

## Performance

- Top speed of 3.5 mph and enough power to get you where you want to go
- Twin motors for traction and precision handling

## Maneuverability

- Compact design for a tight turning radius

## Colors



Red



Blue

(US) 800-800-8586 • (Canada) 888-570-1113 • [www.pridemobility.com](http://www.pridemobility.com)

Range and speed vary with user weight, terrain type, battery charge, battery condition and tire pressure. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. This product is not coded by DME PDAC or does not meet criteria: K0899.

**PRIDE®**  
PERFORMANCE MOBILITY  
— LIVE YOUR BEST™ —