

# GO GO<sup>®</sup>

**LX** WITH CTS SUSPENSION



Comfort-Trac Suspension




## Features, Performance and Travel Ready

The Go-Go<sup>®</sup> LX with CTS (Comfort-Trac Suspension) brings a new level of performance features and value to travel mobility. Featuring CTS independent front and rear suspension, sleek style, feather-touch disassembly and standard lighting, the Go-Go LX with CTS Suspension is the most advanced travel product available.



**Pride<sup>®</sup>**  
Mobility Products Corp.  
— LIVE YOUR BEST<sup>™</sup> —



## SPECIFICATIONS

|                               |               |
|-------------------------------|---------------|
| Model number                  | S50LX         |
| Max. weight capacity          | 300 lbs.      |
| Maximum speed <sup>1</sup>    | Up to 4.4 mph |
| Ground clearance <sup>2</sup> | 1.5" at motor |
| Turning radius <sup>2</sup>   | 33.5"         |
| Overall length <sup>2</sup>   | 39.75"        |
| Overall width <sup>2</sup>    | 20.5"         |

## TIRES:

|           |                 |
|-----------|-----------------|
| Front     | 2" x 8" solid   |
| Rear      | 2.5" x 8" solid |
| Anti-tips | Rear            |

|                  |   |
|------------------|---|
| Standard seating | Foldable, compact, vinyl-covered molded plastic |
|------------------|---|

|            |                                    |
|------------|------------------------------------|
| Drivetrain | 24-volt DC motor, sealed transaxle |
|------------|------------------------------------|

|                     |                                    |
|---------------------|------------------------------------|
| Dual braking system | Regenerative and electromechanical |
|---------------------|------------------------------------|

|                 |           |
|-----------------|-----------|
| Battery charger | Off-board |
|-----------------|-----------|

|                                 |                 |
|---------------------------------|-----------------|
| Per-charge range <sup>1,4</sup> | Up to 5.9 miles |
|---------------------------------|-----------------|


|                                   |                               |
|-----------------------------------|-------------------------------|
| Battery requirements <sup>3</sup> | (2) 12V deep-cycle (standard) |
|-----------------------------------|-------------------------------|

|                      |       |
|----------------------|-------|
| Compatible batteries | 12 AH |
|----------------------|-------|

## WEIGHT:

|                                     |                         |
|-------------------------------------|-------------------------|
| Batteries (2 required) <sup>6</sup> | 12 AH (19.6 lbs. each)  |
| Weight w/o battery pack             | 85.4 lbs.               |
| Heaviest piece                      | 34 lbs. (front section) |

|                               |  |
|-------------------------------|--|
| Interchangeable shroud panels |  |
|-------------------------------|--|



Red



Blue

## FEATURES

- CTS (Comfort-Trac front and rear independent suspension)
- Feather-touch disassembly permits simple frame separation with only one hand
- Pride's exclusive black, non-scuffing tires
- Auto-connecting front-to-rear cable
- Includes 2 sets of easily changeable colored panels in Red and Blue
- Front frame-mounted seat post offers maximum stability
- All-new frame design easily disassembles into 5 super lightweight pieces for convenient transport and storage
- Modular design for easy serviceability
- Microprocessor-based controller offers optimal power management and added safety features
- Standard front tiller-mounted basket included
- Convenient off-board dual voltage charger can charge battery pack on-board or off-board
- Charger port in tiller

## Pride® Lifts & Ramps

These Pride® Lifts are compatible with the Go-Go® LX with CTS Suspension.

- Outlander Exterior Lift
- Power Tote Mini



Power Tote Mini



Outlander

<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition.

This specification can be subject to a variance of (+ or -) 10% - 5%.

<sup>2</sup>Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.

<sup>3</sup>AGM or gel cell type required.

<sup>4</sup>Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4, and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>5</sup>Includes standard seat.

<sup>6</sup>Battery weight may vary based on manufacturer.

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.