



Jazzy[®]

ELITE ES

The Jazzy[®] Elite ES delivers a reliable blend of power, performance and style. In-line, front-wheel drive technology gives the Elite ES excellent stability and maneuverability for solid performance indoors and out.



Pride[®]
Mobility Products Corp.
LIVE YOUR BEST[®]

SPECIFICATIONS

Weight capacity	300 lbs.
Maximum speed	Up to 3.7 mph
Ground clearance	1.5"
Turning radius	24.75"
Overall length	35" without front riggings
Overall width	23"
Drive wheels	9" solid
Front wheels	3" solid anti-tips
Rear wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	40 amp, PG GC 3 controller
Specialty controls	N/A
Battery charger	Off-board, 3 amp
Per-charge range	Up to 8.95 miles
Battery requirements	(2) 12 volt, deep cycle
Battery size	U-1
Battery weight	18 lbs. each
Base weight	68.5 lbs.
Standard seat weight	43 lbs. (comfort high-back)
Max. seat dimensions	20" x 20"
Warranty	5-year limited warranty on frame; 13-month limited warranty on electronics; 13-month limited warranty on drive motors

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

MODELS

Jazzy Elite ES-2S-SS	(Code: K0822)
Jazzy Elite ES-2S-C	(Code: K0823)

FEATURES

- In-line, front-wheel drive technology
- 40 amp, PG GC 3 controller
- Simple main frame design for easy serviceability
- Comfort high-back seat with headrest
- New Jazzy® Armor™ controller guard to protect from daily wear and tear
- Basic Jazzy® Armor™ shroud decals

OPTIONS

- High-back seat with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

RECOMMENDED LIFTS

Pride® Lifts & Ramps


- Backpacker® Plus
- Outlander
- Backpacker® MV



COLOR CHOICES



Jazzy Red



Jazzy Blue