



## KM-8520-22W

*Aluminum high strength, lightweight  
flip-back-arm & backrest wheelchair*

- 6061 T-6 Aircraft-grade aluminum frame in black powder coat finish
- Padded flip-back full length armrests with side panels
- Double crossbraces
- Swing-away detachable footrests with Velcro mount leg strap
- 24" x 1 3/4" spoke wheels with polyurethane tires
- 8" x 2" heavy duty flat free polyurethane casters
- Adjustable toggle wheel locks
- Adjustable length anti-slip footplates
- Flip-back backrest.
- High-grade flame retardant, breathable nylon upholstery

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
KM-8520-22W	22	18	9	20	16	36	30	14	41	35	350

Common Options: \*Not available on all models \*\* Price & specifications are subject to change without notice.

- Elevating Legrests • Anti-Tippers • Brake Extension • Oxygen Tank Holder • Seat Belt • Leg Straps •



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335