



Quickie LXI/LX December 2014

EILX/EILX

Mark For:

Submitting for:

 Quote

 Order

Date: _____

PO#: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: Fax Email

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.
The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME STYLE

<input type="checkbox"/>	065WM31	K0005	Quickie LXI (Adjustable Axle Plate) 265 lbs Weight Capacity	\$1,895.00	<input type="checkbox"/>	065WM30	K0004	Quickie LX (Non Adjustable Axle Plate) 250 lbs Weight Capacity	\$1,805.00
--------------------------	---------	-------	--	------------	--------------------------	---------	-------	---	------------

2. FRAME SET UP

Transit Option

Includes stabilizer bar. Only on LXI. NA w/ padded swing away armrest

<input type="checkbox"/>	065TO1	Transit System	\$315.00
--------------------------	--------	----------------	----------

Frame Width

<input type="checkbox"/>	065W3	12"	NC
<input type="checkbox"/>	065W4	13"	NC
<input type="checkbox"/>	065W5	14"	NC
<input type="checkbox"/>	065W6	15"	NC
<input type="checkbox"/>	065W7	16"	NC
<input type="checkbox"/>	065W8	17"	NC
<input type="checkbox"/>	065W9	18"	NC
<input type="checkbox"/>	065W41	19"	\$125.00
<input type="checkbox"/>	065W11	E2201 20"	\$500.00

Seat Depth

12"-14" NA w/ 8" casters or 24" wheels

<input type="checkbox"/>	065SD5	12"	NC
<input type="checkbox"/>	065SD6	13"	NC
<input type="checkbox"/>	065SD7	14"	NC
<input type="checkbox"/>	065SD8	15"	NC
<input type="checkbox"/>	065SD9	16"	NC
<input type="checkbox"/>	065SD10	17"	NC
<input type="checkbox"/>	065SD11	18"	NC
<input type="checkbox"/>	065SD74	19"	\$125.00
<input type="checkbox"/>	065SD13	E2203 20"	\$550.00

3. FRONT SEAT HEIGHT SET UP

Front Seat Height

Rear STFH ranges from 0" to 1 1/2" lower than front STFH

<input type="checkbox"/>	065SH12	15"	NC
<input type="checkbox"/>	065SH13	15.5"	NC
<input type="checkbox"/>	065SH14	16"	NC
<input type="checkbox"/>	065SH15	16.5"	NC
<input type="checkbox"/>	065SH16	17"	NC
<input type="checkbox"/>	065SH17	17.5"	NC
<input type="checkbox"/>	065SH18	18"	NC
<input type="checkbox"/>	065SH20	18.5"	NC
<input type="checkbox"/>	065SH21	19"	NC
<input type="checkbox"/>	065SH23	19.5"	NC
<input type="checkbox"/>	065SH24	20"	NC
<input type="checkbox"/>	065SH26	20.5" (Only available on swing away)	NC

Wheel	FRONT SEAT HEIGHT ORDERING GUIDE				
	Caster				
	4"	5"	6"	8" POLY	8"x2" PNU
20"	15.0	15.0	16.0		18.0
	15.5	16.0	16.5		
	16.0	16.5	17.0		
	16.5	17.0	17.5		
	17.0	17.5	18.0		
	17.5	18.0			
22"	18.0				
	16.0	16.0	16.0	18.5	19.0
	16.5	16.5	16.5	19.0	
	17.0	17.0	17.0		
	17.5	17.5	17.5		
	18.0	18.0	18.0		
24"	18.5	18.5	18.5		
	17.00	17.0	17.0	18.5	18.0
	17.50	17.5	17.5	19.0	19.0
	18.00	18.0	18.0	19.5	19.5
	18.50	18.5	18.5	20.0	
	19.00	19.0	19.0		
25"	19.50	19.5	19.5		
	20.0	20.0	20.0		
	17.00	17.0	17.0	18.5	19.0
	17.50	17.5	17.5	19.0	19.5
	18.00	18.0	18.0	19.5	
	18.50	18.5	18.5	20.0	
26"	19.00	19.0	19.0		
	19.50	19.5	19.5		
	20.0	20.0	20.0		
	18.00	18.0	18.0	18.5	18.0
	18.50	18.5	18.5	19.0	19.5
	19.00	19.0	19.0	19.5	20.5
	19.50	19.5	19.5	20.0	
	20.0	20.0	20.5		

3. FRONT SEAT HEIGHT SET UP (CONT)

Caster Wheels			
<input type="checkbox"/>	065CW3		4" Low Profile Polyurethane NC
<input type="checkbox"/>	065CW31		4" Micro Lighted \$70.00
<input type="checkbox"/>	065CW72		4" x 1" Aluminum Polyurethane \$100.00
<input type="checkbox"/>	065CW44		4" x 1.25" Semi Pneumatic \$55.00
<input type="checkbox"/>	065CW64		4" x 1.5" Aluminum Soft Roll \$175.00
<input type="checkbox"/>	065CW5		5" Low Profile Polyurethane NC
<input type="checkbox"/>	065CW51		5" Micro Lighted \$85.00
<input type="checkbox"/>	065CW65		5" x 1" Aluminum Polyurethane \$110.00
<input type="checkbox"/>	065CW7		5" x 1.5" Semi-Pneumatic \$85.00
<input type="checkbox"/>	065CW54		5" Soft Roll (Grey) \$85.00 <i>NA w/ Frog Legs</i>
<input type="checkbox"/>	065CW66		5" x 1.5" Aluminum Soft Roll \$185.00
<input type="checkbox"/>	065CW34		6" Polyurethane NC
<input type="checkbox"/>	065CW33	E2214 E2215	6" Pneumatic \$75.00
<input type="checkbox"/>	065CW35		6" x 1.5" Semi Pneumatic \$85.00
<input type="checkbox"/>	065CW68		6" x 1.5" Aluminum Soft Roll \$195.00
<input type="checkbox"/>	065CW16		7" Polyurethane NC
<input type="checkbox"/>	065CW32		8" Polyurethane NC
<input type="checkbox"/>	065CW21	E2214 E2215	8" Pneumatic \$75.00
<input type="checkbox"/>	065CW39		8" x 1.5" Semi Pneumatic \$85.00
<input type="checkbox"/>	065CW36	E2214 E2215	8" x 2" Pneumatic \$145.00
<input type="checkbox"/>	065CW37	E2214 K0108	8" x 2" Pneumatic w/Airless Insert \$225.00
<input type="checkbox"/>	065CW38		Omit Caster Wheel NC

Caster Forks			
<input type="checkbox"/>	065CF62		3" NC
<input type="checkbox"/>	065CF9		4 3/4" NC
<input type="checkbox"/>	065CF12		5 1/4" NC
<input type="checkbox"/>	065CF17		6" NC
<input type="checkbox"/>	065CF23		7" NC
<input type="checkbox"/>	065CF29		8" x 2" NC
<input type="checkbox"/>	065CF30		Omit Caster Fork and Wheel NC

Fork Stems			
<i>Caster pin locks N/A w/ 8"x2" casters</i>			
<input type="checkbox"/>	065CS1		Standard Stem NC
<input type="checkbox"/>	065CS5		3/4" Stem NC
<input type="checkbox"/>	065CS9		1 1/2" Stem NC
<input type="checkbox"/>	065CS7	K0073	Standard Stem w/Pin Locks \$195.00
<input type="checkbox"/>	065CS11	K0073	3/4" Stem w/Pin Locks \$195.00
<input type="checkbox"/>	065CS4	K0073	1 1/2" Stem w/Pin Locks \$195.00

4. HANGER & FOOTPLATE OPTIONS

Hangers			
<i>Front mount hangers only available w/ aluminum angle adj footplates</i>			
<input type="checkbox"/>	065H85		60° Swing In/Out Hanger NC
<input type="checkbox"/>	065H86		70° Swing In/Out Hanger NC
<input type="checkbox"/>	065H91		90° Swing In/Out Hanger NC <i>NA w/ 8" casters. Pull to lock wheel locks suggested</i>
<input type="checkbox"/>	065H24	K0053	Articulating Hanger (15"-20"w) \$350.00 <i>NA w/ Push To Lock Wheel Locks</i>
<input type="checkbox"/>	065H14	E0990	Elevating Hanger (14"-20"w) \$295.00 <i>Pull to lock wheel locks suggested. Only available w/alum and angle adj footplates</i>
<input type="checkbox"/>	065H32		Omit Hanger -\$15.00

Extension Tube			
<input type="checkbox"/>	065E1		Short (See Matrix for Lengths) NA w/ 90° hangers and ELR NC
<input type="checkbox"/>	065E2		Medium (See Matrix for Length) NC
<input type="checkbox"/>	065E5		Long (See Matrix for Length) NA w/ 90o hangers and ELR NC
<input type="checkbox"/>	065E14		Regular ELR NC
<input type="checkbox"/>	065E15		Long ELR NC
<input type="checkbox"/>	065E16	K0037	High Mount Only available with angle adj footplates or 60° and 70° hangers \$55.00

Footplate & Hanger Options (Ext Tube & Direct Mount)			
<input type="checkbox"/>	065HL1	E0951	Heel Loop \$30.00
<input type="checkbox"/>	065HL5	K0038	Leg Strap \$35.00
<input type="checkbox"/>	065HL4		Omit Heel Loop NC
<input type="checkbox"/>	065HL6		Omit Leg Strap NC
<input type="checkbox"/>	065IG		Impact Guard Neoprene \$85.00

Seat Sling to Top of Footrest

60° Swing- Away	Short Ext.	11 ^{1/2} - 15 ^{1/4}
60° Swing- Away	Medium Ext.	14 - 17 ^{3/4}
60° Swing- Away	Long Ext.	16 ^{1/2} - 20 ^{1/4}
60° Platform Flip-Up	Short Ext.	11 ^{1/2} - 15 ^{7/8}
60° Platform Flip-Up	Medium Ext.	14 - 18
60° Platform Flip-Up	Long Ext.	16 ^{1/2} - 20 ^{1/2}
70° Swing- Away	Short Ext.	11 ^{1/2} - 15 ^{1/4}
70° Swing- Away	Medium Ext.	14 - 17 ^{3/4}
70° Swing- Away	Long Ext.	16 ^{1/2} - 20 ^{1/4}
70° Platform Flip-Up	Short Ext.	11 ^{1/8} - 14 ^{7/8}
70° Platform Flip-Up	Medium Ext.	15 ^{3/4} - 18 ^{1/2}
70° Platform Flip-Up	Long Ext.	17 ^{3/4} - 20 ^{1/2}

*All adjustment ranges may vary due to seat height restrictions. Recommended 2^{1/2}" footrest clearance to ground surface.

Footplates			
<i>Heel loop or leg strap not included, See Footplate Options Section</i>			
<input type="checkbox"/>	065F1		Composite (14"-20"w) Only available w/ 60°, 70°, and ALR NC
<input type="checkbox"/>	065F48	K0040	Platform Flip Up \$190.00 <i>Available w/ 60°, 70°, or hemi hangers</i>
<input type="checkbox"/>	065F5		Aluminum Footplate NC
<input type="checkbox"/>	065F3	K0040	Aluminum Angle Adj Adult \$185.00
<input type="checkbox"/>	065F48	K0040	Aluminum Angle Adj Kids \$190.00 <i>NA w/ ELR. Only on 12"-16" wide</i>
<input type="checkbox"/>	065F7		Omit Footplate & Ext Tube -\$15.00

5. REAR SEAT HEIGHT SET UP

Axle			
<input type="checkbox"/>	065AX19	Threaded	NC
<input type="checkbox"/>	065AX3	Quick Release Steel	NC
<input type="checkbox"/>	065AX13	Quad Quick Release	\$185.00

One Arm Drive			
<i>Only available w/ Mag Wheels. Select Standard Rim and either right or left mount</i>			
<input type="checkbox"/>	065AL16	E0958 One Arm Drive Adj. Axle Plate	\$850.00
<input type="checkbox"/>	065OA2	Standard Handrim	NC
<input type="checkbox"/>	065HM1	Right Hand Mount	NC
<input type="checkbox"/>	065HM2	Left Hand Mount	NC

Rear Wheel Type			
<i>Contact Built 4 Me for Spinergy Color Options</i>			
<input type="checkbox"/>	065RW3	20" Mag Wheel	NC
<input type="checkbox"/>	065RW8	22" Mag Wheel	NC
<input type="checkbox"/>	065RW1	24" Mag Wheel	NC
<input type="checkbox"/>	065RW4	20" Spoke Wheel	\$200.00
<input type="checkbox"/>	065RW9	22" Spoke Wheel	\$200.00
<input type="checkbox"/>	065RW2	24" Spoke Wheel	NC
<input type="checkbox"/>	065RW11	26" Spoke Wheel	\$65.00
<input type="checkbox"/>	065RW95	24" Spinergy SPOX Black	\$795.00
<input type="checkbox"/>	065RW96	25" Spinergy SPOX Black	\$795.00
<input type="checkbox"/>	065RW97	26" Spinergy SPOX Black	\$795.00
<input type="checkbox"/>	065RW98	24" Spinergy LX Black	\$995.00
<input type="checkbox"/>	065RW99	25" Spinergy LX Black	\$995.00
<input type="checkbox"/>	065RW50	24" Mountain	\$250.00
<input type="checkbox"/>	065RW25	Omit Rear Wheel and Axle (Intended Size for Set Up _____)	-\$55.00

Rear Wheel Tire			
<input type="checkbox"/>	065RT1	E2211 E2212	Pneumatic NC
<input type="checkbox"/>	065RT2	E2211 E2213	Pneumatic Airless Insert \$80.00
<input type="checkbox"/>	065RT3	E2218	Full Poly NC
<input type="checkbox"/>	065RT6	E2211 E2212	High-Pressure Clincher \$95.00
<input type="checkbox"/>	065RT21	E2211 E2212	Schwalbe Marathon Plus Evolution \$125.00
<input type="checkbox"/>	065RT38	E2218	SHOX Performance Poly \$95.00
<input type="checkbox"/>	065RT12	E2211 E2212	Mountain Tire <i>Only available on 24" Mountain Rear Wheel</i> \$95.00

Tire	Wheel	Mag	Lite Spoke	Ultra Lightweight	Spinergy
Pneumatic		20,22,24	20,22,24,25, 26		24,25,26
Pneumatic w/Airless Insert		20,22,24	20,22,24,25,26		24,25,26
Full Poly Tire		20,22,24	20,22,24		24
High Pressure Clincher			22,24,25,26	22,24	24,25,26
Schwalbe Marathon			24,25,26	22,24	24,25,26
SHOX Performance Poly		24	24,25		24,25

22" = 489 mm 24" = 540 mm 25" = 559 mm 26" = 590 mm

Handrim			
<i>Natural Fits available in 22", 24", and 25".</i>			
<input type="checkbox"/>	065HR1	Aluminum Anodized	NC
<input type="checkbox"/>	065HR4	Plastic Coated -NA on all rims	\$95.00
<input type="checkbox"/>	065HR22	Natural Fit - Standard Grip	\$375.00
<input type="checkbox"/>	065HR23	Natural Fit - Super Grip	\$410.00
<input type="checkbox"/>	065HR7	Omit Hand Rims	NC
<input type="checkbox"/>	065HR8	8 Oblique <i>Only available w/ aluminum handrims. NA w/ Kik</i>	\$325.00

6. BACKREST

Back Type			
<i>Angle Adj std w/ push handles and NA w/ all flip back armrests</i>			
<input type="checkbox"/>	065BU1	Standard	NC
<input type="checkbox"/>	065BT20	Angle Adj 3-5" Pivot (<i>regular</i>)	\$185.00
<input type="checkbox"/>	065BT22	Angle Adj 5-7" Pivot (<i>tall</i>)	\$185.00

Angle Adjustable Backrest Options					
Angle Adjustable Pivot Height					
<input type="checkbox"/>	065AH1	3" NC	<input type="checkbox"/>	065AH4	6" NC
<input type="checkbox"/>	065AH2	4" NC	<input type="checkbox"/>	065AH5	7" NC
<input type="checkbox"/>	065AH3	5" NC			
Angle Adjustable Pre-Set Angle					
<input type="checkbox"/>	065BD5	Back Angle - 15°		NC	
<input type="checkbox"/>	065BD4	Back Angle - 12°		NC	
<input type="checkbox"/>	065BD3	Back Angle - 6°		NC	
<input type="checkbox"/>	065BD2	Back Angle - 0°		NC	
<input type="checkbox"/>	065BD1	Back Angle +3°		NC	
Angle Adjustable Back Tube Height (See Chart)					
<input type="checkbox"/>	065B1	Low		NC	
<input type="checkbox"/>	065B5	Medium		NC	
<input type="checkbox"/>	065B9	Tall		NC	
Angle Adjustable Back Height (Optional)					
<input type="checkbox"/>	065UX10	12" NC	<input type="checkbox"/>	065UX20	17" NC
<input type="checkbox"/>	065UX12	13" NC	<input type="checkbox"/>	065UX22	18" NC
<input type="checkbox"/>	065UX14	14" NC	<input type="checkbox"/>	065UX24	19" NC
<input type="checkbox"/>	065UX16	15" NC	<input type="checkbox"/>	065UX26	20" NC
<input type="checkbox"/>	065UX18	16" NC	<input type="checkbox"/>	065UX28	21" NC

JAY Back (Optional)			
<input type="checkbox"/>	065JB01	J3 Back (See J3 Order Form)	Price List

Angle Adjustable Backrest Set Up

TO DETERMINE THE PROPER PIVOT-POINT HEIGHT, MEASURE BETWEEN THE PSIS AND THE TOP OF THE SEATING SURFACE (A), THEN ADD THE THICKNESS OF THE CUSHION (B) IF ONE IS USED. IF A DROP-SEAT FALLS BELOW THE SEAT RAILS, SUBTRACT THIS DISTANCE FROM THE OVERALL VALUE.

Angle Adjustable Backrest Chart			
Back Tube Height	Low	Medium	Tall
Pivot regular set at 3"	12" 13" 14"	15" 16" 17"	17" 18" 19"
Pivot regular set at 4"	13" 14" 15"	16" 17" 18"	18" 19" 20"
Pivot regular set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot tall set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot tall set at 6"	15" 16" 17"	18" 19" 20"	N/A
Pivot tall set at 7"	16" 17" 18"	19" 20" 21"	N/A

Adjustable pivot point present measurements from seat sling. If back height option is not chosen, middle option will be present.

6. BACKREST (CONT)

Backrest Upholstery

<input type="checkbox"/>	065BU1	Standard Upholstery	NC
<input type="checkbox"/>	065BU60	Tension Adjustable Upholstery	\$175.00
<input type="checkbox"/>	065BU14	Omit Back Upholstery	-\$15.00

Backrest Options

<input type="checkbox"/>	065BU16	Reinforced Back <i>Only available w/ Std Upholstery</i>	\$30.00
<input type="checkbox"/>	065SZ1	Stabilizer Bar	\$105.00

Push Handles

<input type="checkbox"/>	065BP1	Push Handles Straight	NC
<input type="checkbox"/>	065BP3	Push Handles 8" (Std)	\$25.00

Back Height

<input type="checkbox"/>	065B92	16"-19"	NC
<input type="checkbox"/>	065B86	14"-17"	\$50.00

7. ARMREST

Armrest

<input type="checkbox"/>	065AA16	E0973	Padded Swing-Away Adult	\$210.00
<input type="checkbox"/>	065AA17	E0973	Padded Swing-Away Kids	\$210.00
<input type="checkbox"/>	065AA5		Flip-Back	NC
<input type="checkbox"/>	065AA8	E0973	Single Post Ht. Adjustable Adult	\$285.00
<input type="checkbox"/>	065AA9	E0973	Single Post Ht. Adjustable Kids	\$285.00
<input type="checkbox"/>	065AA6	E0973	Flip-Back Height Adjustable	\$295.00
<input type="checkbox"/>	065AA14		Omit Armrest	NC

Armpad

<i>NA w/ padded swing away armrests</i>			
<input type="checkbox"/>	065AM1	Desk Length - (10")	NC
<input type="checkbox"/>	065AM2	Full Length - (14")	NC

8. KEY OPTIONS & ACCESSORIES

Anti-Tip

<input type="checkbox"/>	065AT4	E0971	Rear Anti-Tip	\$140.00
--------------------------	--------	-------	---------------	----------

Wheel Lock

Additional no gos may exist depending on wheel size, seat depth and axle position

<input type="checkbox"/>	065WL1		High Mount Push	NC
<input type="checkbox"/>	065WL2		High Mount Pull	NC
<input type="checkbox"/>	065WL30		Compact	\$125.00
<input type="checkbox"/>	065WL10		Omit	NC
<input type="checkbox"/>	065WL12		Do Not Mount	NC

Grade Aid

<input type="checkbox"/>	065GA1	E0974	Grade Aid <i>NA w/ poly, high-pressure clincher tires</i>	\$150.00
--------------------------	--------	-------	---	----------

Extension Handle

<input type="checkbox"/>	065EH1	E0961	6" Pair	\$75.00
--------------------------	--------	-------	---------	---------

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

<input type="checkbox"/>	065JC20	JAY Bundled Cushion		
		Cushion Name or Part #		Price List
		Width: _____	Depth: _____	

Cushion Options

<input type="checkbox"/>	065CU1		2" Cushion	\$120.00
<input type="checkbox"/>	065CU4		Omit Seat Sling	-\$15.00

Positioning Belts

<input type="checkbox"/>	065A142	E0978	Belt 1.5" Auto Buckle	\$70.00
<input type="checkbox"/>	065A12	E0978	Belt 2" Aircraft Buckle	\$110.00
<input type="checkbox"/>	065A13	E0978	Belt 2" Aircraft Buckle Pad	\$120.00

Accessories

<input type="checkbox"/>	065A1		Quickie Backpack Black	\$55.00
<input type="checkbox"/>	065A2		Quickie Kids Backpack Black	\$50.00
<input type="checkbox"/>	065A8		Quickie Seat Pouch Black	\$45.00
<input type="checkbox"/>	065A29	E0950	Tray Table	\$250.00
<input type="checkbox"/>	065A23	K0065	Black Spoke Guards	\$115.00
<input type="checkbox"/>	065A24	K0065	Quickie Clear Spoke Guards	\$115.00
<input type="checkbox"/>	065A54	K0038	Extra Leg Strap Black	\$30.00
<input type="checkbox"/>	065SG4	K0108	Side Guard Plastic Kids	\$125.00
<input type="checkbox"/>	065SG3	K0108	Side Guard Plastic Regular	\$125.00
<input type="checkbox"/>	065A82		Caddy	\$45.00
<input type="checkbox"/>	065A49	E2207	Crutch Holder	\$195.00
<input type="checkbox"/>	065A83		Tool Kit	\$45.00
<input type="checkbox"/>	065A500		Velcro Gel Cover for Hanger/Frame	\$105.00

Military Service Patches (Optional)

Only available on tension adjustable upholstery and J3 Back

<input type="checkbox"/>	065VU1		Marines Patch	\$90.00
<input type="checkbox"/>	065VU2		Army Patch	\$90.00
<input type="checkbox"/>	065VU3		Navy Patch	\$90.00
<input type="checkbox"/>	065VU4		Air Force Patch	\$90.00
<input type="checkbox"/>	065VU12		Coast Guard Patch	\$90.00
<input type="checkbox"/>	065VU5		Purple Heart Patch	\$90.00
<input type="checkbox"/>	065VU8		Purple Heart with Marines Patch	\$145.00
<input type="checkbox"/>	065VU9		Purple Heart with Army Patch	\$145.00
<input type="checkbox"/>	065VU10		Purple Heart with Navy Patch	\$145.00
<input type="checkbox"/>	065VU11		Purple Heart with Air Force Patch	\$145.00
<input type="checkbox"/>	065VU13		Purple Heart with Coast Guard Patch	\$145.00

8. KEY OPTIONS & ACCESSORIES (CONT)

Standard Frame Colors			
<input type="checkbox"/>	065C140	Aztec Gold	NC
<input type="checkbox"/>	065C2	Black <i>(Not avail as base w/Zebra)</i>	NC
<input type="checkbox"/>	065C141	Black Cherry	NC
<input type="checkbox"/>	065C10	Black Opal	NC
<input type="checkbox"/>	065C38	Blue Opal	NC
<input type="checkbox"/>	065C36	Candy Blue	NC
<input type="checkbox"/>	065C14	Candy Purple	NC
<input type="checkbox"/>	065C35	Candy Red	NC
<input type="checkbox"/>	065C142	Electric Blue	NC
<input type="checkbox"/>	065C143	Evergreen	NC
<input type="checkbox"/>	065C126	Glow	NC
<input type="checkbox"/>	065C144	Green Apple	NC
<input type="checkbox"/>	065C147	Hot Sparkle Pink	NC
<input type="checkbox"/>	065C148	Mauve Pink	NC
<input type="checkbox"/>	065C103	Rootbeer	NC
<input type="checkbox"/>	065C114	Sparkle Silver	NC
<input type="checkbox"/>	065C145	Sunrise Orange	NC
<input type="checkbox"/>	065C113	Titanium Color Paint	NC
<input type="checkbox"/>	065C8	Yellow	NC

Matte Frame Colors			
<input type="checkbox"/>	065C116	Matte Black	NC
<input type="checkbox"/>	065C151	Matte Black Cherry	NC
<input type="checkbox"/>	065C150	Matte Electric Blue	NC
<input type="checkbox"/>	065C149	Matte Evergreen	NC
<input type="checkbox"/>	065C152	Matte Purple	NC

Custom Frame Color			
<input type="checkbox"/>	065C112	Desert Camouflage	\$295.00
<input type="checkbox"/>	065C43	Mossy Oak Camouflage	\$295.00
<input type="checkbox"/>	065C42	Stars & Stripes	\$295.00
<input type="checkbox"/>	065C108	Zebra <i>(Choose above base color) Std with White unless above base paint is chosen NA w/ Black</i>	\$295.00

Built for Me

To inquire about any unique customization for Quickie products, contact: Phone: (800) 333 - 4000, Fax: (800) 333 9011

