



Yes, you can.®

### Invacare eForms offer new functionality!

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

\*Saving and emailing the eForms is available with Acrobat Professional software

#### If you are using Adobe Reader

- **View** and download all forms
- **Save** blank versions of the forms to your desktop
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service

**Please note:** The 'Email' button at the end of the form will not work with Adobe Reader. The 'Email' button will only work while using Adobe Acrobat Professional.

#### If you are using Adobe Acrobat Professional

- **View** and download all forms
- **Save** blank AND completed versions of the forms to maintain electronic records
- **Complete** forms electronically
- **Print** blank or completed forms
- **Fax** completed forms to Customer Service
- **Email** completed forms to Customer Service by using the 'Email' button at the end of the form.

For information on downloading Adobe Reader for free, or purchasing Adobe Acrobat Professional, please visit [www.adobe.com](http://www.adobe.com).

# Invacare® MyOn™ Ultralightweight Folding Wheelchair

## Price List and Order Form

www.invacare.com

Customer Service: 800-333-6900

Fax: 800-678-4682

To order literature: 800-828-6282



Yes, you can.®

### TRANSACTION TYPE:

- Order       Quote

Date \_\_\_\_\_ Account # \_\_\_\_\_ P.O. # \_\_\_\_\_

Company Name \_\_\_\_\_ Phone # \_\_\_\_\_

Address \_\_\_\_\_ City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_ Country \_\_\_\_\_

Contact Name \_\_\_\_\_ Comments \_\_\_\_\_

### Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

## FRAME SELECTIONS

### FRAME TYPE (275 lbs weight capacity)

- MYONADLT Myon ..... \$1,799.00 K0005  
 MYONTS Myon w/Transport Brackets ..... \$2,049.00 K0005

### TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets ..... \$250.00 K0108

#### TRANSPORT TIE-DOWN\*

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

\* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

### HEAVY DUTY OPTION (350 lbs weight capacity)

HD frame includes double crossbraces and reinforced frame. See below chart for HD frame size availability.

- HDPG Heavy Duty Option ..... \$100.00

#### HEAVY DUTY - Frame Size

| Width | Depth |      |      |      |      |
|-------|-------|------|------|------|------|
|       | 16"   | 17"  | 18"  | 19"  | 20"  |
| 15"   | X     | X    | X    | X    | X    |
| 16"   | X     | X    | X    | X    | X    |
| 17"   | X     | X    | X    | X    | X    |
| 18"   | X     | X    | X    | X    | X    |
| 19"   | X     | X    | X    | X    | X    |
| 20"   | X     | X    | X    | N/A  | X    |
| 22"   | HDPG  | HDPG | HDPG | HDPG | HDPG |
| 24"   | N/A   | N/A  | HDPG | HDPG | HDPG |

X: HDPG not available; N/A: Configuration not available

### SEAT WIDTH

- W15 15" Wide ..... No Charge  
 W16 16" Wide ..... No Charge  
 W17 17" Wide ..... No Charge  
 W18 18" Wide ..... No Charge  
 W19 19" Wide ..... No Charge  
 W20 20" Wide ..... \$300.00 E2201  
 W22 22" Wide ..... \$300.00 E2201  
     W22 requires HDPG.  
 W24 24" Wide ..... \$300.00 E2202  
     W24 requires HDPG. W24 not available with D16 and D17.

### SEAT DEPTH

- D16 16" Seat Depth ..... No Charge  
 D17 17" Seat Depth ..... No Charge  
 D18 18" Seat Depth ..... No Charge  
 D19 19" Seat Depth ..... No Charge  
 D20 20" Seat Depth ..... No Charge E2203

### FRONT SEAT-TO-FLOOR HEIGHT

Refer to FSTF, Casters and Forks Compatibility table.

- FSTF14 14" Seat-to-Floor Height ..... No Charge  
     Not available with HDPG.  
     Not available with ELRSS. Not available with 1PCFTPT.  
 FSTF15 15" Seat-to-Floor Height ..... No Charge  
     Not available with HDPG. Not available with ELRSS.  
 FSTF16 16" Seat-to-Floor Height ..... No Charge  
     Not available with HDPG.  
 FSTF17 17" Seat-to-Floor Height ..... No Charge  
     Not available with HDPG.  
 FSTF18 18" Seat-to-Floor Height ..... No Charge  
 FSTF19 19" Seat-to-Floor Height ..... No Charge  
 FSTF20 20" Seat-to-Floor Height ..... No Charge

NOTE: All specifications and dimensions are approximate.

For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).

# Invacare® MyOn™ Ultralight Folding Wheelchair

## REAR SEAT-TO-FLOOR HEIGHT

Maximum 2" difference from Front Seat-to-Floor Height

- RSTF13 13" Seat-to-Floor Height ..... No Charge
  - RSTF14 14" Seat-to-Floor Height ..... No Charge
  - RSTF15 15" Seat-to-Floor Height ..... No Charge
  - RSTF16 16" Seat-to-Floor Height ..... No Charge
  - RSTF17 17" Seat-to-Floor Height ..... No Charge
  - RSTF18 18" Seat-to-Floor Height ..... No Charge
  - RSTF19 19" Seat-to-Floor Height ..... No Charge
  - RSTF20 20" Seat-to-Floor Height ..... No Charge
- Cannot select FSTF14 and RSTF15

## CENTER OF GRAVITY SETTING

Please select a center of gravity setting. Center of gravity is measured from the front of the back post to the center of the axle. The recommended setting is ACG0.

- ACG2 Center Gravity +2 (0 is back canes) ..... No Charge  
Not available with HDPG.
- ACG1 Center Gravity +1 (0 is back canes) ..... No Charge  
Not available with HDPG.
- ACG0 Center Gravity 0 (0 is back canes) ..... No Charge  
Not available with TRBKTS.
- ACG-1 Center Gravity -1 (0 is back canes) ..... No Charge  
Not available with TRBKTS.

## FRAME FINISH SELECTIONS

- 24P Wet Black ..... No Charge
- 60P Silver Metallic ..... No Charge
- 115P Black Prism ..... No Charge
- 105P Midnight Blue ..... No Charge
- 30P Sunny Yellow ..... No Charge
- 61P Electric Red ..... No Charge
- 119P Bubblegum Pink ..... No Charge
- 120P Grape Madness ..... No Charge
- 62P Electric Blue ..... No Charge
- 127P Grasshopper Green ..... No Charge
- 117P Tangerine Red ..... No Charge
- 116P Sunkissed Orange ..... No Charge
- 29P White Snow ..... No Charge

## CASTER SELECTIONS

Refer to FSTF, Casters and Forks Compatibility table.

### CASTER OPTIONS (Pair)

- 5ALUM 5" Aluminum Urethane Caster ..... No Charge  
Not available with HDPG. Not available with FRKS3 fork option.
- AT901 5" Composite Urethane ..... \$64.00  
Not available with HDPG.
- AT5SOFT 5"x 1 1/2" Wide Soft Roll ..... No Charge  
Not available with HDPG.
- AT902 6" Composite Urethane 3 Spokes ..... No Charge  
Not available with HDPG.
- 6PNEU 6" x 1 1/4" Composite Pneumatic ..... \$64.00  
Not available with HDPG.
- AT903 8" x 1 1/4" Composite Urethane Casters ..... \$78.00
- AT1292 8" x 1 1/4" Composite Pneu Casters ..... \$86.00
- 1492 8" x 1 1/4" Pneu w/ Flat Free Insert ..... \$182.00

## FORK SELECTIONS

Refer to FSTF, Casters and Forks Compatibility table.

### FORK OPTIONS (Pair)

- FRKL1 Super Low, 1 Hole Fork ..... No Charge
- FRKS1 Normal Range 1 Hole Fork ..... No Charge
- FRKS2 2 Hole Standard Fork ..... No Charge
- FRKS3 3 Hole Standard Fork ..... No Charge

## FSTF, Casters, and Forks Compatibility

| FSTF   | CASTER OPTIONS   | FORK OPTIONS |
|--------|--|--------------|
| FSTF14 | 5ALUM, AT901   | FRKL1        |
| FSTF15 | 5ALUM, AT901   | FRKL1, FRKS1 |
| FSTF16 | 5ALUM, AT901   | FRKS1        |
|        | 5ALUM, AT901, AT5SOFT, AT902, 6PNEU                      | FRKS2, FRKS3 |
| FSTF17 | 5ALUM, AT901, AT5SOFT, AT902, 6PNEU                      | FRKS2, FRKS3 |
| FSTF18 | 5ALUM, AT901, AT5SOFT, AT902, 6PNEU, AT903, AT1292, 1492 | FRKS2, FRKS3 |
| FSTF19 | 5ALUM, AT901, AT5SOFT, AT902, 6PNEU, AT903, AT1292, 1492 | FRKS2, FRKS3 |
| FSTF20 | AT903, AT1292, 1492                                      | FRKS3        |

## REAR WHEEL SELECTIONS

Please note, chair must be shipped with a pair of wheels - if a choice is not made we will select the most common No Charge option. Refer to Rear Wheels, Tires and Handrims Compatibility Table.

### WHEEL OPTIONS (Pair)

- 20COM 20" Composite Wheels ..... No Charge
- 22COM 22" Composite Wheels ..... No Charge  
Not available with FSTF14 or RSTF14
- 24COM7 24" Composite - 7 Spoke Wheels ..... No Charge
- SPIN24 24" Spinergy Spoke Wheels ..... \$695.00  
Not available with HDPG.
- SPIN25 25" Spinergy Spoke Wheels ..... \$695.00  
Not available with HDPG.
- SPIN26 26" Spinergy Spoke Wheels ..... \$695.00  
Not available with HDPG.
- 22ITESPK 22" Lightweight Performance Stylized Wheels .... No Charge  
Not available with HDPG.  
Not available with FSTF14 or RSTF14
- 24ITESPK 24" Lightweight Performance Stylized Wheels .... No Charge  
Not available with HDPG.
- 25ITESPK 25" Lightweight Performance Stylized Wheels .... No Charge  
Not available with HDPG.
- 26ITESPK 26" Lightweight Performance Stylized Wheels .... No Charge  
Not available with HDPG.

## FSTF & Rear Wheels Compatibility

| FSTF   | REAR WHEELS OPTIONS   |
|--------|---|
| FSTF14 | 20COM, 22COM  |
| FSTF15 | 20COM, 22COM, 24COM7, FN1624, 22ITESPK, 24ITESPK                              |
| FSTF16 | ALL REAR WHEELS AVAILABLE   |
| FSTF17 | ALL REAR WHEELS AVAILABLE   |
| FSTF18 | ALL REAR WHEELS AVAILABLE   |
| FSTF19 | 22COM, 24COM7, SPIN24, SPIN25, SPIN26, 22ITESPK, 24ITESPK, 25ITESPK, 26ITESPK |
| FSTF20 | 24COM7, SPIN24, SPIN25, SPIN26, 24ITESPK, 25ITESPK, 26ITESPK                  |

## SPINERGY SPOKE COLOR

Clear

If color is not chosen, spokes will be black.

- SPCBK Black ..... Standard
- SPCY Yellow ..... No Charge
- SPCW White ..... No Charge
- SPCR Red ..... No Charge
- SPCBL Blue ..... No Charge

Please note for all colors below there will be a 15 day lead time due to special order

**NOTE: All specifications and dimensions are approximate.**

For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).

# Invacare® MyOn™ Ultralight Folding Wheelchair

## TIRE OPTIONS (Pair)

Refer to Rear Wheels, Tires and Handrims Compatibility Table.

- PNEU Pneumatic Tires ..... \$60.00
- PNEUFF Flat Free Inserts ..... \$65.00 E2213
- URETH Urethane (treaded) Tires ..... No Charge
- PRIMO Primo Tires ..... \$60.00
- KIKB KIK® Black Tires..... No Charge
- CRF123 Schwalbe Marathon Plus ..... \$90.00 E2213

## HANDRIM OPTIONS (Pair)

Refer to Rear Wheels, Tires and Handrims Compatibility Table.

- ALUMPC Aluminum Anodized Handrim ..... No Charge
- A1512 Black Plastic Coated Handrims ..... \$72.00
- NFSDLT Natural-Fit® LT Handrims Standard Grip ..... \$325.00
- NFITSD Natural-Fit® Handrims Standard Grip ..... \$295.00
- AT1224-8 Plastic Coated w/8 Oblique Projections ..... \$240.00
- AT1225-8 Plastic Coated w/8 Vertical Projections ..... \$240.00
- AT1225-12 Plastic Coated w/12 Vertical Projections ..... \$240.00

## MOUNTING HANDRIMS (Spoke Wheels Only)

- OUT Mount Outside Handrim ..... Standard
- IN Mount Inside Handrim ..... No Charge

### Rear Wheels, Tires & Handrims Compatibility

| REAR WHEEL | TIRE OPTIONS                      | HANDRIM OPTIONS  |
|------------|-----------------------------------|--|
| 20COM      | PNEU, PNEUFF, URETH               | ALUMPC, A1512  |
| 22COM      | PNEU, PNEUFF, URETH               | ALUMPC, A1512, AT1225-8                                      |
| 24COM7     | PNEU, PNEUFF, URETH               | ALUMPC, A1512, NFITSD, NFSDLT, AT1224-8, AT1225-8, AT1225-12 |
| SPIN24     | PRIMO, CRF123                     | ALUMPC, A1512, NFITSD, NFSDLT                                |
| SPIN25     | PRIMO, CRF123                     | ALUMPC, A1512, NFITSD, NFSDLT                                |
| SPIN26     | PRIMO, CRF123                     | ALUMPC, A1512, NFITSD, NFSDLT                                |
| 22ITESPK   | PNEU, PNEUFF                      | ALUMPC, A1512  |
| 24ITESPK   | PNEU, PNEUFF, PRIMO, KIKB, CRF123 | ALUMPC, A1512, NFITSD, NFSDLT                                |
| 25ITESPK   | PRIMO, KIKB, CRF123               | ALUMPC, A1512, NFITSD, NFSDLT                                |
| 26ITESPK   | PNEU, PRIMO, CRF123               | ALUMPC, A1512, NFITSD, NFSDLT                                |

22COM will be excluded if TRBKTS with RSTF14 or RSTF15 is selected

## ONE-ARM DRIVE

Dual Handrim, 24" Spoke, Urethane/Pneu Tire only.

- DUHARU Right Side One Arm Drive 24" Ureth ..... \$790.00 E0958
- DUHALU Left Side One Arm Drive 24" Ureth ..... \$790.00 E0958
- DUHARP Right Side One Arm Drive 24" Pneu ..... \$790.00 E0958
- DUHALP Left Side One Arm Drive 24" Pneu ..... \$790.00 E0958



One Arm Drive

## CAMBER ANGLE

- C00 0° Camber ..... No Charge
- C03 3° Camber ..... No Charge

## REAR AXLE POSITION (Pair)

- A1577 Quick Release Axles with Dome ..... No Charge
- PAXLE Permanent Axles ..... No Charge  
Not available with PNEU pneumatic tires.
- A1573 Quad Release Axles ..... \$89.00  
Not available with 20COM, 22COM, 24COM7 rear wheels options.

## WHEEL LOCKS (Clamp On)

Please note: Chair must be shipped with a pair of wheel locks - if a choice is not made we will select the most common No Charge option

- WLPULL Pull to Lock Fold Down Wheel Locks ..... No Charge
- WLPUSH Push to Lock Fold Down Wheel Locks ..... No Charge
- WLPULLX Pull to Lock Wheel Lock w/ Extension ..... \$60.00 E0961
- WLPUSHX Push to Lock Wheel Lock w/ Extension ..... \$60.00 E0961

**NOTE: All specifications and dimensions are approximate.**  
For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).

# Invacare® MyOn™ Ultralight Folding Wheelchair

## BACK SELECTIONS

### FIXED HEIGHT

- FBH14 14" w/ Push Handles, 8 Degree Bend ..... No Charge
- FBH16 16" w/ Push Handles, 0 Degree Bend ..... No Charge
- FBH20 20" w/ Push Handles, 0 Degree Bend ..... No Charge
- 8BEND Add 8 Degree Bend ..... No Charge

### FIXED HEIGHT FOLDING BACK CANES

- FBH16F 16" Folding Back Canes, 8 Degree Bent ..... \$65.00 K0108
- FBH17F 17" Folding Back Canes, 8 Degree Bent ..... \$65.00 K0108
- FBH20F 20" Folding Back Canes, 8 Degree Bent ..... \$65.00 K0108

### ADJUSTABLE HEIGHT WITH PUSH HANDLES

- ABH1519PH 15"-19" Range, Straight Cane ..... No Charge

### ADJUSTABLE HEIGHT WITHOUT PUSH HANDLES

- ABH1519 15"-19" Range, Straight Cane ..... No Charge
- 8BEND Add 8 Degree Bend ..... No Charge

### ADJUSTABLE ANGLE BACK WITH PUSH HANDLES

- AABH1519 15"-19" Range, Straight Cane ..... \$165.00 K0108
- 8BEND Add 8 Degree Bend ..... No Charge

## UPHOLSTERY SELECTIONS

### BACK UPHOLSTERY TYPE

- SLIPBACK Slip on Back Upholstery ..... Standard
- DAJSTUPH Adjustable Back Upholstery ..... \$125.00
- L8601 Omit Back Upholstery ..... Less \$10.00

### SEAT UPHOLSTERY

- SAILCL Non-Adjustable ..... No Charge
- ADJTENS Adjustable Tension Seat Upholstery..... \$100.00
- L8651 Omit Seat Upholstery ..... Less \$10.00

| SEAT UPHOLSTERY Availability |                  |                  |                  |                  |                  |
|------------------------------|------------------|------------------|------------------|------------------|------------------|
| Width                        | Depth            |                  |                  |                  |                  |
|                              | 16"              | 17"              | 18"              | 19"              | 20"              |
| 15"                          | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | X                |
| 16"                          | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | X                |
| 17"                          | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj |
| 18"                          | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj |
| 19"                          | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj |
| 20"                          | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | X                | SAILCL/Depth Adj |
| 22"                          | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj |
| 24"                          | X                | X                | SAILCL/Depth Adj | SAILCL/Depth Adj | SAILCL/Depth Adj |

X: Not available

## ARM SELECTIONS

### ARM TYPE (Pair)



CAD, CAF, CNAD and CNAF  
Conventional Flip Back Arms Type



TAD and TAF  
T Arm Adjustable Height Arms Type

- CAD Conventional Flip Back Adj Hgt Desk Arms ..... \$195.00 E0973  
Not available with HDPG.  
Not available with seat depth D19 or D20.  
Not available with AABH1216 and AABH1519 back options.
- CAF Conventional Flip Back Adj Hgt Full Arms ..... \$195.00 E0973  
Not available with HDPG.  
Not available with seat depth D19 or D20.  
Not available with AABH1216 and AABH1519 back options.
- CNAD Convent Flip Back Non-Adj Hgt Desk Arms ..... \$145.00  
Not available with HDPG.  
Not available with seat depth D19 or D20.  
Not available with AABH1216 and AABH1519 back options.
- CNAF Convent Flip Back Non-Adj Hgt Full Arms ..... \$145.00  
Not available with HDPG.  
Not available with seat depth D19 or D20.  
Not available with AABH1216 and AABH1519 back options.
- TAD T Arm Adj Hgt Desk Length w/ Locking Pins .. \$225.00 E0973
- TAF T Arm Adj Hgt Full Length w/ Locking Pins .... \$225.00 E0973
- LESSARMS Less Arms ..... Less \$24.00

### ARM ACCESSORIES (Pair)

- A1571 Fendered Side Guard (Mudguard) Removable. \$122.00 K0108

## POSITIONING STRAPS/ACCESSORIES

### SEAT POSITIONING STRAPS (Each)

- A1515 Push Button Buckle Seat Pos. Seat ..... \$50.00 E0978  
Available on W11 and wider.
- A1321BK Hook & Loop Seat Pos. Seat ..... \$40.00 E0978  
Available on W11 and wider.

### BACK ACCESSORIES

- BP02 Back Pack ..... \$50.00

### SEAT ACCESSORIES

- ACTRY Acrylic Lap Tray ..... \$199.00 E0950  
Not available with any of the Arm Accessories.
- SP01 Seat Pouch ..... \$35.00

### MISCELLANEOUS ACCESSORIES

- A1307 Spoke Protector (For 20", 22", 24", 26") ..... \$90.00 K0065
- 1209 Tire Pump ..... \$35.00
- A1580 Rear Anti-Tippers (pair) ..... \$100.00 E0971
- A1580P Rear Anti-Tippers w/ Obstacle Climb Aid ..... \$100.00 E0971
- A1569 Cane Holder ..... \$45.00 E2207

**NOTE: All specifications and dimensions are approximate.**

*For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).*

# Invacare® MyOn™ Ultralightweight Folding Wheelchair

## FRONT RIGGING SELECTIONS

### FRONT RIGGING (Pair)

- 70NTAPER 70° Swingaway Non Tapered Footrests ..... No Charge
- 80NTAPER 80° Swingaway Non Tapered Footrests ..... No Charge
- 90STRGHT 90° Swing In/Out Non Tapered Footrests .. No Charge  
**90STRGHT front rigging is always colored black.**  
**Not available with FLIPFTS footplate option.**  
**Not available with FSTF18 and 8" casters (AT903, AT1292, 1492).**
- ELRSS Elevating Legrests ..... \$255.00 E0990  
**ELRSS Elevating Legrests are always colored black.**  
**Not available with W22 and W24.**



80NTAPER 80° Swingaway  
Non Tapered Footrests



ELRSS Elevating Legrests

### INDIVIDUAL FOOTPLATE TYPE

- FLIPFTS Flip-Up Composite Footplates ..... No Charge  
**Not available with 90STRGHT front rigging option or ELRSS on HD option.**
- FLDPADJ Flip-Up Adjustable Depth/Angle Footplates ..... \$165.00 K0040



FLDPADJ Flip-Up Adjustable Depth/Angle Footplates

### ONE PIECE FOOTPLATE TYPE

- IPCFTPT One Piece Adjustable Angle Flip Up Footplate . \$165.00 K0040  
**Not available with FSTF14**



IPCFTPT One Piece Adjustable Angle Flip Up  
Footplate

### FOOTREST/LEGREST ACCESSORIES

- A1337 Calf Strap (each) ..... \$24.00 K0038
- A9342 "H" Strap (each) ..... \$51.00 K0039
- HLOOP Heel Loop w/o Ankle Straps (pair) ..... \$51.00 E0951

| FOOTPLATE | FSTF            | KNEE-TO-HEEL RANGE |
|-----------|-----------------|--------------------|
| FLIPFTS   | FSTF14          | 13" - 16.25"       |
|           | FSTF15 - FSTF16 | 13.25" - 18.25"    |
|           | FSTF17 - FSTF20 | 13.25" - 18.25"    |
| FLPDADJ   | FSTF14          | 13.25" - 16.25"    |
|           | FSTF15 - FSTF16 | 13" - 18"          |
|           | FSTF17 - FSTF20 | 13" - 18"          |
| IPCFTPT   | FSTF15 - FSTF20 | 13.5" - 18.25"     |

## SEATING

### SEATING CUSHION SELECTION

Clear

- EC Invacare Absolute Cushion ..... \$104.00 E2601  
**Width sizes: 16", 18" and 20". Depth sizes: 16" and 18".**
- SD Invacare Single Density Cushion ..... \$170.00 E2601  
**Width sizes: 10" through 24". Depth sizes: 7" through 20".**
- ITSG Matrx Stabilite Cushion ..... \$320.00 E2605  
**Width sizes: 14" through 24". Depth sizes: 14" through 22".**
- ITFG Matrx Flovair Gentle Contour Cushion ..... \$390.00 E2607  
**Width sizes: 15" through 24". Depth sizes: 15" through 22".**
- ITFM Matrx Flovair Max Contour Cushion ..... \$390.00 E2607  
**Width sizes: 15" through 24". Depth sizes: 15" through 22".**

For additional seating needs please refer to the following:  
 Invacare Matrx Seating and Positioning Products Price List Form  
 #11-036.

Please contact Invacare Customer Service to ship Invacare Matrx  
 Seating options with your chair order.

Comments

**NOTE: All specifications and dimensions are approximate.**  
**For the most current pricing information, go to [www.invacare.com/homecare](http://www.invacare.com/homecare).**



Invacare Corporation

USA

One Invacare Way  
 Elyria, OH USA  
 44035-4190  
 800-333-6900

[www.invacare.com](http://www.invacare.com)

©2010 Invacare Corporation  
 All rights reserved. Trademarks are identified by the symbols™ and ®.  
 All trademarks are owned by or licensed to Invacare Corporation unless  
 otherwise noted. Freedom Designs is a registered trademark of Freedom  
 Design, Inc.