



Quickie QXi November 2012

EIF6

Mark For:

Submitting for:

Quote

Order

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

PO#: _____
ADDITIONAL SHIPPING INFORMATION
 Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.
 The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓ Option # HCPCS Definition Price ✓ Option # HCPCS Definition Price

1. FRAME STYLE & SETUP

Base Model			
<input type="checkbox"/>	196WM78	K0005 Quickie QXi (265 lbs) For 14" sizes use E1236	\$1,845.00

Frame Style			
<i>See Seat Height Matrix on page 2</i>			
<input type="checkbox"/>	196FS9	Standard Frame	NC
<input type="checkbox"/>	196FS30	Hemi Frame	NC

Transit			
<input type="checkbox"/>	196TO38	Transit System	\$315.00

Frame Width			
<input type="checkbox"/>	196W5	14"	NC
<input type="checkbox"/>	196W6	15"	NC
<input type="checkbox"/>	196W7	16"	NC
<input type="checkbox"/>	196W8	17"	NC
<input type="checkbox"/>	196W9	18"	NC
<input type="checkbox"/>	196W10	19"	NC
<input type="checkbox"/>	196W11	E2201 20"	\$450.00

Seat Depth			
<input type="checkbox"/>	196SD7	14"	NC
<input type="checkbox"/>	196SD8	15"	NC
<input type="checkbox"/>	196SD9	16"	NC
<input type="checkbox"/>	196SD10	17"	NC
<input type="checkbox"/>	196SD11	18"	NC
<input type="checkbox"/>	196SD71	19"	NC
<input type="checkbox"/>	196SD13	E2203 20"	\$475.00



Quickie QXi, Ultra Lightweight Adjustable Folding Wheelchair



Quickie QXi featuring the Jay ION™ cushion

2. FRONT SEAT HEIGHT SET UP

Front Seat Height

See Front Seat Height Matrix. 12.5 & 13.0 only available w/ 60° hanger
12.5" & 13" only available w/ Alum Angle Adj Footplate

<input type="checkbox"/>	196SH34	12.5"	NC	<input type="checkbox"/>	196SH14	16.0"	NC
<input type="checkbox"/>	196SH35	13.0"	NC	<input type="checkbox"/>	196SH15	16.5"	NC
<input type="checkbox"/>	196SH40	13.5"	NC	<input type="checkbox"/>	196SH16	17.0"	NC
<input type="checkbox"/>	196SH41	14.0"	NC	<input type="checkbox"/>	196SH17	17.5"	NC
<input type="checkbox"/>	196SH42	14.5"	NC	<input type="checkbox"/>	196SH18	18.0"	NC
<input type="checkbox"/>	196SH12	15.0"	NC	<input type="checkbox"/>	196SH20	18.5"	NC
<input type="checkbox"/>	196SH13	15.5"	NC	<input type="checkbox"/>	196SH21	19.0"	NC
				<input type="checkbox"/>	196SH23	19.5"	NC

Caster Wheels

<input type="checkbox"/>	196CW3	4" Low Profile Polyurethane	NC
<input type="checkbox"/>	196CW44	4" x 1.25" Semi Pneumatic	\$55.00
<input type="checkbox"/>	196CW5	5" Low Profile Polyurethane	NC
<input type="checkbox"/>	196CW7	5" x 1.5" Semi Pneumatic	\$85.00
<input type="checkbox"/>	196CW34	6" Polyurethane	NC
<input type="checkbox"/>	196CW35	6" x 1.5" Semi Pneumatic	\$85.00
<input type="checkbox"/>	196CW32	8" Polyurethane	NC
<input type="checkbox"/>	196CW39	8" x 1.5" Semi Pneumatic	\$85.00

Caster Fork

<input type="checkbox"/>	196CF71	3.5" (standard stem bolt only)	NC
<input type="checkbox"/>	196CF12	5.25"	NC
<input type="checkbox"/>	196CF17	6"	NC
<input type="checkbox"/>	196CF23	7"	NC

Caster Stem Bolt

<input type="checkbox"/>	196CS1	Standard	NC
<input type="checkbox"/>	196CS5	3/4"	NC
<input type="checkbox"/>	196CS9	1 1/2"	NC

Front Seat Height Matrix

Caster	Standard Stem Bolt							
	3.5" FORK				5 1/4" FORK			
	STD		HEMI		STD		HEMI	
4 Low Poly	14.5	15.0	12.5	13.0	16.0	17.0	14.0	15.0
4 x 1.25 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
5 Low Poly	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
5 x 1.5 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
6 Poly	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
8 Poly	N/A	N/A	N/A	N/A	19.0	19.0	17.0	17.0
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.0	16.5	17.0
Caster	6" FORK				7" FORK			
	STD		HEMI		STD		HEMI	
	STD		HEMI		STD		HEMI	
4 Low Poly	17.5	18.0	15.5	16.0	18.5	19.0	16.5	17.0
4 x 1.25 Semi Pneu	17.5	18.0	15.5	16.0	18.5	19.0	16.5	17.0
5 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
5 x 1.5 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
6 Poly	17.5	19.0	15.5	17.0	18.5	N/A	16.5	18.0
6 x 1.5 Semi Pneu	17.5	19.0	15.5	17.0	18.5	N/A	16.5	18.0
8 Poly	19.5	N/A	17.5	18.0	19.5	N/A	17.5	19.0
8 x 1.5 Semi Pneu	18.5	19.5	16.5	17.5	19.0	N/A	17.0	18.5

*Seat squeeze and seat depth can affect FSH availability
For 3/4" and 1 1/2" Stem bolt matrix see page 5*

3. HANGER & FOOTPLATE OPTIONS

Front Hangers

Footplate mounts with Extension Tube using clamping hardware

<input type="checkbox"/>	196H145	60° Swing In/Out Hanger	NC
<input type="checkbox"/>	196H139	70° Swing In/Out Hanger	NC
<input type="checkbox"/>	196H141	80° Swing In/Out Hanger	NC
<input type="checkbox"/>	196H143	E0990 Elevating Legrest (Front Latch) <i>All ELR's will be Black paint</i>	\$295.00
<input type="checkbox"/>	196H144	Omit Hanger & Footplate	NC

Extension Tubes

See Lower Leg Length Matrix for range

<input type="checkbox"/>	196E1	Short	NC
<input type="checkbox"/>	196E2	Medium	NC
<input type="checkbox"/>	196E5	Long	NC

Footplate

Heel loop not included, see Footplate Option section below

<input type="checkbox"/>	196F1	Composite Footplate	NC
<input type="checkbox"/>	196F3	K0040 Aluminum Angle Adj. Footplate	\$185.00
<input type="checkbox"/>	196F5	Aluminum <i>Only available w/Elevating Legrest</i>	NC

Hanger & Footplate Options


<input type="checkbox"/>	196HL12	E0951 Adjustable Heel Loop	\$35.00
<input type="checkbox"/>	1968A02	Gel Cover For Hangers (Pair)	\$110.00

UPDATED

Lower Leg Length Range Matrix

Hanger	Ext Tube	Footplate		
		Composite	Alum Angle Adj	Aluminum
60°	Short	11.25" - 14.75"	9.25" - 13.25"	NA
	Medium	11.25" - 16.00"	9.25" - 14.75"	NA
	Long	12.50" - 16.00"	10.75" - 16.25"	NA
70°	Short	12.00" - 15.75"	10.00" - 14.50"	NA
	Medium	12.00" - 16.50"	10.75" - 16.50"	NA
	Long	13.75" - 16.50"	12.25" - 16.50"	NA
80°	Short	10.50" - 15.25"	10.50" - 15.25"	NA
	Medium	11.75" - 16.50"	11.75" - 16.50"	NA
	Long	13.50" - 16.50"	11.75" - 16.50"	NA
ELR	Short	11.25" - 15.50"	11.00" - 14.50"	11.50" - 16.50"
	Medium	12.50" - 16.50"	11.25" - 16.25"	13.50" - 16.50"
	Long	14.50" - 16.50"	13.00" - 16.50"	15.25" - 16.50"

All values are ± 0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full ELR values based on 70° hanger angle



4. REAR SEAT HEIGHT SET UP

Rear Seat Height Matrix

Rear Wheel Size	Std Frame	Hemi Frame
20" Rear Wheel	13.0" - 16.5"	12.0" - 15.0"
22" Rear Wheel	14.0" - 17.5"	12.5" - 16.0"
24" Rear Wheel	15.0" - 18.5"	13.5" - 17.0"
25" Rear Wheel	15.5" - 19.0"	14.0" - 17.5"
26" Rear Wheel	16.0" - 19.5"	14.5" - 18.0"

Rear Seat Height

See Rear Seat Height Matrix

<input type="checkbox"/> 196RH9	12"	NC
<input type="checkbox"/> 196RH10	12.5"	NC
<input type="checkbox"/> 196RH11	13"	NC
<input type="checkbox"/> 196RH12	13.5"	NC
<input type="checkbox"/> 196RH13	14"	NC
<input type="checkbox"/> 196RH14	14.5"	NC
<input type="checkbox"/> 196RH15	15"	NC
<input type="checkbox"/> 196RH16	15.5"	NC
<input type="checkbox"/> 196RH17	16"	NC
<input type="checkbox"/> 196RH18	16.5"	NC
<input type="checkbox"/> 196RH19	17"	NC
<input type="checkbox"/> 196RH26	17.5"	NC
<input type="checkbox"/> 196RH20	18"	NC
<input type="checkbox"/> 196RH27	18.5"	NC
<input type="checkbox"/> 196RH21	19"	NC

Seat Squeeze Matrix

Frame Depth	Max Seat Squeeze
14"	2.0"
15" - 16"	2.5"
17" - 18"	2.5"
19" - 20"	3.0"

Rear Tire & Wheel Matrix

Tire	Wheel	Mag	Spoke
Pneumatic		20,22,24,26	20,22,24,26
Pneumatic w/ Airless Insert		20,22,24,26	20,22,24,26
Full Poly Tire		20,22,24	20,22,24
High Pressure Clincher		N/A	22,24,25,26
		22" = 489 mm	24" = 540 mm

Rear Wheel Size

See Rear Tire & Wheel Matrix for availability

<input type="checkbox"/> 196RW3	20" Mag Wheel	NC
<input type="checkbox"/> 196RW8	22" Mag Wheel	NC
<input type="checkbox"/> 196RW1	24" Mag Wheel	NC
<input type="checkbox"/> 196RW90	20" Lite Spoke Wheel	\$200.00
<input type="checkbox"/> 196RW91	22" Lite Spoke Wheel	\$200.00
<input type="checkbox"/> 196RW92	24" Lite Spoke Wheel	NC
<input type="checkbox"/> 196RW100	25" Lite Spoke Wheel	\$200.00
<input type="checkbox"/> 196RW101	26" Lite Spoke Wheel	\$200.00
<input type="checkbox"/> 196RW25	Omit	NC

Rear Wheel Tire

See Rear Tire & Wheel Matrix for availability

<input type="checkbox"/> 196RT1	E2211 E2212	Pneumatic	NC
<input type="checkbox"/> 196RT2	E2211 E2213	Pneumatic w/ Airless Insert	\$80.00
<input type="checkbox"/> 196RT3	E2218	Full Poly	NC
<input type="checkbox"/> 196RT6	E2211 E2212	High-Pressure Clincher	\$95.00

Handrim

Natural Fit only available w/ 22, 24, 25 Lite Spoke. Plastic NA w/ 25, 26"

<input type="checkbox"/> 196HR1	Aluminum Anodized	NC
<input type="checkbox"/> 196HR4	Plastic	\$95.00
<input type="checkbox"/> 196HR22	Natural Fit - Standard Grip	\$375.00
<input type="checkbox"/> 196HR23	Natural Fit - Super Grip	\$395.00

Wheel Lock

<input type="checkbox"/> 196WL1	High Mount Push	NC
<input type="checkbox"/> 196WL2	High Mount Pull	NC
<input type="checkbox"/> 196WL30	Compact	\$125.00
<input type="checkbox"/> 196WL10	Omit Wheel Lock	NC
<input type="checkbox"/> 196WL12	Do Not Mount	NC

Extension Handle

Ext. handles only available with high mount push or pull wheel locks

<input type="checkbox"/> 196EH7	E0961	Extension Handle (Pair)	\$60.00
---------------------------------	-------	-------------------------	---------

Axle

<input type="checkbox"/> 196AX1	Quick Release Steel	NC
<input type="checkbox"/> 196AX11	Quad Quick Release Steel NA on Spoke	\$90.00

5. ARMREST

Armrest - Padded Swing Away

Some limitations when used with the Angle Adjustable Back in the rear most setting

Please circle a setting, will default to 11" if no selection is made

<input type="checkbox"/> 196AA18	E0973	Padded Swing-Away STD (9" to 13")	\$175.00
		Setting <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13"	

Armrest - Height Adjustable

<input type="checkbox"/> 196AA104	E0973	Single Post Ht. Adj - Std (9.75"-14.75")	\$260.00
<input type="checkbox"/> 196AA105	E0973	Single Post Ht Adj - Low (8.25"-12.25")	\$260.00
<input type="checkbox"/> 196AA65	E0973	Flip Back Ht Adj LITE (10"-15")	\$275.00
<input type="checkbox"/> 196AA14		Omit Armrest	NC
<input type="checkbox"/> 196AA103	E0973	Dual Post Flip Back Ht Adj NEW!	\$275.00

Armpad - For Height Adjustable

<input type="checkbox"/> 196AM23	Desk Length (10")	NC
<input type="checkbox"/> 196AM24	Full Length - Full (14")	NC

Armrest Options

<input type="checkbox"/> 1968A01	Gel Cover For Swing Away Arms (pair)	\$110.00
----------------------------------	--------------------------------------	----------

6. BACKREST

Back Type			
<input type="checkbox"/>	196BT10	Standard	NC
<input type="checkbox"/>	196BT100	Angle Adjustable Back	\$165.00

Backrest Height			
<input type="checkbox"/>	196B5	Med (13"-17")	NC
<input type="checkbox"/>	196B9	Tall (16"-20")	NC

Push Handles			
<input type="checkbox"/>	196PH6	Integral Push Handles 8° for Std Backrest, Straight for Angle Adj. Backrest	NC
<input type="checkbox"/>	196PH5	Fold-Down Push Handles Only available with Adjustable Upholstery. Used in combination w/Angle Adjust only avail in 16"-20"	\$100.00
<input type="checkbox"/>	196PH2	Omit Push Handles 8° for Std Backrest, Straight for Angle Adj. Backrest	NC

Back Upholstery			
<input type="checkbox"/>	196BU1	Standard Upholstery	NC
<input type="checkbox"/>	196BU70	E2611 Backrest Tension Adj Ballistic Nylon *Reference 186BU70 for billing	\$315.00
<input type="checkbox"/>	196BU71	E2611 Backrest Tension Adj 3DX Vented *Reference 186BU71 for billing	\$330.00
<input type="checkbox"/>	196BU14	Omit Back Upholstery	-\$15.00

JAY J3 Back			
<input type="checkbox"/>	196JB01	JAY J3 Back	See OF

Military Service Patches			
<i>Only available on tension adjustable upholstery and Jay 3 Back</i>			
<input type="checkbox"/>	196VU1	US Marines Patch	\$90.00
<input type="checkbox"/>	196VU2	US Army Patch	\$90.00
<input type="checkbox"/>	196VU3	US Navy Patch	\$90.00
<input type="checkbox"/>	196VU4	US Air Force Patch	\$90.00
<input type="checkbox"/>	196VU12	US Coast Guard	\$90.00
<input type="checkbox"/>	196VU5	Purple Heart	\$90.00
<input type="checkbox"/>	196VU8	Purple Heart with Marines Patch	\$145.00
<input type="checkbox"/>	196VU9	Purple Heart with Army Patch	\$145.00
<input type="checkbox"/>	196VU10	Purple Heart with Navy Patch	\$145.00
<input type="checkbox"/>	196VU11	Purple Heart with Air Force Patch	\$145.00
<input type="checkbox"/>	196VU13	Purple Heart with Coast Guard	\$145.00

7. SEATING

Seat Sling			
<input type="checkbox"/>	196SL23	Seat Sling Depth Adjustable	NC
<input type="checkbox"/>	196SL6	Omit Seat Sling	-\$15.00

NEW!		JAY Cushions (Optional)	
-------------	--	-------------------------	--

To order any JAY cushions on-chair, please see the JAY price list.
Cushion name or part number must be filled out with desired width and depth

<input type="checkbox"/>	196JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width: _____ Depth: _____	

8. Key Options & Accessories

Anti Tips			
<input type="checkbox"/>	196AT1	E0971 Anti Tips (Pair)	\$135.00

Side Guards			
<input type="checkbox"/>	196SG33	K0108 Plastic - Standard (8.5" High)	\$125.00
<input type="checkbox"/>	196SG29	K0108 Plastic - Low (6.5" High)	\$125.00

Pack Accessories			
<input type="checkbox"/>	196A1	Backpack - Black/Charcoal Adult	\$55.00
<input type="checkbox"/>	1961A01	Seat Pouch Black	\$45.00

Positioning Belt			
<input type="checkbox"/>	196A141	E0978 1.5" Auto Buckle w/ Clamps	\$70.00
<input type="checkbox"/>	196A804	K0108 Positioning Belt Slide Clamps	\$50.00

Spoke Guards			
<input type="checkbox"/>	1964A01	K0065 Quickie Black Spoke Guard	\$115.00

Miscellaneous			
<input type="checkbox"/>	1968A03	Caddy	\$45.00
<input type="checkbox"/>	196A2100	E2207 Crutch Holder	\$195.00

Frame Color			
<input type="checkbox"/>	196C2	Black	NC
<input type="checkbox"/>	196C10	Black Opal	NC
<input type="checkbox"/>	196C38	Blue Opal	NC
<input type="checkbox"/>	196C36	Candy Blue	NC
<input type="checkbox"/>	196C35	Candy Red	NC

Front Seat Height Matrix Cont.

3/4" STEM								
Caster	3.5" FORK				5 1/4" FORK			
	STD	HEMI	STD	HEMI	STD	HEMI	STD	HEMI
4 Low Poly	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
5 Low Poly	N/A	N/A	N/A	N/A	17.0	18.0	14.5	15.5
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	15.5	16.5
6 Poly	N/A	N/A	N/A	N/A	17.5	18.5	15.5	16.5
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	15.5	16.5
8 Poly	N/A	N/A	N/A	N/A	19.5	19.5	17.5	17.5
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	19.0	19.5	17.0	17.5
Caster	6" FORK				7" FORK			
	STD	HEMI	STD	HEMI	STD	HEMI	STD	HEMI
4 Low Poly	18.0	18.5	16.0	16.5	19.0	19.5	17.0	17.5
4 x 1.25 Semi Pneu	18.0	18.5	16.0	16.5	19.0	19.5	17.0	17.5
5 Low Poly	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0
5 x 1.5 Semi Pneu	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0
6 Poly	18.0	19.5	16.0	17.5	19.0	19.5	17.0	18.5
6 x 1.5 Semi Pneu	18.0	19.5	16.0	17.5	19.0	19.5	17.0	18.5
8 Poly	N/A	N/A	18.0	18.5	19.5	19.5	18.0	19.5
8 x 1.5 Semi Pneu	19.0	19.5	17.0	18.0	19.5	19.5	17.5	19.0
<i>Seat squeeze and seat depth can affect FSH availability</i>								

1 1/2" STEM								
Caster	3.5" FORK				5 1/4" FORK			
	STD	HEMI	STD	HEMI	STD	HEMI	STD	HEMI
4 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	15.5	16.5
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	15.5	16.5
5 Low Poly	N/A	N/A	N/A	N/A	18.0	19.0	16.0	17.0
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.0	19.0	16.0	17.0
6 Poly	N/A	N/A	N/A	N/A	18.5	19.5	16.5	17.5
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.5	16.5	17.5
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	18.5	18.5
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	19.5	19.5	17.5	18.0
Caster	6" FORK				7" FORK			
	STD	HEMI	STD	HEMI	STD	HEMI	STD	HEMI
4 Low Poly	19.0	19.5	17.0	17.5	N/A	N/A	18.0	18.5
4 x 1.25 Semi Pneu	19.0	19.5	17.0	17.5	N/A	N/A	17.5	18.5
5 Low Poly	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5
5 x 1.5 Semi Pneu	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0
6 Poly	18.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0
6 x 1.5 Semi Pneu	18.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0
8 Poly	N/A	N/A	18.5	19.0	N/A	N/A	18.5	19.5
8 x 1.5 Semi Pneu	N/A	N/A	18.0	19.0	N/A	N/A	18.5	19.5
<i>Seat squeeze and seat depth can affect FSH availability</i>								

