

SPECIFICATIONS

	Victory® 9	Victory® 9PS	Victory® 10
Max. weight capacity:	300 lbs.	300 lbs.	400 lbs.
Maximum speed: ¹	Up to 5.25 mph	Up to 5.5 mph	Up to 5.25 mph
Ground clearance: ²	3W: 2" / 4W: 2.25"	2" at midframe	3" at midframe
Turning radius: ²	3W: 40.25" / 4W: 51.5"	41"	3W: 42" / 4W: 54"
Overall length: ²	3W: 43.5" / 4W: 45.5"	43.5"	3W: 45.25" / 4W: 47"
Overall width: ²	3W: 22.25" 4W: 22.5"	23"	22.25"
Battery charger:	Off-board, 3 amp	Off-board, 3 amp	Off-board, 3 amp
Per-charge range: ⁵	Up to 13 miles	Up to 13 miles	Up to 15.5 miles
Tire size:	3W: Front: 8" solid/Rear: 9" solid 4W: Front: 9" solid/Rear: 9" solid	3W: Front: 8" solid/Rear: 9" solid	3W: Front: 10" solid 4W: Rear: 10" solid
Total weight without batteries:	3W: 117.5 lbs./4W: 131 lbs.	131 lbs.	3W: 131.5 lbs./4W: 132.5 lbs.
Heaviest piece: ² when disassembled:	3W: 48.5 lbs. (front section) 4W: 45 lbs. (front & rear section)	59.5 (front section)	3-wheel: 54.5 lbs. (front section) 4-wheel: 61 lbs. (front section)

¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).


³AGM or Gel-Cell type required.

⁴Battery weight may vary based on manufacturer

⁵Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.

⁶Includes standard seat weight

NOTE: All specifications subject to change without notice.



Dealer Information

Pride[®]
 Mobility Products Corp.
 — LIVE YOUR BEST —

VICTORYTRI_05-10

VICTORY[®] SERIES



Pride[®]
 Mobility Products Corp.
 — LIVE YOUR BEST —

VICTORY® 9



FEATURES

- Exclusive, one-hand Feather-touch disassembly provides the ultimate quick transport system
- Pride®'s exclusive stylish, lightweight, low profile, non-scuffing wheels
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Seat includes sliders to position oneself closer or further from tiller
- Available in 3 & 4-wheel models

VICTORY® 9



FEATURES

- Standard power-elevating seat
- Exclusive, one-hand Feather-touch disassembly provides the ultimate quick transport system
- Pride®'s exclusive stylish, lightweight, low profile, non-scuffing wheels
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Available in Candy Apple Red or Viper Blue
- Only available as a 3-wheel model

VICTORY® 10



FEATURES

- Exclusive, one-hand Feather-touch disassembly provides the ultimate quick transport system
- Crush proof basket
- Delta tiller with wraparound handles
- Pride®'s exclusive stylish, lightweight, non-scuffing, black, low profile wheels
- Backlit battery gauge
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Easy access tiller-mounted charger port
- Available as 3 & 4 wheel model