

ALL COURT



Quickie All Court – A Revolution for Court Sports

Already the world's most successful sports wheelchair the All Court can be used for any position in all court or impact sports. The chair adjusts with a patented center-of-mass adjustment feature that eliminates the need to align the wheels – a feature unique to Quickie. New for 2001, the All Court includes a six-point attachment to tightly secure the bumper at 11 cm and a new quick-release fifth wheel option for added convenience. This new lighter weight frame design offers optimal adjustability and performance for recreational players and serious athletes alike.



ALL COURT

YOUR ABILITY

