

Item #	Seat	Color	UOM
SFSCOUT3	16.5" Folding Seat	Red/Blue 	1/cs



HCPCS
K0800



Spitfire Scout 3

Compact Travel Scooter

- Quick and easy disassembly.
- Interchangeable color panels hide scuffs and scratches.
- Armrests are padded and adjustable.
- Quick connect batteries.
- Ergonomic throttle control.
- Height-adjustable swivel seat with fold-down backrest.
- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Delta tiller.
- Easy freewheel operation

specifications

HCPCS	K0080	Charger	2A Offboard
Weight Capacity	300 lbs.	Brakes	Electromagnetic
Top Speed	4.25 mph	Front Wheels	8" x 2" Flat-free
Maximum Range	9 miles	Rear Wheels	8" x 2" Flat-free
Turning Radius	45.5"	Base Weight	49 lbs.
Climbing Angle	6°	Battery Weight (pair)	20 lbs.
Ground Clearance	2.5"	Seat Weight	16 lbs.
Dimensions	42.5" (L) x 21.25" (W)	Total Weight	85 lbs.
Floor-to-Seat Height	17" - 19"	Heaviest Piece Weight (Front Section)	29 lbs.
Seat Dimensions	16.5" (W) x 16" (D)	Warranty on Frame	Lifetime
Controller	CTE 45A	Warranty on Electronic Controller/Drive Train Components	14 months
Motor	24V x 270W x 4700 rpm	Warranty on Batteries	6 months
Batteries (pair)	12V x 12AH		