

About EZ-ACCESS

EZ-ACCESS is the world's leading manufacturer of high quality aluminum ramp systems. We are a family-owned and operated company with over 30 years of experience in designing and manufacturing access solutions. Our products are made in the USA and are known for their safety, durability, and workmanship.

In addition to the versatile Suitcase® ramp line, we manufacture the Passport® Vertical Platform Lift and a variety of portable ramps and other mobility accessories for both commercial and residential use.


EZ-ACCESS can make any type of home, regardless of age or style, accessible. Our products provide flexible short- and long-term solutions for a variety of needs.

Through an extensive network of local authorized dealers, EZ-ACCESS serves the entire United States with professional products and quick installation.

For additional information, please contact us today or visit ezaccess.com.



customerservice@ezaccess.com | 1.800.451.1903 | www.ezaccess.com



Scan the QR code with your smartphone to see product specifications.



SUITCASE®

Anytime. Anywhere.
Portable ramps for your every need.



EZ-ACCESS® SUITCASE® ramps allow durable, dependable solutions to a variety of mobility and accessibility needs. From the SUITCASE Singlefold to the Singlefold AS and Trifold AS, you are sure to find a convenient, portable, and easy-to-use ramp perfect for wheelchair and scooter use.

Visit ezaccess.com for product specifications


SUITCASE® Ramps

- Ideal for use at home or on the go
- Designed for wheelchairs and scooters
- Available in 3 styles, each with their own unique features
- Made in the USA



SUITCASE® Singlefold

This portable ramp, with its single-fold design, sets up quickly and can be carried like a suitcase when not in use. Features a full-length hinge for additional strength and stability, and an extruded slip-resistant surface for superior traction. This durable ramp was designed for high traffic use and is simply the best in the industry. 800-pound weight capacity.



SUITCASE® Singlefold AS

Can be used as one unit or easily separated into two lightweight sections, making it even more convenient to use, carry, and store. This ultra-strong ramp provides the safety and strength required while remaining light and portable. Features an applied, slip-resistant surface. 800-pound weight capacity.



SUITCASE® Trifold® AS

Unique 3-fold design allows for ultimate portability when longer ramp lengths are needed. The ramp is designed to be used and carried as one unit, or can be quickly separated into two individual sections for easier transport and reduced handling weight. Features an applied, slip-resistant surface. 800-pound weight capacity.

Quick Specs

		Singlefold	Singlefold AS	Trifold AS
2'	Usable Size	2'x30"	2'x29-1/2"	—
	Folded Size	2'x15"x4"	2'x15"x4-1/2"	—
	Weight (lbs.)*	9	6 ea. (2 sections)	—
	Application	Curbs, 1 step	Curbs, 1 step	—
3'	Usable Size	3'x30"	3'x29-1/2"	—
	Folded Size	3'x15"x4"	3'x15"x4-1/2"	—
	Weight (lbs.)*	14	8 ea. (2 sections)	—
	Application	Curbs, 1 step	Curbs, 1 step	—
4'	Usable Size	4'x30"	4'x29-1/2"	—
	Folded Size	4'x15"x4"	4'x15"x4-1/2"	—
	Weight (lbs.)*	18	10 ea. (2 sections)	—
	Application	Curbs, 1-2 steps	Curbs, 1-2 steps	—
5'	Usable Size	5'x30"	5'x29-1/2"	5'x29-1/2"
	Folded Size	5'x15"x4"	5'x15"x4-1/2"	30-3/4"x15"x11-1/2"
	Weight (lbs.)*	25	13 ea. (2 sections)	16 ea. (2 sections)
	Application	Curbs, 1-2 steps, some minivans	Curbs, 1-2 steps, some minivans	Curbs, 1-2 steps, some minivans
6'	Usable Size	6'x30"	6'x29-1/2"	6'x29-1/2"
	Folded Size	6'x15"x4"	6'x15"x4-1/2"	36-3/4"x15"x11-1/2"
	Weight (lbs.)*	30	15 ea. (2 sections)	19 ea. (2 sections)
	Application	1-3 steps, some minivans & SUVs	1-3 steps, some minivans & SUVs	1-3 steps, some minivans & SUVs
7'	Usable Size	7'x30"	—	7'x29-1/2"
	Folded Size	7'x15"x4"	—	42-3/4"x15"x11-1/2"
	Weight (lbs.)*	36	—	22 ea. (2 sections)
	Application	1-3 steps, some minivans & SUVs	—	1-3 steps, some minivans & SUVs
8'	Usable Size	8'x30"	—	8'x29-1/2"
	Folded Size	8'x15"x4"	—	48-3/4"x15"x11-1/2"
	Weight (lbs.)*	42	—	25 ea. (2 sections)
	Application	2-3 steps, some minivans & SUVs	—	2-3 steps, some minivans & SUVs
10'	Usable Size	—	—	10'x29-1/2"
	Folded Size	—	—	60-3/4"x15"x11-1/2"
	Weight (lbs.)*	—	—	30 ea. (2 sections)
	Application	—	—	2-4 steps, some minivans & SUVs

Top Lip Extension (TLE)



Our optional TLE increases the reach of EZ-ACCESS SUITCASE ramps by extending the top lip from 3" to 9", allowing the ramp to clear the rear bumper of most SUVs and vans.

Be Prepared

Like a trusty umbrella or roadside emergency kit, your vehicle should always have a small wheelchair ramp tucked away for emergencies. A 2-ft. SUITCASE® Singlefold ramp is the ideal choice and its compact size allows for easy, take-anywhere portability. The lightweight ramp makes a great traveling companion because it fits perfectly in the trunk of a car, under the seat of a van, or is easily checked at the airport.

Call today for more information

customerservice@ezaccess.com | 1.800.451.1903 | www.ezaccess.com

* Weights are approximate and may vary. Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.