

  
 **TOP END**®



## Invacare® Top End® T-5™ 7000 Series Tennis Wheelchair

The **Invacare® Top End® T-5™ 7000 Series Tennis Wheelchair** is faster, lighter and stronger.\* The most recent enhancement included new material, enhanced headtube design and axle securing system, reducing overall weight by two pounds and improving performance. Each T-5 7000 Series Tennis wheelchair is custom built to each individual athlete's specifications and is overseen by Paul Schulte who is the design engineer for Top End and three-time Paralympian as a member of the USA National Men's Wheelchair Basketball Team.

\*As compared to previous Invacare® Top End® model wheelchairs.

## Features

Custom-built fixed frame in 7005 aluminum  
Adjustable height fixed swivel anti-tip  
59mm fifth wheel  
Fixed back with fixed or adjustable height  
Adjustable tension seat and back upholstery  
Fixed camber: 15°, 18° or 20°  
High performance Sun® wheels:

24", 25", 26" or 700C  
Precision fork caster system  
72mm casters  
Powder coat paint  
Rear frame protector  
250 lb. weight capacity

## Model #: T57



1. Custom front frame angle and tubular footrest
2. Adjustable height fixed swivel anti-tip
3. 7005 aluminum heat treating process



4. Precision fork caster system
5. Optional welded fixed tube sideguards
6. Adjustable center-of-gravity option

## Options

Spinergy® SLX wheels  
Adjustable Center-of-Gravity  
Welded seat rigidizer bars with solid seat  
Telescoping removable anti-tip bar  
Click straps  
Bolt-on fixed side guards  
Fixed side guards with welded tube  
Fixed aluminum or carbon fiber side guards  
Coated pushrims

Footrest covers  
Frame protectors  
Custom strapping  
Alignment gauge  
Adjustable angle footplate  
Recessed platform footplate with toe stops/footrest strap  
Close mount/Cut pushrim tab  
O-ring positioning belt

See order form for current options and complete details

**Invacare Corporation**  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

### Top End

www.topendwheelchair.com  
Sports and Recreation  
Products  
4501 63rd Circle North  
Pinellas Park, Florida  
33781-5914  
(800) 532-8677  
Fax: (727) 522-1007

©2018

Invacare Corporation.  
All rights reserved.  
Trademarks are identified  
by the symbols™ and®.  
All trademarks are owned  
by or licensed to Invacare  
Corporation unless  
otherwise noted. Spinergy  
is a registered trademark  
of Spinergy Inc. Sun and  
Equalizer are registered  
trademarks of Sun  
Components. Kenda is a  
registered trademark of  
Kenda Rubber Industrial  
Co., Ltd. Specifications are  
subject to change without  
notification.  
Form No. 18-180  
180738

 **TOP END®**