





HCPCS
E1238



Wallaby Pediatric Folding Wheelchair

Item #	Product Seat Dimensions	Product Overall Width Open	Product Cane Height	Color	UOM
WB 1200	12" (W) x 12" (D)	20"	14.5"	Red 	1/ea
WB 1400	14" (W) x 12" (D)	22"	15"	Blue 	1/ea

specifications	
Overall Width (Closed)	12"
Seat to Floor (Front)	19"
Seat to Floor (Rear)	18"
Seat to Footrest	10" - 15"
Push Handle Height	33" - 39"
Armrest Height	7"
Overall Length w/ Footrests	37"
Overall Length without Footrests	28.5"
Back Angle	90°
Product Weight w/ Footrests	37 lbs.
Product Weight w/out Footrests	33 lbs.
Weight Capacity	150 lbs.
Weight Capacity when Used in Bus Transport	80 lbs.
Warranty	Limited Lifetime

- WC19 Bus Transit Approved when used in conjunction with WB 8021 and WB 8030.
- Powder-coated colored steel frame with black cross brace.
- Easy-to-clean lightweight durable nylon upholstery.
- Height-adjustable push handles.
- Flip-back padded desk length arms.
- Pelvic belt with reversible cover.
- 22" rear and 6" front swivel solid rubber tires mounted on composite wheels.
- Push-to-lock brakes.
- Swing-away footrests with plastic plates and calf straps.

Headrest Extension



HCPCS
E0966

Item #	For Use With	UOM
WB 8021	WB 1200 & WB 1400	1/bx

- Width and height adjustable.
- Provides additional support for taller users.
- Used on bus-transit models to prevent head whiplash during rear-end collision.

Five Point Harness



HCPCS
E0960

Item #	For Use With	UOM
WB 8030	WB 1200 & WB 1400	1/bx

- Converts three-point harness to a five-point harness.
- Padded reversible strap covers feature imitation lamb's wool on one side and velour on the other side.
- Padded covers increase user's comfort and reduce friction.

Elevating Legrests



HCPCS
E0990

Item #	For Use With	UOM
WB ELR	WB 1200 & WB 1400	1 pr/bx

- Seating with raised legs assists with decreasing dependent edema and is essential for patients with below-the-knee amputation, knee extension contractures or other joint abnormalities.
- Includes a calf support for the lower leg.

Anti-Tippers



Item #	For Use With	UOM
STDS 833N	WB 1200 and WB 1400	1 pr/bx