

Zippie 2[™]





The Zippie 2 is a lightweight folding wheelchair that folds small enough to fit in most compact cars. The Zippie 2 is easily adjustable and has the ability to grow as the child grows

Folding frame allows for easier portability. **Locking crossbraces** gives the Zippie 2 an easy-to-push, solid feel.

Seat rail is eliminated to provide an **open seat frame** which easily interfaces with seating systems. Elimination of the upper seat rails allows for lower seat heights and enables children to interact at eye level.

Interchangeability in the components allows the chair to grow as the child grows. The Zippie 2 offers a wide variety of swing in/out hangers and other adjustable, quick-release components.







ZIPPIE 2 SPECIFICATIONS	
HCPCS Code:	E1238
Overall Weight:*	25 lbs. w/o footrest
Weight Capacity:	165 lbs.
Overall Product Width:	20" -24"
Transit Option:	Yes
Front Seat Height:	15" - 20"
Rear Seat Height:	13.5" - 19"
Hanger Options:	70°, 80°, 90° (Swing in/out & HD lift off) ELR, ALR
Warranty	Lifetime on frame and cross braces, 1 year on parts

^{*} Actual measurements may vary according to the wheelchair's configuration and accessories

