



Full-Electric Bed

Item #	Product Description	UOM
15005	Full-Electric Bed	1/cs
15005BV-FR	Full-Electric Bed with Full-Length Side Rails (item #15001ABV)	1/cs
15005BV-HR	Full-Electric Bed with Half Rails (item #15201BV)	1/cs
15005BV-PKG	Full-Electric Bed with Full-Length Side Rails (item #15001ABV) and 80" Innerspring Mattress (item #15006)	1/cs
15005BV-PKG-T	Full-Electric Bed with Full-Length Side Rails (item #15001ABV) and 80" Therapeutic Support Mattress (item #15019)	1/cs
15005BV-PKG-1	Full-Electric Bed with Half Rails (item #15201BV) and 80" Innerspring Mattress (item #15006)	1/cs
15005BV-PKG-1-T	Full-Electric Bed with Half Rails (item #15201BV) and 80" Therapeutic Support Mattress (item #15019)	1/cs
15005BV-PKG-2	Full-Electric Bed with Full-Length Side Rails (item #15001ABV) and 80" Foam Mattress (item #15007)	1/cs

- Quiet, smooth operation.
- The easy-to-use hand pendant electrically controls the positioning of the upper body, knees and overall height of the bed frame to provide maximum comfort and versatility for the patient or caregiver.
- Side clips store shaft when bed is unassembled.

specifications

Bed Height	16" - 24.5"
Unit Dimensions	36" (W) x 88" (L)
Weight	181 lbs.
Carton Shipping Weight	194 lbs.
Weight Limit	450 lbs.
Warranty	(Original Purchaser) Limited Lifetime on Welds, 5 Years on Frame, 1 Year on All Other Parts and Components

- Bed ships in two cartons. Easy to set up.
- Head and foot adjustment allow for an anatomically correct sleep surface.
- Channel frame construction provides superior strength and reduced weight.
- Reinforced frame resists twisting and bending.
- Zinc-coated spring deck.
- Heavy-duty bed ends are attractive in any decor and easy to maintain.
- Head and foot ends are interchangeable with Invacare® and Sunrise®.
- Half rails and full rails available. (Sold separately) (Figure A - B)
- UL 73 approved. CSA approved.
- Manufactured in ISO9001 facility.