



# Invacare® Top End® Pro™-2 All Sport Wheelchair

The Invacare® Top End® Pro™-2 All Sport Wheelchair is all about sports. If you want to get into the game, this is the chair for you. Just name it: basketball, tennis, softball, table tennis, badmitton, floor hockey, racketball, lacrosse or any other sport that is played on a court and the Top End Pro-2 Wheelchair can get you started. The Pro-2 All Sport Wheelchair is offered in five seat widths with lots of adjustments which is perfect for clubs and aspiring athletes. Get in the game.

# Model TE10014



10<sup>11</sup>

Invacare Top End Pro-2 All Sport Wheelchair

10" Seat Width Neon Green "Mini"

# **Invacare Top End Pro-2 All Sport Wheelchair Features**

- Choice of 5 seat widths 10",12",14", 16" or 18"
- Sturdy, lightweight aluminum frame
- Adjustable seat depth:10" & 12" wheelchair seat depth 10" & 12"14", 16", 17" wheelchair seat depth 13", 15", 17"
- Screw on seat upholstery
- Back Height 11"-14"
- Adjustable front and rear seat heights
- Adjustable clamp-on tubular recessed footrest
- Adjustable center of gravity
- Removable single swivel anti-tip bar
- Camber: 15°
- 25" High performance wheels with high pressure clincher tires

- Frame Adjustable axle position 2.5"-5"
- Quick release axles
- Non-folding adjustable back: angle and height
- Adjustable tension back upholstery
- Aluminum side guards
- Cushion
- Positioning strap
- Rear frame protectors
- 250 lb. weight capacity



Invacare Top End Pro-2 All Sport Wheelchair

12" Seat Width Safety Yellow "Prep"

121

Invacare Top End Pro-2 All Sport Wheelchair

14" Seat Width White Snow "Junior Varsity"

14!





Invacare Top End Pro-2
All Sport Wheelchair
18" Seat Width

18" Seat Width Simple Blue "Adult" Size: 15021



Invacare Top End Pro-2 All Sport Wheelchair

16" Seat Width Simple Red "Varsity" Size: 15017











- 1. Welded bumper/recessed footrest
- 2. Adjustable seat depth: 13" 17"
- 3. Removable single swivel anti-tip bar
- 4. Adjustable seat height: front and rear
- 5. Optional competition basketball wing

# **Options**

Basketball carrying bag Click straps

Note: Some assembly is required.

# **Invacare Corporation**

www.invacare.com

#### USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

#### Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

### Top End

www.topendwheelchair.com Sports and Recreation Products 4501 63rd Circle North Pinellas Park, Florida 33781-5914 (800) 532-8677 Fax: (727) 522-1007

## ©2018

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®.
All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted.
Specifications are subject to change without notification.
Form No. 18-324
181126

