



Tote

Made in the U.S.A!

Model 003



WheelChair Carrier
7325 Douglas Road
Lambertville, MI 48144
www.wheelchaircarrier.com
(800) 541.3213
admin@wheelchaircarrier.com

- Low cost wheelchair carrier
- Lifting the wheelchair onto the carrier is required

Specifications and Features:

- For folding wheelchairs up to 15 1/2 inches wide
- Does NOT tilt for loading
- Does NOT fold for storage
- Attaches to any ball mount

Capacity: 75 lbs Weight: 16 lbs