



# Introducing New Additions to ONTRAC®




Without Handrail

- Ramps available in lengths of 3', 4', 5', 6', 7', 8', 9' & 10'
- Choice of ramps with or without handrail
- Installs in minutes
- Durable and weather resistant
- Features 2" curbing
- Slip resistant grooved aluminum pathway surface with 4-WayTraction
- All aluminum construction
- Little or no maintenance
- ADA extended handrail available upon request
- Optional Powder Coated Handrail Available



Now ALL Lengths are Deliverable by UPS or FedEx!



4-Way Traction



Shown with Handrail

Item Number (with handrails)	Length	Width	Weight	Capacity*
XPS336	3'	36"	49 lbs.	850 lbs.
XPS436	4'	36"	59 lbs.	850 lbs.
XPS536	5'	36"	69 lbs.	850 lbs.
XPS636	6'	36"	88 lbs.	850 lbs.
XPS736	7'	36"	98 lbs.	850 lbs.
XPS836	8'	36"	108 lbs.	850 lbs.
XPS936	9'	36"	117 lbs.	850 lbs.
XPS1036	10'	36"	127 lbs.	850 lbs.

\* Maximum ramp capacity: 425 lbs. (one axle) – 850 lbs. (two axles).

Item Number	Length	Width	Weight	Capacity*
XPS336NR	3'	36"	25 lbs.	850 lbs.
XPS436NR	4'	36"	33 lbs.	850 lbs.
XPS536NR	5'	36"	42 lbs.	850 lbs.
XPS636NR	6'	36"	53 lbs.	850 lbs.
XPS736NR	7'	36"	61 lbs.	850 lbs.
XPS836NR	8'	36"	69 lbs.	850 lbs.
XPS936NR	9'	36"	76 lbs.	850 lbs.
XPS1036NR	10'	36"	84 lbs.	850 lbs.

\* Maximum ramp capacity: 425 lbs. (one axle) – 850 lbs. (two axles).



Prairie View Industries

• Phone 800.554.7267(RAMP) • info@pviramps.com • www.pviramps.com •