



NEW & IMPROVED

The Kanga 2G

Growable Tilt-in-Space Wheelchair

FEATURES

- Newly designed, reinforced frame
- New rehab style seating components with multiple adjustments (**sold separately**)
- Newly designed push button height and depth adjustable, removable "T" style armrests
- Growable frame width - no added components needed
- Folding frame tilt-in-space wheelchair with 45° of tilt
- Height & angle adjustable back canes
- 8" height and angle adjustable front casters and 12" quick-release rear wheels
- Independently depth adjustable footplates
- Removable foot hangers that swing inward and outward for safe transfer
- Removable anti tippers included
- Crash tested in accordance with ISO 7176-19 for use in motor vehicles (transit brackets included)


For children and young adults in need of minimal to moderate support, The Kanga wheelchair offers tilt-in-space adjustability to maintain a proper sitting posture. The Kanga wheelchair grows from 10" to 14" width and offers 45° of tilt for pressure relief and gravity assisted positioning. For use with Inspired by Drive newly available complex rehab seating option. Also accommodates other manufacturer's seating.



Tilt & Recline



Folds for Easy Transport



Adjustable Width Crossbrace



Standard Transit Brackets & Anti-Tipper



Seating system sold separately

ACCESSORIES

See our Poziform brochure for more positioning product offering!



Rehab Style Seating System



Angle Adjustable Footplates



Elevating Legrests



20" Rear Wheel Assembly



Poziform Contoured Chest Harness



Poziform Padded Pelvic Belt

ORDER FORM

EMAIL FORM

CLEAR FORM

PRINT FORM

Dealer: _____ Marked For: _____ Account #: _____

Date: _____ PO#: _____ Ordered By: _____

Phone: _____ Fax: _____ E-mail: _____

Ship to address: _____

KANGA WHEELCHAIR FRAME	\$2,845	<input type="radio"/> KGP-8010-2GR	<input type="radio"/> KGP-8012-2GR	<input type="radio"/> KGP-8014-2GR
8" casters, 12" wheels, 3.5" x 8" fixed angle footplates				
Weight Capacity		200 lb	200 lb	200 lb
Seat Width		10"	12"	14"
Seat Rail Depth (from front of back cane)		14" - 18"	14" - 18"	14" - 18"
Seat Rail Height		17.25"	17.25"	17.25"
Seat Rail to Footplate Length (leg length)		5.5" - 11.5"	5.5" - 11.5"	5.5" - 11.5"
Back Cane Height (from seat rail)		22" - 26"	22" - 26"	22" - 26"
Back Angle		85° - 113°	85° - 113°	85° - 113°
Armrest Height (seat rail to top of armrest)		8" - 12"	8" - 12"	8" - 12"
Tilt Angle		0° - 45°	0° - 45°	0° - 45°
Overall Length (w/o footrests)		29"	29"	29"
Overall Height		38.5" - 42.5"	38.5" - 42.5"	38.5" - 42.5"
Overall Width (outside rear wheels)		17.5"	19.5"	21.5"
Folded Dimensions		29" x 27" x 13.5"	29" x 27" x 13.5"	29" x 27" x 13.5"
Product Weight		40 lb	40 lb	40 lb
HCPCS		E1234	E1234	E1234

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims. Referenced HCPCS codes are also not guarantee of payment.

SEATING & ACCESSORIES (all accessories come with mounting hardware)				
Seat & Back Package w/Mounting Hardware Contoured I Back / Contoured Seat (E2609/E2617)	\$1,595	<input type="radio"/> 1300-921110 seat depth 10"-12"	<input type="radio"/> 1300-921112 seat depth 12"-14"	<input type="radio"/> 1300-921114 seat depth 14"-16"
Comfort Max Headrest w/Offset Hardware (E0955/E1028)	\$315	<input type="radio"/> 4008-401010	<input type="radio"/> 4008-401012	<input type="radio"/> 4008-401014
Abductor Pommel, Flip Down (E0957/E1028)	\$325	<input type="radio"/> 11060-10110	<input type="radio"/> 11060-10112	<input type="radio"/> 11060-10114
Thigh Supports, QR, Straight & Offset (pr) (E0957/E1028)	\$650	<input type="radio"/> 11072-34410	<input type="radio"/> 11072-34412	<input type="radio"/> 11072-34414
Trunk Support, QR, Adjustable (pr) (E0957/E1028)	\$595	<input type="radio"/> 3040-021010	<input type="radio"/> 3040-021012	<input type="radio"/> 3040-021014
Poziform Pelvic Belt, Padded (K0108)	\$85	<input type="radio"/> PPB-100PBXS	<input type="radio"/> PPB-100PBS	<input type="radio"/> PPB-100PBM
Poziform Contoured Chest Harness w/Cam Locks (K0108)	\$145	<input type="radio"/> PCH-340XS	<input type="radio"/> PCH-340S	<input type="radio"/> PCH-340M
Poziform Calf Panel, Soft Padded Nylon (K0108)	\$85	<input type="radio"/> PCP-280S	<input type="radio"/> PCP-280M	<input type="radio"/> PCP-280L
Poziform Ankle Positioners, Side Release (K0108)	\$130	<input type="radio"/> PAP-200XS	<input type="radio"/> PAP-200S	<input type="radio"/> PAP-200M
Upper Extremity Support, Polycarbonate (K0950/K018)	\$430	<input type="radio"/> 12011-10	<input type="radio"/> 12011-12	<input type="radio"/> 12011-14
Angle Adjustable Footplates (pair) (K0052/K0040)	\$205	<input type="radio"/> KGP-FR-AA	<input type="radio"/> KGP-FR-AA	<input type="radio"/> KGP-FR-AA
Elevating Legrests (pair) (E0990)	\$235	<input type="radio"/> KGP-ELR-2G	<input type="radio"/> KGP-ELR-2G	<input type="radio"/> KGP-ELR-2G
20" Rear Wheel Assembly (pair) (K0108)	\$205	<input type="radio"/> KG-RW20	<input type="radio"/> KG-RW20	<input type="radio"/> KG-RW20