



Item #	Seat	Color	UOM
SFSCOUT4	16.5" Folding Seat	Red ■	1/cs



HCPCS
K0800



Spitfire Scout 4 Compact Travel Scooter

- Quick and easy disassembly. (Figure A)
- Color-through panels to hide scuffs and scratches.
- Armrests are padded and adjustable.
- Quick connect batteries.
- Ergonomic throttle control.
- Height-adjustable swivel seat with fold-down backrest.
- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Delta tiller.
- Easy freewheel operation

specifications

HCPCS	K0800	Charger	2A Offboard
Weight Capacity	300 lbs.	Brakes	Electromagnetic
Top Speed	4.25 mph	Front Wheels	8" x 2" Flat-free
Maximum Range	9 miles	Rear Wheels	8" x 2" Flat-free
Turning Radius	53.75"	Base Weight	58 lbs.
Climbing Angle	6°	Battery Weight (pair)	20 lbs.
Ground Clearance	2.5"	Seat Weight	16 lbs.
Dimensions	42.5" (L) x 20.5" (W)	Total Weight	94 lbs.
Floor-to-Seat Height	17" - 19"	Heaviest Piece Weight (Front Section)	39 lbs.
Seat Dimensions	16.5" (W) x 16" (D)	Warranty on Frame	Lifetime
Controller	CTE 45A	Warranty on Electronic Controller/Drive Train Components	14 months
Motor	24V x 270W x 4700 rpm	Warranty on Batteries	6 months
Batteries (pair)	12V x 12AH		