

SENTRA FULL RECLINING



**THE 22" SEAT WIDTH
FEATURES A WEIGHT
CAPACITY OF 450 LBS.**



Deluxe Sentra Full Reclining (K0007)

**COMES STANDARD
WITH SWING-AWAY
ELEVATING LEGRESTS.**

FEATURES:

- New state-of-the-art hydraulic reclining mechanism allows for infinite adjustments up to 180°.
- Carbon Steel frame with triple coated chrome for an attractive, chip-proof, maintainable finish.
- Durable, heavy gauge naugahyde upholstery.
- Padded armrests provide added patient comfort.
- Composite Mag-style wheels with chrome hand rims are lightweight and maintenance free.
- Comes standard with swing-away elevating legrests.
- Wheels set back on frame prevents tipping.
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability.
- Standard dual axle provides easy transition of seat to hemi-level (excluding 22").
- Aluminum footplates are attractive, crack-proof and lightweight.
- Front caster forks are adjustable in 3 positions.
- **Carry pocket standard.**
- **Rear anti-tippers standard.**
- 8" front casters.
- Headrest extension with cushioned head immobilizer standard.
- Comes with push to lock wheel locks.
- Black upholstery, chrome frame
- Limited Lifetime Warranty.

Due to product improvements and changes, all images in this catalog are subject to change.

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*