

Blue Streak Wheelchair
Single Axle



HCPCS
K0001

COMPETITIVE
EDGE Line

Item #	Seat Width	Armrest	Product Front Riggings	UOM
BLS18FBD-SF	18"	Flip-Back Desk Arm	Swing-away Footrests	1/cs
BLS18FBD-ELR	18"	Flip-Back Desk Arm	Swing-away, Elevating Leg rests	1/cs
BLS20FBD-SF	20"	Flip-Back Desk Arm	Swing-away Footrests	1/cs
BLS20FBD-ELR	20"	Flip-Back Desk Arm	Swing-away, Elevating Leg rests	1/cs

- Blue Powder-coated frame, black cross brace.
- Flip-back, padded desk-length arm.
- Nylon upholstery is durable, lightweight, attractive and easy to clean.
- Solid rubber tires mounted on composite wheels provide durability and low maintenance.
- Comes with push-to-lock wheel locks.
- Comes with calf strap.
- Comes with swing-away footrests or elevating leg rests.
- Weight Capacity: 250 lbs.

specifications

Width Open with 18" Seat	24"	Seat Dimensions with 18" Seat	18" (W) x 16" (D)
Width Open with 20" Seat	26"	Seat Dimensions with 20" Seat	20" (W) x 16" (D)
Width Closed	12.5"	Length with Riggings	42"
Height	36"	Product Weight	41.25 lbs.
Seat-To-floor Height	19.5"	Carton Shipping Weight	47 lbs.
Back Height	16"	Warranty	Limited 3 Year
Weight Limit	250 lbs.		

accessories

Anti-Tippers

STDS833	Anti-Tippers With Wheels	1/pr
STDS833N	Anti-Tippers With Wheels	1/pr

Seat Belts

STDS850	Seat Belts, Auto-Clasp Type Closure	1/ea
STDS851	Seat Belts, Velcro® Type Closure	1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea
STDS856	Seat Belts, Velcro® Type Closure, Bariatric	1/ea

Wheel Lock Extension

STDS801GT	Wheel Lock Extension	1/ea
-----------	----------------------	------

Anti-Fold Bar

STDS806	Anti-Fold Bar	1/ea
---------	---------------	------

Wheelchair Cane/Crutch Holder

STDS1034	Wheelchair Cane/Crutch Holder	1/ea
----------	-------------------------------	------

Heel Loops

STDS831	Heel Loops	1/pr
---------	------------	------

Limb Support

WASR	Limb Support, Right	1/ea
WASL	Limb Support, Left	1/ea