Jazzy[®] EVO 613 SERIES

THE MOST **ADVANCED** JAZZY[®] EVER

You will love the all-new Medicarereimbursable Jazzy® EVO 613 series! You can choose between the EVO 613 or the EVO 613 Li, which uses a lithium iron phosphate battery for longer life, more convenient charging and fewer service calls.



13"
WHEELS



Get the lithium advantage with the Jazzy® EVO 613 Li.

> FDA Class II Medical Device¹



Scan QR code to learn more.













PRIDE

EVO 613 SERIES

AVAILABLE MODELS

Jazzy® EVO 613 / Jazzy® EVO 613 Li Group 2, Standard (K0822/K0823)

FDA Class II Medical Device¹





Lithium Iron Phosphate Battery

(Available on Jazzy® EVO 613 Li)

- Faster charging time
- Sustained power performance
- Lighter weight
- More useable energy
- Longer life

FEATURES

- Maximum speed up to 5.1 mph
- 13" mid-wheel drive tires
- Narrow base width of only 22"
- Group 2, Standard (K0822/K0823)
- Choice of six contemporary removable, matte-finish color-through shrouds
- Depth adjustable, limited-recline, memory foam, high-back with single seat post allows the user to swivel
- Pride's patented Active-Trac® (ATX) Suspension
- Mid-Wheel 6® Technology

ACCESSORIES

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag kit
- Oxygen tank holder
- Rear basket
- Saddlebag
- Walker holder
- Weather cover
- 2.1A XLR USB charger

OPTIONS

- Elevated leg rests
- Swing-away leg rests
- Residual limb support
- Heel loops
- Unoccupied transit kit
- Depth-adjustable solid seat pan, limited recline, high-back with swivel

SPECIFICATIONS

WEIGHT CAPACITY 300 lbs.

MAXIMUM SPEED² Up to 5.1 mph

DIMENSIONS³ Width: 22" Length: 41" with footrest

TOTAL WEIGHT EVO 613: 189 lbs. EVO 613 Li: 160.5 lbs.

TURNING RADIUS 22"

RANGE PER CHARGE³ EVO 613: Up to 14 miles @ 300 lbs.

EVO 613: Up to 17.1 miles @ 200 lbs. EVO 613 Li: Up to 17.6 miles @ 300 lbs. EVO 613 Li: Up to 21 miles @ 200 lbs.

WARRANTY

FRAME: 5-year limited
DRIVETRAIN: 2-year limited
ELECTRONICS: 2-year limited

SLA BATTERY: 13-month limited LITHIUM BATTERY: 3-year limited

- $1.\ Pride^{@}\ FDA\ Class\ II\ Medical\ Devices\ are\ designed\ to\ aid\ individuals\ with\ mobility\ impairments.$
- 2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or 3%).



