

Cirrus Plus LT

HCPCS: K0822



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

- Foldable, lightweight frame for easy portability
- Attractive steel cross-brace frame
- Armrests are padded, adjustable and removable
- Adjustable-length controller mount
- Adjustable-tension back upholstery
- Swing-away footrests with heel loops
- Positioning belt included
- Padded back upholstery with storage pouch and cushioned seat
- Programmable controller
- Flat-free, non-marking tires
- Easy free-wheel operation
- Manual wheel locks
- General foam seat cushion included
- Anti-tip wheels

PRODUCTS

Description

CP18FBAN	CirrusPlusPowerChair,18 FlpBkArm,1/cs
CP20FBAN	CirrusPlusPowerChair,20 FlpBkArm,1/cs
CP22FBAN	CirrusPlusPowerChair,22 FlpBkArm,1/cs

SPECIFICATIONS

Description

Anti-Tip Wheels	Yes
Base Weight	90 lbs
Batteries (pair)	12V x 33AH
Battery Weight (pair)	47 lbs
Brakes	Electromagnetic
Casters	8" x 2" solid
Charger	4A Offboard
Climbing Angle	6 degrees
Color	Black
Controller	PG Vsi 50A
Drive Wheels	12.5 x 2.25 solid PU tire
Freewheel Mode	Yes
Ground Clearance	4"
Heaviest Piece	Frame: 90 lbs
Material	Steel

SPECIFICATIONS	Description
Max Speed	5 mph
Maximum Range	12.5 miles
Motor	2 x (24V 320W x 4600 rpm)
Number of Wheels	4
Product Weight	Without Battery: 90 lbs
Seat Depth	16.5"
Seat To Floor Height	21.5"
Seat Width	18"
Turning Radius	30"
Warranty	Batteries: 12 Months Drive Train Components: 24 Months Electronics: 24 Months Frame: Limited Lifetime
Weight Capacity	300 lbs