

## AL825 Pickup 225 Speed



### FEATURES

- 2x the lifting speed
- Delivers chair close to driver or passenger side door
- Powered lifting and rotation
- Only two motors to control with remote hand control
- Mounts on either side of vehicle
- Docking device included
- New Dash Mounted Control Panel

### SPECIFICATIONS

Lifting Capacity:	225 lbs.
Operation:	Fully Powered
Adjustable Boom Height:	34"-41"
Boom Clearance:	19"-26"
Adjustable Boom Length:	28.75"-46.5"
Adjustable Boom Pitch:	0-70 degrees
Adjustable Rotation Limit:	Yes
Warranty:	3-Year Transferable

### DESCRIPTION

The AL825 is the perfect solution for pickup truck owners that don't need the 350 lbs. capacity of the AL835 and want a faster cycle time. The AL825 provides ample clearance and reach with one less motor than the competition, meaning fewer buttons and simpler operation.

This model is equipped with a newly designed anti-rotate strap system and compact base for a better fit in the truck bed. This convenient lift brings the manual wheelchair, scooter or power chair right to the driver's or passenger's side door. The AL825 Pickup 225 Speed Lift easily loads virtually any mobility device available into the bed of your pickup truck and unload it to the driver or passenger side door.

© 2012 Harmar MKT-SS-0512