







For new players, adjustability is key. The **Invacare® Top End® Pro™ Tennis Wheelchair** has just that, featuring adjustable front and rear seat heights, center-of-gravity, footrest angle/height, backrest angle/height and anti-tip. Plus, the Pro Tennis Wheelchair is available in two frame height choices: short or tall. The short frame is ideal for children and the tall frame is for adults, which makes this chair a great choice for programs that have multiple users. Adjusts in minutes, super easy to use, lightweight and durable.

Features

Adjustable front and rear seat heights Anti-tip/5th wheel: adjustable height Fixed back Adjustable back height Adjustable tension back upholstery Adjustable center-of-gravity Fixed camber 15°, 18° or 20° High performance Sun® wheels: 24", 25", 26" Anodized pushrims High pressure clincher tires Smart fork caster system 72mm casters Adjustable platform footrest Powder coat paint 250 lb. weight capacity

Model #: SPTN



- 1. Front seat height and center-of-gravity adjustments
- 2. Rear seat height adjustment
- 3. Quick release lever
- 4. Single swivel fixed fifth wheel
- 5. Precision fork caster system
- 6. Adjustable footrest system

Options

- Clamp-on fixed aluminum side guards Spinergy[®] SLX wheels Sun[®] Equalizer[™] wheels 12 Tab pushrims Kenda[®] tires Click straps O-ring positioning belt
- Spoke guards Frame protectors Alignment gauge Telescoping anti-tip bar Cushion Folding Back Fold-in side guards

Invacare Corporation www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

Top End

www.topendwheelchair.com Sports and Recreation Products 4501 63rd Circle North Pinellas Park, Florida 33781-5914 (727) 522-8677 (800) 532-8677 Fax: (727) 522-1007

©2016

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols [™] and [®]. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Spinergy is a registered trademark of Spinergy Inc. Sun and Equalizer are registered trademarks of Sun Components. Kenda is a registered trademark of Kenda Rubber Industrial Co., Ltd. Specifications are subject to change without notification. Form No. 12-507 Rev. 03/16 151185

