



## General Use Back Cushion with Lumbar Support

Item #	Product Description	Product Dimensions	Product HCPCS	UOM
14906 (Replaces Mason # 8030)	Wheelchair Back Cushion	16" (W) x 17" (D) x 2.5" (H)	E2611	1/ea
14889 (Replaces Mason # 8031)	Wheelchair Back Cushion	18" (W) x 17" (D) x 2.5" (H)	E2611	1/ea
14920 (Replaces Mason # 8032)	Wheelchair Back Cushion	20" (W) x 17" (D) x 2.5" (H)	E2611	1/ea
8033	Wheelchair Back Cushion	22" (W) x 19" (D) x 2.5" (H)	E2612	1/ea

- Memory-foam back cushion provides support, comfort and a tactile curve to improve posture.
- Durable Masonite board secured on back of foam core prevents hammocking in sling-back wheelchairs.
- Conforms to lumbar curve to support lower back.
- Sealed liner protects foam core from moisture.
- Non-slip vinyl back prevents the cushion from sliding.

### specifications

Cover Material	Half P/U, Half Vinyl
Weight Limit	N/A
Warranty	18 Months Limited