



Made in the U.S.A!

Lift 'n' Go Electric Lift

Model 210



Rated "BEST BUY" in Vehicle Lifts according to Consumers Digest! Our standard platform design is a popular choice which universally accommodates both scooters and power chairs.

Specifications and Features:

- 28" X 48" of usable platform
- 24" of possible ground clearance for smooth traveling
- Automatically folds up to stowaway position
- Switch box allows for easy operation of the lift
- Includes 2 retractable ratchets for securing the chair
- Includes wire harness with 30 amp resettable circuit breaker
- Ask about our Power Booster Box if you are transporting a heavier mobility device OR your battery is weaker



WheelChair Carrier
7325 Douglas Road
Lambertville, MI 48144
www.wheelchaircarrier.com
(800) 541.3213
admin@wheelchaircarrier.com