



# PATHWAY® Classic Series RAMP

AVAILABLE WITH OPTIONAL HANDRAILS



Pathway shown with optional handrails

Designed for safety and engineered for reliability, the PATHWAY® Classic Series portable wheelchair ramp is a smart alternative to costly wooden ramps. When mobility needs extend beyond the scope of a conventional portable ramp, this ramp is your answer. Can be easily moved or relocated as needs change. Designed for the rigors of repeated, everyday use in all weather conditions. Made of aluminum and features an extruded non-skid surface. Select sizes are available with optional handrails.

Additional sizes available.



Pathway without handrails



SPECIFICATIONS	PATHWAY CS4	PATHWAY CS6	PATHWAY CS8	PATHWAY CS10
Made of aluminum	•	•	•	•
Non-skid surface	•	•	•	•
Usable size (L x W)	4' x 36"	6' x 36"	8' x 36"	10' x 36"
Outside width	37"	37"	37"	37"
Ramp weight (lbs.) *	45	60	75	90
Weight capacity (lbs.)	850	850	850	850

SPECIFICATIONS (with handrails)	PATHWAY CSHR4	PATHWAY CSHR6	PATHWAY CSHR8	PATHWAY CSHR10
Made of aluminum	•	•	•	•
Non-skid surface	•	•	•	•
Usable size (L x W x H)	4' x 36" x 35¾"	6' x 36" x 35¾"	8' x 36" x 35¾"	10' x 36" x 35¾"
Outside width	40¼"	40¼"	40¼"	40¼"
Ramp weight (lbs.) *	70	85	105	130
Weight capacity (lbs.)	850	850	850	850

Additional sizes are available upon request. Please call for details.

\* Weights are approximate and may vary.

Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.



Bariatric weight capacity = 800 lbs. or greater



Bariatric width = 36" or greater