

Cruiser III Wheelchair

Lightweight, Dual-Axle Wheelchair



HCPCS
K0003

All the features of the original
Cruiser III and the New Built-in
Seat Rail System and
Extendable Upholstery


STEP 1



STEP 2



STEP 3



Item #	Seat Width	Armrest	Product Front Riggings	UOM
K316DDA-SF	16"	Flip Back Detachable Desk Arm	Swing-away Footrest	1/cs
K316DDA-ELR	16"	Flip Back Detachable Desk Arm	Elevating Leg rest	1/cs
K316DFA-SF	16"	Flip Back Detachable Full Arm	Swing-away Footrest	1/cs
K316DFA-ELR	16"	Flip Back Detachable Full Arm	Elevating Leg rest	1/cs
K316ADDA-SF	16"	Flip Back Detachable and Adjustable Height Desk Arm	Swing-away Footrest	1/cs
K316ADDA-ELR	16"	Flip Back Detachable and Adjustable Height Desk Arm	Elevating Leg rest	1/cs
K318DDA-SF	18"	Flip Back Detachable Desk Arm	Swing-away Footrest	1/cs
K318DDA-ELR	18"	Flip Back Detachable Desk Arm	Elevating Leg rest	1/cs
K318DFA-SF	18"	Flip Back Detachable Full Arm	Swing-away Footrest	1/cs
K318DFA-ELR	18"	Flip Back Detachable Full Arm	Elevating Leg rest	1/cs
K318ADDA-SF	18"	Flip Back Detachable and Adjustable Height Desk Arm	Swing-away Footrest	1/cs
K318ADDA-ELR	18"	Flip Back Detachable and Adjustable Height Desk Arm	Elevating Leg rest	1/cs
K320DDA-SF	20"	Flip Back Detachable Desk Arm	Swing-away Footrest	1/cs
K320DDA-ELR	20"	Flip Back Detachable Desk Arm	Elevating Leg rest	1/cs
K320DFA-SF	20"	Flip Back Detachable Full Arm	Swing-away Footrest	1/cs
K320DFA-ELR	20"	Flip Back Detachable Full Arm	Elevating Leg rest	1/cs
K320ADDA-SF	20"	Flip Back Detachable and Adjustable Height Desk Arm	Swing-away Footrest	1/cs
K320ADDA-ELR	20"	Flip Back Detachable and Adjustable Height Desk Arm	Elevating Leg rest	1/cs

- Built-in seat rail extensions and extendable upholstery easily adjust seat depth from 16" to 18" (Figure A - C)
- Weighs under 36 lbs. (excluding front riggings).
- Carbon steel frame with silver vein finish.
- Removable flip-back arms allow for easy transfer. (Figure D)
- New frame style eliminates seat guides and allows for custom back inserts and accessories.
- Nylon upholstery is durable, lightweight, attractive and easy to clean.
- Composite, Mag-style wheels are lightweight and maintenance free.
- 8" front casters are adjustable in three positions.
- Padded armrests provide additional comfort.
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings. (Figure E)
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability.
- Dual axle provides easy transition of seat height to hemi-level.
- Comes with push-to-lock wheel locks.

specifications

Width Closed	12"	Product Weight (with 18" Seat)	36 lbs.
Height	35"	Product Weight (with 20" Seat)	38 lbs.
Seat--To-Floor Height	17.5" - 19.5"	Weight Limit (with 16" and 18" Seat)	300 lbs.
Back Height	16"	Weight Limit (with 20" seat)	350 lbs.
Seat Depth	16" - 18"	Carton Shipping Weight	43 lbs. or less
Length with Riggings	42"	Warranty	Limited Lifetime
Product Weight (with 16" Seat)	35 lbs.		

Universal Forward Stabilizer

STDS860	Universal Forward Stabilizer	1/pr
---------	------------------------------	------

Anti-Tippers

STDS819	Anti-Tippers With wheels	1/pr
---------	--------------------------	------

Seat Belts

STDS850	Seat Belts, Auto-Clasp Type Closure	
STDS851	Seat Belts, Hook and Loop Fastener Type Closure	1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric	1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric	1/ea

Telescoping I.V. Pole Attachment

STDS820	Telescoping IV Pole Attachment	1/ea
---------	--------------------------------	------

Universal Oxygen "E" Cylinder Holder

STDS803	Universal Oxygen "E" Cylinder Holder	1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment	1/ea

6" Wheel Lock Extensions

STDS801	Wheel Lock Extension	1/ea
---------	----------------------	------

Overhead Anti Theft Device

STDS823	Overhead Anti-Theft Device	1/ea
STDS821	Overhead Anti-Theft Device With IV Hooks	1/ea
STDS834	Overhead Anti-Theft Device Single Pole	1/ea

Anti-Fold Bar

STDS806	Anti-Fold Bar	1/ea
---------	---------------	------

Wheelchair Cane/Crutch Holder

STDS1034	Wheelchair Cane/Crutch Holder	1/ea
----------	-------------------------------	------

Heel Loops

STDS831	Heel Loops	
---------	------------	--

Limb Support

WASR	Limb Support, Right	1/ea
------	---------------------	------

Chart Carry Pocket

STDS835	16", 18" and 20" Wheelchairs w/Nylon Upholstery	1/ea
---------	---	------

Super Hemi-Kit

STDSSHK-K3	Super Hemi-Kit	1/ea
------------	----------------	------