Splash!™ Series



This fixed lifting system is designed for swimming facilities seeking to provide user-friendly access to their swimming pools.

- · Third-party tested & verified ADA compliant
- 400 lb/181.5 kg lifting capacity*
- LiftOperator® Intelligent Control includes a 24 Volt rechargeable battery
- · 344° rotation helps ensure a safe transfer area
- Powder-coated stainless steel and aluminum construction
- * Lifting capacity of the Extended Reach model is 300 lb/136 kg



Lift Color



Due to printing technology actual color may differ.

Splash! variations include:

- Hi/Lo Useful for facilities that have both in-ground pools and raised spas
- Spa Designed for raised spas
- Extended Reach For use on pools with unusual gutter configurations requiring a longer reach
- ER Hi/Lo For facilities that have both in-ground pools and raised spas, requiring a longer reach than the Splash Hi/Lo


			Shipping – (Class 85)			
Model No.	Description	Weight	Length	Width	Height	
300-0000	Splash!	274 lbs	40"	48"	57"	
350-0000	Splash! Hi/Lo	280 lbs	40"	48"	57"	
370-0000	Splash! Extended Reach	282 lbs	40"	48"	75"	
375-0000	Splash! Spa	268 lbs	40"	48"	57"	
385-0000	Splash! ER Hi/Lo	281 lbs	40"	48"	88"	
390-0000	Splash! California Package*	285 lbs	40"	48"	75"	
390-2000	Splash! ER California Package*	293 lbs	40"	48"	75"	

Notes:

- Includes battery, charger, battery console cover, water-resistant hand control, stainless steel anchor socket with cover, footrest and seat belt assembly
- To order lifts with arm rests, drop the 0 on the end of the part number and add a 5 example: 300-0000 becomes 300-0005 for a Splash! with arm rests
- To order a lift with an Activation Key (for LiftOperator), add a K to the end of the part number
- To order without an anchor jig add an N to end of part number
- * California package includes arm rests, mast cover, spineboard attachment and caddy

Parts & Accessories:

- 400-0000 Caddy
- 170-1000 Arm Rest Assembly
- 920-5000 Splash! Cover
- 900-4000 Seat Pad
- 970-0000 Seat Saver Cover
- 500-1000 Spineboard Attachment
- 900-2000 Stability Vest
- 300-6200L Locking Anchor & Cap
- 1001495 Battery



Seat is 18" wide

